

Hand-bikes

EASYBIKE

Technical features

- Steel frame.
- 24" front wheel equipped with high pressure tyre.
- Tilt and width adjustable fixing arms.
- Adjustable stand.
- Front light and bell.

EASYBIKE 8CR

Revo 8 gear hub shift with coaster brake.
Available also for kids.

EASYBIKE 8F

8 gear hub shift with lever brake.
Available also for kids.

EASYBIKE 27

27 gear shift and twin brake with controls on handles.

Tetraplegic handle



Front light



Cyclo computer



For safety reasons, Easybike can be fitted only to wheelchairs with not detachable foot-rests.

Order form:

- EASYBIKE KID 8F
 EASYBIKE KID 8CR

- EASYBIKE 27
 EASYBIKE 8CR
 EASYBIKE 8F

-Model of wheelchair to be fitted to Easybike :

-Outer width of the frame front tube: (L) cm

-User height: cm



- Frame colour: metallic grey
 glossy black

metallic blue

SET-UP NOTES

To find the right unicycle inclination, the fixing arms must be turned in a symmetrical way considering the tooth mark (A) referred to the groove (B).

Easybike should be adjusted setting the crankset at about 5 cm from legs and 15 cm from chest.



Options

- Cyclo computer:
Rear-view mirror:
Tetraplegic handles:

Other frame colour:

(see "Frame colour" leaflet)



Via dell' Artigianato II, 29
35010 Villa del Conte (PD) Italy
Tel. +39 049 9325733 - Fax +39 049 9325734
www.offcarr.com - offcarr@offcarr.com

Dealer address

Date: _____
Order number: _____
Customer ref.: _____
Notes: _____