

Order N° _____ Quotation
 Date _____ Signature _____
 Customer Ref. _____ User weight Kg _____

PAGE 1/2


2020, January

Dealer

ZODIAC

Key standard ○ option

- L** Seat width 34 36 38 40 42 44
- P** Seat depth 36 38 40 42 44
- HS** Back-rest height 27 29 31 33 35 37 39
- HA** Front seat/floor height 46 48 50
- HP** Rear seat/floor height 40 42 44 46
- DP** Foot-rest/seat distance 35 37 39 41 43 45 47 49 cm




Frame colour ①

- Metallic grey Pearl white Matt black
- Glossy black Cyclamen light purple
- (see "Colours" catalogue)
- Neon green Neon orange

Upholstery ②


- Tension adjustable nylon back-rest
- Tension adj. breathable back-rest
- Without back-rest for seating system:
- SAMEKH**
- other seating system:.....

Single tubular foot-plate ③




- Open Closed

Pushing handles ④




- Without pushing handles Standard pushing handles Height adjustable pushing handles

Distance between the hub centre and back-rest ⑤




- 7 cm 13 cm
- 9 cm (standard) 15 cm
- 11 cm

Back-rest tilt ⑥



- 89°
- 91° (standard)
- 93°

Front angle ⑦



- 83° (standard) (A=16cm)
- 75° (A=23cm)

Rear wheel camber ⑧

- 1° 2°
- 3° 4°

Side guards ⑨



- Aluminium straight Aluminium with carbon fibre wing - detachable
- Aluminium with carbon fibre wing

"L" arm-rest, detachable, h 22cm ⑩




- Tip-up

Elbow-rest


- Tubular
- Straight (26 cm long)
- Classic (23 cm long)
- Classic (30 cm long)

Brakes ⑪



- Standard
- Scissor
- Composite

Revolving anti-tip device ⑫



- Right
- Left

Notes

