


Dealer _____

ALHENA



Key	<input type="checkbox"/> standard	<input type="radio"/> option
L Seat width	<input type="checkbox"/> 34 <input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46	
P Seat depth	<input type="checkbox"/> 34 <input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46	
HS Back-rest height	<input type="checkbox"/> 29 <input type="checkbox"/> 31 <input type="checkbox"/> 33 <input type="checkbox"/> 35 <input type="checkbox"/> 37 <input type="checkbox"/> 39 <input type="checkbox"/> 41	
HA Front seat/floor height	<input type="checkbox"/> 42 <input type="checkbox"/> 43 <input type="checkbox"/> 44 <input type="checkbox"/> 46 <input type="checkbox"/> 48 <input type="checkbox"/> 50	
HP Rear seat/floor height	<input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46	
DP Foot-rest/seat distance	<input type="checkbox"/> 33 <input type="checkbox"/> 35 <input type="checkbox"/> 37 <input type="checkbox"/> 39 <input type="checkbox"/> 41 <input type="checkbox"/> 43 <input type="checkbox"/> 45 <input type="checkbox"/> 47	cm

- Frame colour** ①
- Shiny titanium
 - Glossy black
 - Pearl white
 - Avio blue
 - Antique gold
 - Potassium light purple
 - Pearl grey
 - Sage green
 - (see "Colours" catalogue)

- Upholstery** ②
- Tension adjustable nylon back-rest
 - Tension adj. breathable back-rest
 - Without back-rest
 - Without seat

- Upholstery colour** (nylon only) ③
- Black
 - Red Orange
 - Green Light blue
 - Cornflower blue

Not detachable foot-rest ④
(with DP=43 cm and P=40 cm)
(standard)

Open hinged titanium separate tubular foot-plates ⑤

Foot-holding belts

Closed titanium single tubular foot-plate ⑥

Foot-holding belts

Automatic closure foot-plate ⑦

Foot-holding belts

Not detachable foot-rest
(with DP=43 cm and P=40 cm)

(approx. measurement)

Separate aluminium profile tilt adjustable foot-plates ⑧

forward backward

Foot-holding belts

Single aluminium profile tilt adjustable foot-plate ⑨

forward backward

Foot-holding belts

Side guards ⑩

- Aluminium with carbon fibre wing
- Plastic
- Curved aluminium
- Straight aluminium

Detachable sport arm-rest h. 22 cm (Tip-up) ⑪

Elbow-rest

- Tubular
- Classic (23 cm long)
- Classic (30 cm long)

Tip-up height adjustable from 22 to 28 cm "L" arm-rest ⑫
(not detachable, open with sideguards)

Elbow-rest

- Tubular
- Straight (26 cm long)
- Classic (23 cm long)
- Classic (30 cm long)

Rear wheels ⑬

- 24" Ultralight (24 spokes)
- 24" Light (24 spokes)
- 24" Star 6x3 (18 spokes)
- 24" (36 spokes)
- 24" Spinergy
- 24" Mountain-bike
- 22" Light (24 spokes)

Tyres

- Pneumatic Marathon Plus
- Pneum. high press. plugged
- Pneum. high press. smooth
- Solid (not avail. with Ultralight)

Handrims ⑭

- Aluminium silver
- Aluminium dark
- Chromium-plated steel
- Steering (alum.dark) (*)
- Rubber coated
- Tween Grepp
- Max Grepp
- Max Grepp Ergo (*)
- Titanium
- Silicon covers

(*) Only 36-spoke wheels

Dealer _____


Spoke guards (15)

- Transparent
- White
- Black
- Multicolour
- Black plumed red
- Black plumed blue
- Black plumed green
- Alphabet

Front castors (16)


- 80 (3") solid
- 80 (3") solid with lights
- 100 (4"x1.1/4) solid plastic rim
- 100 (4"x1) solid with lights
- 100 (4"x1) solid with alum. rim
- 100 (4"x1.1/4) solid with alum. rim
- 125 (5"x1.1/4) solid plastic rim
- 125 (5"x1) solid with alum. rim
- 125 (5"x1.1/4) solid with alum. rim
- 150 (6") pneumatic
- 150 (6") solid

Revolving anti-tip device (17)



- Right
- Left


Extended wheel plate (18)



(Not available with scissor brakes)


Back-rest pushing handles (19)

- Standard
- Without



- Height adjustable

Distance between hub centre and back-rest (20)



- 3,5 cm
 - 6 cm (standard)**
 - 8,5 cm
- (with HP=44 cm)

Anodized set (21)

(Rear wheel rims, rear wheel hubs, forks and fittings)

- Black
- Silver
- Green
- Light blue
- Red
- Orange
- Cyclamen light purple


(Light wheels only)

Lined cushion in foam polyurethane (22)




- 3 cm high
- 5 cm high

Pockets (24)



- Small under the seat



- On the back-rest

Wheel camber (25)

- 1°
- 2°
- 3°
- 4°

Brakes (23)

- Classic - short handle
- Classic - pull to brake short handle
- Classic - long handle
- Scissor

Notes

Options

- 26 "M" velcro strip (hook) on seat
- 27 "F" velcro strip (loop) on seat
- 28 Seat/back-rest connection flap
- 29 Padded rigid seat
- 30 Chest harness
- 31 Chest belt
- 32 Abductor belt
- 33 Catheter holder pocket
- 34 4 cm extended frame
- 35 Back-rest stretched bar
- 36 Back-rest anatomic board
- 37 Rx Lx Crutch holder
- 38 Pump for pneumatics
- 39 Suitcase holder
- 40 45° pelvic belt
- 41 45° padded pelvic belt
- 42 Black holdall
- 44 Umbrella holder

9,7 kg configuration

- Wheelchair 38x40x31 cm
- Nylon upholstery
- Frame colour: shiny titanium
- Automatic closure foot-rest
- Fixed back-rest pushing handles
- Aluminium sideguards with carbon fibre wings
- Scissor brakes
- Ultralight wheels
- Pneumatic high pressure smooth tyres
- Ultralight handrims
- 100 mm front castors
- Without pockets
- Max load: 100 kg