

Dealer

Key	<input type="checkbox"/> standard	<input type="radio"/> option	<h2>CHILDREN 3000</h2>
CHILDREN 3000 L Seat width <input type="checkbox"/> 24 <input type="checkbox"/> 27 <input type="checkbox"/> 30 <input type="checkbox"/> 33 <input type="checkbox"/> 36 CHILDREN 3000 PLUS L Sead width (adjustable) <input type="checkbox"/> 24 <input type="checkbox"/> 27 <input type="checkbox"/> 30 <input type="checkbox"/> 33 <input type="checkbox"/> 36 cm			
P Seat depth <input type="checkbox"/> 28 <input type="checkbox"/> 30 (adjustable from 28 to 34 cm) <input type="checkbox"/> 32 <input type="checkbox"/> 34 <input type="checkbox"/> 36 <input type="checkbox"/> 38 (adjustable from 32 to 38 cm) HS Back-rest height <input type="checkbox"/> 25 <input type="checkbox"/> 27 <input type="checkbox"/> 29 <input type="checkbox"/> 31 <input type="checkbox"/> 33 <input type="checkbox"/> 35 <input type="checkbox"/> 37 <input type="checkbox"/> 39 <input type="checkbox"/> 41 HA Front seat/floor height <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46 <input type="checkbox"/> 48 HP Rear seat/floor height <input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46 DP Foot-rest/seat distance <input type="checkbox"/> 17* <input type="checkbox"/> 19* <input type="checkbox"/> 21* <input type="checkbox"/> 23 <input type="checkbox"/> 25 <input type="checkbox"/> 27 <input type="checkbox"/> 29 <input type="checkbox"/> 31 <input type="checkbox"/> 33 <input type="checkbox"/> 35 <input type="checkbox"/> 37 <input type="checkbox"/> 39 <small>(*) Not available with detachable foot-rest</small>			
cm			
Frame colour ① <input type="checkbox"/> Metallic light blue <input type="checkbox"/> Orange <input type="checkbox"/> Pearl red <input type="checkbox"/> Metallic pink <input type="checkbox"/> Emerald green <input type="radio"/> Shiny aluminium <input type="radio"/> (see "Colour" catalogue)	Upholstery ② <input type="checkbox"/> Tension adjustable nylon back-rest <input type="radio"/> Tension adj. breathable b. <input type="radio"/> Without back-rest <input type="radio"/> Without seat	Upholstery colour ③ <small>(only nylon)</small> <input type="checkbox"/> Black <input type="checkbox"/> Light blue <input type="checkbox"/> Red <input type="checkbox"/> Green <input type="checkbox"/> Cornflower blue <input type="checkbox"/> Orange	
Front frame ④ <input type="checkbox"/> Standard <small>(Not suggested for L < 32cm)</small> <input type="checkbox"/> Abducted 	Foot-rest ⑤ <input type="checkbox"/> Detachable <input type="checkbox"/> Not detachable <input type="radio"/> Elevating (detachable)	Separate aluminium profile foot-plates ⑥ <input type="radio"/> Foot-holding belts	Single aluminium profile foot-plate ⑦ <input type="radio"/> Foot-holding belts
Side guards ⑧ <input type="checkbox"/> Plastic (standard) <input type="radio"/> Aluminium - curved <input type="radio"/> Aluminium - straight 	Height adjustable arm-rest ⑨ <small>(only with plastic side guards)</small> Elbow-rest <input type="checkbox"/> Classic (23 cm long) <input type="checkbox"/> Classic (30 cm long) <input type="checkbox"/> Straight (26 cm long)	Standard pushing handles ⑩ <input type="checkbox"/> Straight height adjustable pushing handles <input type="checkbox"/> Abducted height adjustable pushing handles 	
Rear wheels ⑪ <input type="checkbox"/> 24" Light (24 spokes) (*) <input type="checkbox"/> 22" Light (24 spokes) <input type="checkbox"/> 20" Light (20 spokes) Tyres <input type="checkbox"/> Pneum. high pressure smooth <input type="checkbox"/> Pneum. high pressure plugged <input type="radio"/> Solid <small>(*) Only with P = 34 or 36 cm and HP from 40 to 46 cm</small>	Handrims ⑫ <input type="checkbox"/> Aluminium silver <input type="checkbox"/> Chromium-plated steel <input type="radio"/> Rubber coated <input type="radio"/> Steering (*) <input type="radio"/> Tween Greep (*) <input type="radio"/> Max Grepp <input type="radio"/> Max Grepp Ergo <input type="radio"/> Silicon covers <small>(*) Not available with 20" wheels</small>	Spoke guards ⑬ <input type="radio"/> Transparent <input type="radio"/> White <input type="radio"/> Black <input type="radio"/> Multicolour <input type="radio"/> Black plumed red <input type="radio"/> Black plumed blue <input type="radio"/> Black plumed green <input type="radio"/> Alphabet	Front castors ⑭ <input type="checkbox"/> 80 (3") solid with lights <input type="checkbox"/> 80 (3") solid <input type="checkbox"/> 100 (4") solid with lights <input type="checkbox"/> 100 (4") solid <input type="radio"/> 100 (4") solid with alum. rim <input type="checkbox"/> 125 (5") solid <input type="radio"/> 125 (5") solid with alum. rim <input type="checkbox"/> 150 (6") pneumatic <input type="checkbox"/> 150 (6") solid
Rear wheel camber ⑮ <input type="checkbox"/> 1° <input type="checkbox"/> 3° <input type="checkbox"/> 5°			

