

# Hand-bikes

# COSTA

## Technical features

- Aluminium frame.
- 26" wheels assembled with pneumatic tyres.
- Height, depth, tilt adjustable seat and back-rest.
- Foot-rest adjustable in several positions.
- Rear wheel cambers: 12° or 6°.
- Brake and shift controls on handles.
- Twin front brake.
- Parking brake.
- Tension adjustable back-rest upholstery.
- Leg keeper and foot-holding belts.
- Rear detachable bumper bar.
- Wireless computer.
- Weight: 16 kg (standard set-up).



Shift handle



Disk brake



Bumper bar



Front light

### COSTA 8CR

Revo 8 gear hub shift and coaster brake.

### COSTA 8F

8 gear hub shift with lever brake.

### COSTA 27

Shimano XT 27 gear shift group with 22-32-44 chainwheel teeth cranksset, supplied with a 9 cassette sprockets:

11-12-14-16-18-21-24-28-32 or  
12-13-14-15-16-17-19-21-23 sets.



F




H

Tetraplegic handles

## Order form

COSTA 27       COSTA 8CR       COSTA 8F

- Seat width: 36  38  40  42  cm
- Seat depth: 38  40  ..... cm
- Back-rest height:  ..... cm
- Distance crankset centre/back-rest: (A) ..... cm
- Distance seat upholstery/crankset centre: (B) ..... cm
- Front seat height from the foot-rest: (C) ..... cm (from 33 to 42 cm).
- Rear seat height from the floor:  34,5 cm     38,5 cm
- Frame colour:  Ferrari red       metallic grey
- Cassette sprockets (27 gear only):  11-12-14-16-18-21-24-28-32  
 12-13-14-15-16-17-19-21-23
- Rear wheel camber:  6°     12°



Standard measurements:

- A = 55 cm
- B = 40 cm
- C = 40 cm

## Options

- Spoke guards:  with **OFFCARR** logo     without **OFFCARR** logo
- Bell:     Rear-view mirror:
- Front light:     Disk brake:
- Lined foam cushion:  3 cm     5 cm     7 cm
- Water-bottle holder:     Water-bottle:     Termic water-bottle:
- Padded rigid back-rest:
- Off-road wheels:     Flag:     Bumper bar:
- Other frame colour: ..... (see "Frame colour" leaflet)
- Tetraplegic handles:  model (F)     model (H)



Via dell' Artigianato II, 29  
35010 Villa del Conte (PD) Italy  
Tel. +39 049 9325733 - Fax +39 049 9325734  
www.offcarr.com - offcarr@offcarr.com

Dealer address

Date: \_\_\_\_\_  
Order number: \_\_\_\_\_  
Custume ref.: \_\_\_\_\_  
Notes: \_\_\_\_\_