


Dealer \_\_\_\_\_

## EMERALD



- L** Seat width  24  26  28  30  32  34  36  
**P** Seat depth  26  28  30  32  34  36  
**HS** Back-rest height  25  27  29  31  33  35  37  39  41  
**HA** Front seat/floor height  40  42  44  46  48  
**HP** Rear seat/floor height  36  38  40  42  44  46  
**DP** Foot-rest/seat distance  17\*  19\*  21\*  23  25  27  29  31  33  
 35  37  39  41

(\*) Not available with detachable foot-rest

cm

### Frame colour

- Emerald green  Orange  Metallic light blue  
 Metallic pink  Pearl red  
 Shiny aluminium  
 ..... (see "Colour" catalogue)

①

### Upholstery

- Tension adjustable nylon back-rest  
 Tension adj. breathable b.  
 Without back-rest  
 Without seat

②


### Upholstery colour (only nylon)

- Black  Light blue  
 Red  Green  
 Orange  Cornflower blue


③

### Front frame

- Standard  
 (Not suggested for L < 32cm)



- Abducted




④

### Foot-rest

- Detachable  
 Not detachable

⑤


### Separate aluminium profile foot-plates



- Foot-holding belts

⑥

### Single aluminium profile foot-plate



- Foot-holding belts

⑦

### Side guards

- Plastic (standard)



- Aluminium - curved




- Aluminium - straight



⑧

### Height adjustable arm-rest (with plastic side guards)




#### Elbow-rest


- Classic (23 cm long)  
 Classic (30 cm long)  
 Straight (26 cm long)

⑨


### Standard pushing handles



### Straight height adjustable pushing handles



### Abducted height adjustable pushing handles



⑩

### Rear wheels

- 24" Light (24 spokes) (\*)  
 22" Light (24 spokes)  
 20" Light (20 spokes)

### Tyres

- Pneum. high pressure smooth  
 Pneum. high pressure plugged  
 Solid

(\*) Only with P=36 or 34 cm and HP from 40 to 46 cm

⑪

### Handrims

- Aluminium silver  
 Chromium-plated steel  
 Rubber coated  
 Steering (\*)  
 Tween Greep (\*)  
 Max Grepp  
 Max Grepp Ergo  
 Silicon covers

⑫

### Spoke guards

- Transparent  
 White  
 Black  
 Multicolour  
 Black plumed red  
 Black plumed blue  
 Black plumed green  
 Alphabet

⑬

### Front castors

- 80 (3") solid with lights  
 80 (3") solid  
 100 (4") solid with lights  
 100 (4") solid  
 100 (4") solid with alum. rim  
 125 (5") solid  
 125 (5") solid with alum. rim  
 150 (6") pneumatic  
 150 (6") solid

⑭

### Rear wheel camber

- 1°  3°  5°

⑮

(\*) Not available with 20" wheels

