

Dealer


Key

standard

option

QUASAR KID ADJUSTABLE AXLE

L Seat width	<input type="checkbox"/> 24 <input type="checkbox"/> 26 <input type="checkbox"/> 28 <input type="checkbox"/> 30 <input type="checkbox"/> 32 <input type="checkbox"/> 34 <input type="checkbox"/> 36
P Seat depth	<input type="checkbox"/> 26 <input type="checkbox"/> 28 <input type="checkbox"/> 30 <input type="checkbox"/> 32 <input type="checkbox"/> 34 <input type="checkbox"/> 36
HS Back-rest height	<input type="checkbox"/> 25 <input type="checkbox"/> 27 <input type="checkbox"/> 29 <input type="checkbox"/> 31 <input type="checkbox"/> 33 <input type="checkbox"/> 35 <input type="checkbox"/> 37 <input type="checkbox"/> 39 <input type="checkbox"/> 41
HA Front seat/floor height	<input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46 <input type="checkbox"/> 48
HP Rear seat/floor height	<input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46
DP Foot-rest/seat distance cm (min. 15 cm, max 35 cm for HA=40 cm) (min. 15 cm, max 41 cm for HA=46 cm)



Frame colour

- Shiny titanium Pearl white Metallic light blue
 Sage green Water green Antique gold
 Avio blue Potassium light purple Pearl grey
(see "Colours" catalogue)

①

Upholstery

- Tension adjustable nylon back-rest
 Tension adj. breathable bark-rest
 Without back-rest
 Without seat

②

Upholstery colour


(nylon only)

- Black Light blue
 Red Green
 Cornflower blue Orange

③

Front frame


④



- Standard (not suggested with L < 32)
 Abducted
 Y (only with L= 34-36) (single foot-rest only)

- Fixed handles Height adjustable handles
 Abducted height adjustable handles Height adjustable pushing bar

⑤



Separate aluminium foot-rest

⑥


Only for L>32 cm (standard front frame)
Only for L>28 cm (abducted front frame)



- Foot-holding belts

Single aluminium foot-rest

⑦




- Foot-holding belts

Single closed carbon foot-rest

⑧

(only with DP > 37)



- Open

Side guard

⑨

- Aluminium with carbon fibre wing



- Straight aluminium Plastic



Height adjustable arm-rest

⑩

(only with plastic clothes guard)




Elbow-rest

- Classical (23 cm long)
 Classical (30 cm long)
 Straight (26 cm long)

Distance between the hub centre and back-rest


⑪



- 3 cm 7 cm
 5 cm (standard) 9 cm

Back-rest tilt


⑫



- 89°
 91° (standard)
 93°

Revolving anti-tip device

⑬



- Left
 Right
 Left and Right

