

Dealer \_\_\_\_\_


### Key

standard

option

# RUBY

<b>L</b> Seat width	<input type="checkbox"/> 24 <input type="checkbox"/> 26 <input type="checkbox"/> 28 <input type="checkbox"/> 30 <input type="checkbox"/> 32 <input type="checkbox"/> 34 <input type="checkbox"/> 36
<b>P</b> Seat depth	<input type="checkbox"/> 26 <input type="checkbox"/> 28 <input type="checkbox"/> 30 <input type="checkbox"/> 32 <input type="checkbox"/> 34 <input type="checkbox"/> 36
<b>HS</b> Back-rest height	<input type="checkbox"/> 25 <input type="checkbox"/> 27 <input type="checkbox"/> 29 <input type="checkbox"/> 31 <input type="checkbox"/> 33 <input type="checkbox"/> 35 <input type="checkbox"/> 37 <input type="checkbox"/> 39 <input type="checkbox"/> 41
<b>HA</b> Front seat/floor height	<input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46 <input type="checkbox"/> 48
<b>HP</b> Rear seat/floor height	<input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46
<b>DP</b> Foot-rest/seat distance	..... cm (min. 15 cm, max 35 cm for HA=40 cm) (min. 15 cm, max 41 cm for HA=46 cm)



### Frame colour

- Pearl red       Pearl white       Metallic pink  
 Metallic light blue       Orange  
 Shiny aluminium  
 ..... (see "Colours" catalogue)

①

### Upholstery

- Tension adjustable nylon back-rest  
 Tension adj. breathable b.  
 Without back-rest  
 Without seat

②

### Upholstery colour


(nylon only)

- Black       Light blue  
 Red       Green  
 Cornflower blue       Orange

③

### Front frame


④



- Standard       Abducted  
 (not suggested with L < 32)

### Single closed carbon foot-rest

(only with DP > 37 cm)




- Open

⑤

### Separate aluminium foot-rest


Only for L>32 cm (standard front frame)  
Only for L>28 cm (abducted front frame)



- Foot-holding belts

⑥

### Single aluminium foot-rest



- Foot-holding belts

⑦


### Standard pushing handles

### Straight height adjustable pushing handles

### Abducted height adjustable pushing handles

### Height adjustable pushing bar

⑧



### Side guards

- Plastic  
 Straight aluminium



⑨

### Height adjustable arm-rest

(only with plastic clothes guard)




### Elbow-rest

- Classic (23 cm long)  
 Classic (30 cm long)  
 Straight (26 cm long)

⑩

### Distance between the hub centre and back-rest


⑪



- 3 cm  
 5 cm (standard)  
 7 cm  
 9 cm

### Back-rest tilt


⑫



- 89°  
 91° (standard)  
 93°

### Revolving anti-tip device

⑬



- Right  
 Left  
 Right and Left

### Notes

