

Dealer \_\_\_\_\_

### Key


standard

option

- L** Seat width 34 36 38 40 42 44 46  
**P** Seat depth 36 38 40 42 44  
**HS** Back-rest height 27 29 31 33 35 37 39  
**HA** Front seat/floor height 46 48 50  
**HP** Rear seat/floor height 39 41 43 45 47  
**DP** Foot-rest/seat distance 31 33 35 37 39 41 43 45 47 49

cm

## VENUS ELITE



### Frame colour

- Shiny titanium  Pearl white  Glossy black  
 Sage green  Antique gold  
 Avio blue  Potassium light purple  Pearl grey  
 ..... (see "Colours" catalogue)

①


### Upholstery

②

- Tension adjustable nylon back-rest  
 Tension adj. breathable back-rest  
 Without back-rest

### Single tubular titanium foot plate

③




Open

Closed in carbon fibre

### Pushing handles

④




Without pushing handles

Standard pushing handles

Height adjustable pushing handles

### Distance between the hub centre and back-rest


⑤



- 4 cm  9 cm  
 5 cm  11 cm  
 7 cm (standard)

### Back-rest tilt


⑥



- 89°  
 91° (standard)  
 93°

### Frame length

⑦




- 42,5 (std)  
 46,5  
 .....

cm

### Seat with shock absorbers

⑧



User's weight : ..... kg

It is not possible to have two anti-tip devices

### Side guards

⑨




Aluminium with carbon fibre wing

Straight aluminium

### "L" detachable arm-rest - h 22 cm

⑩



### Elbow rest

- Tubular  
 Straight (26 cm long)  
 Classic (23 cm long)  
 Classic (30 cm long)

### Brakes

⑪

- Push brake  
 Pull brake




Scissor brake



### Revolving anti-tip device

⑫



Right

Left

### Notes

