

ALHENA... LA LEGGEREZZA NEL SUO DNA

*THE LIGHTNESS
ON THE DNA*

Telaio pieghevole
in titanio

Titanium folding
frame



 **CARRER®**

ALHENA

Leggerezza e tecnologia nel DNA di un prodotto unico e raffinato.
Telaio pieghevole in titanio.

State-of-the-art and lightness in the DNA of a unique and refined product.
Titanium folding frame.

Una nuova carrozzina superleggera è entrata a far parte della costellazione **OFFCARR**.

Con soli 9,7 kg **ALHENA**, tra le più leggere della sua categoria, offre agilità, maneggevolezza e versatilità.

L'elegante e dinamica piastra HwP (Helicoid wheel plate) propone un'ampia scelta di posizioni per la ruota posteriore e consente di personalizzare al meglio l'assetto della carrozzina.

A new ultra lightweight wheelchair has joined **OFFCARR**'s constellation.

ALHENA, within 9,7 kg only, offers agility, handling and versatility becoming one of the lightest of its range.

The smart and dynamic rear wheel adjustment, HwP (Helicoid wheel Plate), provides a whole range of rear wheel positions to personalise the wheelchair's set-up.

Fornita di serie con:

- salvaabiti in fibra di carbonio
- ruote e corrimano "Ultralight"
- pedana a chiusura automatica
- doppia crociera.

The basic wheelchair version includes:

- carbon fibre clothes guards
- "Ultralight" wheels and pushrims
- automatic closure foot-rest
- double cross bar



Piastra HwP
HwP plate



Pedana a chiusura automatica
Automatic closure foot-rest



Minimo ingombro
Takes up less space



34	52 57	34	29	A: min. 38 max. 46	media average 87	29	9,7 Kg con seduta with seat 38 x 40 cm
36	54 59	36	31	B: min. 43 max. 50			
38	56 61	38	33	misure indicative approximate measurements			
40	58 63	40	35				
42	60 65	42	37				
44	62 67	44	39				
46	64 69	46	41				



OFF CARR s.r.l.
Via dell'Artigianato II, 29
35010 Villa del Conte
Padova Italy

Tel. +39 049 9325733
Fax +39 049 9325734
offcarr@offcar.com
commerciale@offcarr.com

www.offcarr.com

AZIENDA CON SISTEMA QUALITÀ CERTIFICATO
=UNI EN ISO 9001:2000=

AZIENDA CON SISTEMA DI GESTIONE
PER LA QUALITÀ CERTIFICATO
=UNI EN ISO 13485:2004=

