


OFFCARR Everyday Ultralight Adult wheelchairs — Comparison Chart

 more details can be found in specific catalogues/order forms																					
																					
Model name	Heka	Eos	Eos3	Quasar	Quasar GT	Qpx	Fixed	Venus Adventure	Venus / Venus Elite	Themis Prestige	Themis	Idra 2.0 / Idra 2.E	Fenice	Funky	Zodiac	Alhena	Halley Ti	Halley	Diva	Vega / Vega.E	Althea
Category	Ultralight	Ultralight	Ultralight	Ultralight	Ultralight	Ultralight	Ultralight	Ultralight	Ultralight	Ultralight	Ultralight	Ultralight	Ultralight	Ultralight	Ultralight	Ultralight	Ultralight	Ultralight	Ultralight	Ultralight	Light
Frame	Rigid	Rigid	Hinged	Rigid	Rigid	Hybrid	Rigid / Hinged	Rigid	Rigid	Rigid	Rigid	Rigid	Rigid	Rigid	Rigid	Folding	Folding	Folding	Folding	Folding	Folding
Centre of Gravity	Welded	Welded / Adjustable	Welded / Adjustable	Welded / Adjustable	Welded / Adjustable	Adjustable	Welded	Adjustable	Adjustable	Welded	Adjustable	Welded / Adjustable	Adjustable	Adjustable	Welded	Adjustable	Adjustable	Adjustable	Adjustable	Adjustable	Adjustable
Setup (front and rear height)	Welded	Welded	Welded	Welded	Welded	Welded	Welded	Adjustable	Adjustable	Welded	Adjustable	Welded	Adjustable	Adjustable	Welded	Adjustable	Adjustable	Adjustable	Adjustable	Adjustable	Adjustable
Frame material	Titanium	Titanium	Titanium	Titanium	Titanium	Aluminium / Titanium	Titanium	Titanium + Aluminium	Titanium + Aluminium / Titanium	Aluminium	Aluminium	Aluminium	Aluminium	Aluminium	Aluminium	Titanium	Titanium	Aluminium	Aluminium	Aluminium	Aluminium
Removable Footrests	No	No	No	No	No	No	No	No	No	No	No	Idra 2.E	No	No	No	Possible	Possible	Possible	Possible	Possible / Vega.E standard	Possible
Maximum user weight	120 kg	120 kg	120 kg	120 kg	120 kg	120 kg	120 kg	120 kg	120 kg	120 kg	120 kg	120 kg	120 kg	120 kg	120 kg	150 kg	120 kg	120 kg	120 kg	120 kg	120 kg
Crash tested ISO 7176-19	No	No	No	No	Yes 120 kg	No	No	No	No	No	Yes 120 kg	No	No	No	No	Yes 150 kg	Yes 120 kg	No	No	Yes 120 kg	No
Complete weight (minimum) ¹	7.6	7.5 / 7.8	8.5 / 8.8	7 / 7.4	7.6 / 8	10.6 / 10.3	9.8 / 10.1	11.3	10.4 / 8.9	9	10.5	2.0 9.1 9.6 2.E 10.5 10.9	10.5	9.8	10.5	9.7	9.9	12.4	11.5	12.7	15.8
Transport weight (minimum) ²	4.6	4.4 / 4.7	5.4 / 5.7	4 / 4.4	4.8 / 5.2		6.5 / 6.8	8.1	7.3 / 5.8	6.1	7.5	2.0 6.1 2.E 5.9	7.6	7.0	7.6						
Seat Width (L) ³	34 – 48	34 – 48	34 – 44	34 – 48	34 – 44	34 – 44	34 – 44	34 – 48	34 – 48	34 – 48	34 – 50	34 – 44	34 – 48	34 – 50	34 – 54	34 – 48	34 – 48	34 – 48	34 – 50	34 – 48	34 – 60
Seat Depth (P) ³	34 – 50	36 – 50	36 – 46	36 – 50	36 – 46	36 – 46	36 – 46	36 – 48	36 – 48	34 – 48	34 – 48	36 – 46	36 – 48	36 – 48	36 – 50	34 – 50	36 – 50	36 – 48	34 – 50	34 – 50	36 – 50
Front Height (HA) ³	40 – 52	46 – 50	46 – 50	46 – 50	46 – 50	46 – 50	46 – 52	46 – 50	46 – 50	40 – 52	40 – 52	46 – 50	46 – 50	46 – 50	46 – 50	42 – 50	46 – 50	46 – 50	46 – 50	42 – 54	42 – 54
Rear Height (HP) ³	37 – 50	40 – 48	40 – 48	40 – 48	40 – 48	40 – 48	40 – 48	43 – 47	39 – 47	37 – 50	37 – 50	40 – 48	40 – 46	40 – 46	40 – 46	38 – 46	40 – 46	40 – 46	38 – 50	38 – 50	38 – 50
Frame shapes	Ergo / Y / V / U	Ergo / Y / V / U	Ergo / U	Ergo / Y / V / U / Abducted	Ergo / Y / V / U	Ergo / Y / U		Y	Y	V / U	V / U	Ergo / 2.0 Y / V / U 2.E U	Y / V / U / Abducted	Y / U	Ergo / Y / V / U	Y / U	Y / U	U	U / Abducted	U	U / Abducted
Tube type	Ergonomic 35×30 mm	Round 31.7 mm	Round 31.7 mm	Round 25.4 mm	Round 25.4 mm	Round 25 / 25.4 mm	Oval 34×30 mm	Round 31.7 mm	Round 31.7 mm	Ergonomic 35×32 mm	Ergonomic 35×32 mm	Round 25 mm	Oval 32×28 mm	Round 30 mm	Oval 32×28 mm	Round 25.4 mm	Round 25.4 mm	Round 25 mm	Round 25 mm	Round 25 mm	Round 25 mm
Frame colour	31 available colours	31 available colours	Shiny or natural titanium	31 available colours	31 available colours	29 / 31 colours	31 available colours	30 available colours	29 / 31 colours	29 available colours	29 available colours	29 available colours	29 available colours	Matte black	29 available colours	31 available colours	31 available colours	29 available colours	29 available colours	29 available colours	29 available colours
¹ Weight, usually 38×40 cm, with ultralight wheels, pushrims, brakes and sideguards, more details in catalogues. ² Weight with welded backrest, removing everything possible (e.g. wheels, sideguards, footplates). ³ A wider range of measurements is possible with specific limitations.																					
																			