

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____ User's weight Kg _____

PAG. 1/2
Oct. 2026


Dealer _____

Key ☐ standard ☐ option

L Seat width ☐34 ☐36 ☐38 ☐40 ☐42 ☐44 ☐46 ☐48
P Seat depth ☐34 ☐36 ☐38 ☒40 ☐42 ☐44 ☐46 ☐48 ☐50
HS Back-rest height ☐29 ☐31 ☐33 ☐35 ☐37 ☐39 ☐41
HA Front seat/floor height ☐42 ☐43 ☐44 ☐46 ☐48 ☒50
HP Rear seat/floor height ☐38 ☐40 ☐42 ☒44 ☐46
DP Foot-rest/seat distance ☐33 ☐35 ☐37 ☐39 ☐41 ☐43 ☐45 ☐47 cm

ALHENA

**Titanium folding frame
 crash-tested**




Frame colour
☐ Shiny titanium ☐ Natural titanium
☐ Glossy anthracite ☐ Matte black
☐ Basic colours
☐ Special colours
 (see "Colours" catalogue)


Backrest upholstery
☐ Tension adjustable nylon back-rest
☐ Tens.adjustable breathable backrest
☐ Without any backrest upholstery
☒ **FAMEKH** custom postural backrest
 (fill up the proper order form)

Upholstery colour
☐ Black
☐ Coloured.....
 Available colours:
 NYLON: orange, red, green, cornflower blue, light blue;
 BREATHABLE: red, grey.

☐ **Not detachable foot-rest** ④
 (with DP=43 cm and P=40 cm)
 (standard)




☐ **Not detachable foot-rest**
 (with DP=43 cm and P=40 cm)



☐ **Detachable foot-rest**


☐ 84° (std)
☐ 76°




(approx.measurements)

Front Frame ⑤

☐ STANDARD





☐ Y
 (only in case of
 not detachable foot-rest)



**Separate aluminium profile
 tilt adjustable footplates** ⑧

☐ forward


☐ backward



☐ Foot-holding belts

☐ **Closed carbon single
 footplate** ⑥


(Note: it might interfere with DP)
 (max W 48cm)



☐ Foot-holding belts

☐ **AUTOMATIC
 closure footplate
 in carbon** ⑦
 (up to W46)

for W 48cm
 (Note: it might interfere with DP)
 (max W 48cm)




☐ Foot-holding belts

**Single aluminium profile
 tilt adjustable footplate** ⑨


(Note: it might interfere with DP)

☐ forward

☐ backward




☐ up & outward



☐ Foot-holding belts

Titanium tube structure ⑩

(Note: it might interfere with DP)



☐ open
☐ closed by carbon plate

☐ Forward

☐ Backward
 (default)

Side-guards ⑩

☐ Aluminium with carbon
 fibre wing



☐ Plastic




☐ Curved aluminium
 (same colour or black)



☐ Straight aluminium
 (same colour or black)



☐ **Tip-up,
 detachable sport
 armrest h. 22 cm** ⑪
 (open with sideguards)




Elbow-rest

☐ Tubular
☐ Classic (23 cm long)

☐ **Tip-up height
 adjustable from 22
 to 28 cm "L" armrest** ⑫
 (detachable, open with sideguards)

☐ Not detachable



Elbow-rest

☐ Tubular
☐ Straight (26 cm long)
☐ Classic (23 cm long)
☐ Classic (30 cm long)

Rear wheel ⑬

☐ 24" Ultralight (24 spokes) black
☐ 24" Light (24 black spokes)
☐ 24" Star 6x3 (18 sps) silver or black
☐ 24" (36 spokes)
☐ 24" Mountain-bike CITYBIKE tyre
☐ 24" Mountain-bike OFFROAD tyre
☐ 24" Spineray LX
 Spokes colour.....(black def.)
☐ 25" Spineray LX
 Spokes colour.....(black def.)
☐ 25" Light (24 spokes) black

Tyre

☐ Pneumatic Marathon Plus 24"/25"
☐ Pneum. high press.plugged 24"
☐ Pneum. high press. smooth
☐ Solid black KIK 24"
 (not available with Ultralight)

Pushrim ⑭

☐ Aluminium silver
☐ Aluminium dark
☐ OVAAL alum. silver 24"
☐ OVAAL alum. dark 24"
☐ OVAALG(dark w/silicon ins.)24"
☐ Shiny steel 24"
☐ Rubber coated
☐ Titanium
☐ Carbolife CURVE L24"
☐ CL CUREL TETRAGRIP 24"
☐ CL QUADRO Rubber ctd 24"
☐ Carbolife GEKKO 24"
☐ Silicon covers
 OVAAL/S/G only in size 24"
 Only OVAAL/GEKKO available for Star6x3
 OVAAL & Carbolife not available for MTB

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____

PAG. 2/2
 Oct.2026

ALHENA

Dealer _____

Spokeguards

15

- ☐ Transparent 24"
- ☐ Black 24"
- ☐ Fantasy
- code pattern
- (see catalogue of spokeguards)

Castor wheel

16

- ☐ 80 (3") solid
- ☐ 80 (3") solid with lights
- ☐ 100 (4"x1.1/4) black solid plastic rim
- ☐ 100 (4"x1) solid with lights
- ☐ 100 (4"x1.1/4) sol.w/black alum.rim
- ☐ 100 (4"x1.1/4) sol.w/silver alum.rim
- ☐ 125 (5"x1.1/4) black solid plastic rim
- ☐ 125 (5"x1.1/4) sol.w/black alum.rim
- ☐ 125 (5"x1.1/4) sol.w/silver alum.rim
- ☐ 150 (6") pneumatic
- ☐ 150 (6") black solid

Revolving anti-tip device

17

- ☐ Right
- ☐ Left




Unbalancing Pedal

- ☐ Right
- ☐ Left



Extended wheel plate

18



(Not available with scissor brakes)


Push handles

19

- ☐ Standard
- ☐ Without
- ☐ Height adjustable
- ☐ Foldable push handles

Distance between hub centre and backrest

20



- ☐ 3,5 cm
- ☐ 6 cm (standard)
- ☐ 8,5 cm
- (with HP=44 cm)

Anodized set


(Rear wheel rims, rear wheel hubs, forks and fittings)

21

- ☐ Black
- ☐ Silver
- ☐ Light blue
- ☐ Red
- ☐ Orange
- ☐ Cyclamen light purple
- (Light 24" wheels only)

Lined cushion in foam polyurethane


22




- ☐ 3 cm high
- ☐ 5 cm high

Pockets

24



☐ Small under the seat



☐ On the back-rest

Camber

25

- ☐ 1° ☐ 3°
- ☐ 2° ☐ 4°
- ☐ Quick rel.axles Ø 12 mm with 42mm cambered bushes (default)
- ☐ Quick rel. axles Ø 1/2" (12,7 mm) with 42mm cambered bushes
- ☐ Quick rel. axles Ø 1/2" (12,7 mm) with 50mm cambered bushes

Brakes

23

- ☐ Classic - short handle
- ☐ Classic - pulling brake short handle
- ☐ Classic - long handle
- ☐ Scissor
- ☐ Composite

Options

- 26 ☐ "M" velcro strip (hook) on seat
- 27 ☐ "M" velcro for VICAIR strip (hook) on seat
- 28 ☐ "F" velcro strip (loop) on seat
- 29 ☐ Seat/backrest connection flap
- 30 ☐ Chest harness
- 31 ☐ Chest belt
- 32 ☐ Abductor belt
- 33 ☐ Catheter holder pocket
- 34 ☐ Backrest stretched bar
- 35 ☐ Rx ☐ Lx Crutch hold./single
- 36 ☐ Rx ☐ Lx Crutch hold./double
- 37 ☐ Pump for pneumatics
- 38 ☐ Suitcase holder
- 39 ☐ 45° pelvic belt
- 40 ☐ 45° padded pelvic belt
- 41 ☐ Front protective covers
- 42 ☐ Reinforced frame (max load 100-120Kg)
- 43 ☐ Stickers for car fastenings

crash tested

9,7 kg configuration

- Wheelchair 38x40x31 cm
- Nylon upholstery
- Titanium frame col.: shiny titanium
- Automatic closure footrest
- Fixed backrest pushing handles
- Aluminium sideguards with carbon fibre wings
- Scissor brakes
- Ultralight wheels
- Pneumatic high pressure smooth tyres
- Ultralight handrims
- 100 mm front castors
- Without pockets
- Max load: 100 kg

Notes