

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____ User's weight Kg _____

Key ☐ standard
☐ option


PAGE 1/2
 Oct. 2026

Dealer _____

- L** Seat width ☐34 ☐36 ☐38 ☐40 ☐42 ☐44 ☐46 ☐48 ☐50
☐52* ☐56* ☐60 (including eventual stretched bar, tension adj. upholstery & cover for footplate according to setup)
- P** Seat depth ☐36 ☐38 ☐40 ☐42 ☐44 ☐46
 Adjustable up to 6 cm (max.50cm)
☐48 ☐50
- HS** Back-rest height ☐29 ☐31 ☐33 ☐35 ☐37 ☐39 ☐41 ☐43 ☐45
☐47 ☐49 ☐51 ☐55 ☐60 (included stretched bar)
- HA** Front seat/floor height ☐38(fc80max) ☐39(fc100max) ☐40(fc125max) ☐42(fc150max) ☐45(fc175max)
 (with MAX size of front castors)
 (NOTE: 175 front castors do not allow 82° detachable footrests) ☐44(fc150max) ☐46(fc175max)
- HP** Rear seat/floor height ☐34 ☐36 ☐38 ☐40 ☐42 ☐44 ☐45 ☐46 ☐48
 NOTE: 24" rear wheels (if HP min.38/max 48) 22" (if HP min.36/max 45) 20" (if HP min.34/max 43)
- DP** Foot-rest/seat distance ☐34 ☐35 ☐37 ☐39 ☐41 ☐43 ☐45
 (from 17 to 32 cm footrest directly assembled on the frame) cm

ALTHEA.19 LOWERED SET-UP

alum.foldable frame
 crash tested



Frame colour

- ☐ Metallic blue ☐ Pearl red ☐ Matte black
- ☐ Basic colour..... from:
- ☐ Special colour..... from:
- (see "Colours" catalogue)

Upholstery

- ☐ Tension adj. Nylon. backrest up.
☐ Tension adj. Breath. backrest up.
☐ Without any backrest upholstery
☐ With divided seat (joint by velcro, def.)
☐ With one entire seat (not depth adjustable)

Upholstery colour

★ (Coloured uphol. only Nylon fabric)

- ☐ Black
☐ Red★ ☐ Orange★
☐ Light blue★ ☐ Green★
☐ Cornflower blue★

☐ Not Detachable foot-rest


☐ Detachable foot-rest 82° angle

(DP 41 cm)
 (STANDARD)


☐ Detachable foot-rest 68° angle

(DP 41 cm)

☐ Elevating foot-rest (detachable)




☐ Plastic separate



- ☐ Cover for footplates*
☐ Foot-holding belts (up to W50)
☐ Heel stops (up to W50)
 * Eventual cover f/footplates included for W52-60

Single aluminium profile tilt adjustable foot-plate:


(not available w/detachable or elevating footrests)
 (Note: it might interfere with DP)



- ☐ Foot-holding belts
☐ Heel stops

Separate aluminium profile tilt adjustable foot-plates:


- ☐ forward ☐ backward



- ☐ Foot-holding belts (up to W50)
 * Eventual cover f/footplates included for W52-60

☐ AUTOMATIC closure foot-plate in carbon (up to W46)

(not available w/detachable or elevating footrests)
 (Note: it might interfere with DP)
 (max W 48cm)



- ☐ Foot-holding belts

Side guards

- ☐ Plastic
- ☐ Aluminium
 (same colour)

Detachable DESK armrest - tip-up

- ☐ OPEN H22
 with side-guard
 (only if HP ≥44 cm)
- ☐ OPEN H27
 with side-guard
 (if HP≤42cm)

- ☐ CLOSED H22
 (from HP40/R24"x1
 only if without
 any table fastenings)
- ☐ CLOSED H27 cm
 from seat
 (compulsory for
 HP=42cm&R24")

Elbow-rest

- ☐ Classic (23 cm long)
☐ Straight (26 cm long)
☐ Tubular (only if open armrest)
☐ Hemiplegic RX
☐ Hemiplegic LX

note: see also Hemi modular armrest ord.form

Detachable SPORT armrest - tip-up

- ☐ OPEN H22
 with side-guard
 (only if HP ≥44 cm)
- ☐ OPEN H27
 with side-guard
 (if HP≤42cm)

- ☐ CLOSED H22
 (from HP40/R24"x1
 only if without
 any table fastenings)
- ☐ CLOSED H27 cm
 from seat
 (compulsory for
 HP=42cm&R24")

Elbow-rest

- ☐ Classic (23 cm long)
☐ Straight (26 cm long)
☐ Tubular (only if open armrest)
☐ Hemiplegic RX
☐ Hemiplegic LX

note: see also Hemi mod.armrest ord.form

Detachable "U" armrest - tip-up

- ☐ OPEN H22
 with side-guard
 (only if HP ≥44 cm)
- ☐ OPEN H27
 with side-guard
 (if HP≤42cm)

- ☐ CLOSED H22
 (from HP40/R24"x1
 only if without
 any table fastenings)
- ☐ CLOSED H27 cm
 from seat
 (compulsory for
 HP=42cm&R24")

Elbow-rest

- ☐ Straight (26 cm long)
☐ Straight (36 cm long)
☐ Tubular (only if open armrest)
☐ Hemiplegic RX
☐ Hemiplegic LX

note: see also Hemi mod.armrest ord.form

Height adjustable tip-up armrest (only closed)

- ☐ Desk model
- ☐ Sport model
- ☐ "U" model

Elbow-rest

- ☐ Classic (30 cm long)
☐ Classic (36 cm long)
☐ Straight (26 cm long)
☐ Straight (36 cm long)
☐ Hemiplegic RX
☐ Hemiplegic LX

note: see also Hemi modular armrest ord.form

☐ Order N° _____ ☐ Quotation
Date _____ Signature _____
Customer Ref. _____

PAGE 2/2
Oct.2026

ALTHEA19 LOWERED SET-UP

Dealer _____

Rear wheels ⑭

- ☐ 24" 36 spokes ☐ 12" solid*
☐ 22" 36 spokes ☐ 12" pneumatic*
☐ 20" 36 spokes (N.A. if HP=42 if HP44: only curved rear frame)
☐ 24" MTB citybike tyre
☐ 24" MTB offroad tyre
☐ 24" Light (24 spokes)

- ☐ Pneum. Plugged (x 1"3/8) **Tyres**
Pneum. high press. ☐ Plugged ☐ Smooth
☐ Solid grey 24x1.3/8"/22"x1or1.3/8"
☐ Solid black KIK 24"
☐ 22"/24"Pneum.Marathon Plus

Pushrims ⑮

- ☐ Aluminium silver
☐ Shiny steel 24"
☐ Aluminium dark
☐ Rubber coated
☐ Carbolife CURVE L 24"(*)
☐ OVAAL silver alum. (*)
☐ OVAALS dark alum. (*)
☐ OVAALG dark w/sil. insert (*)
☐ Carbolife GEKKO 24" (*)
☐ Silicon covers 24"
(*) not available with MTB wheels

Spoke guards ⑯


- ☐ Transparent
☐ Black
☐ Fantasy
code pattern:
(see catalogue of spokeguards)

Front castors ⑰


- ☐ 100 (4"x1.1/4) solid - plastic rim
☐ 100 (4"x1.1/4) sol.silver alum. rim
☐ 100 (4"x1.1/4) sol.black alum. rim
☐ 125 (5"x1.1/4) solid - plastic rim
☐ 125 (5"x1.1/4) sol.silver alum. rim
☐ 125 (5"x1.1/4) sol.black alum. rim
☐ 150 (6") pneumatic
☐ 150 (6") solid
☐ 175 (7") solid
☐ 175 (7") pneumatic

Rear frame ⑱


- ☐ STRAIGHT rear frame with small wheels for narrow passages (default)




- ☐ STRAIGHT rear frame WITHOUT any small wheels f/narrow passages



- ☐ STRAIGHT rear frame with LEVER-activated narrow passage wheels



- ☐ CURVED rear frame



(Only if rear seat/floor height HP is between 44 and 48 cm)

(It is not possible to assemble the small wheels for narrow passages)

Max
user weight
100 kg

Unbalancing system ⑲

- ☐ Right ☐ Right
☐ Left ☐ Left
(Only with straight rear frame) (Only with curved rear frame)

Anti-tip device (Only with straight rear frame)

- ☐ Right
☐ Left


Revolving anti-tip device e (Only with curved rear frame)

- ☐ Right
☐ Left


Single drive ⑳

(only w/24" or 22" single drive rear w.)

- ☐ Right
☐ Left




Extended wheel plate ㉑



Push handles ㉒

- ☐ Standard ☐ Height adjustable(*) ☐ Hinged backrest (no stretched bar)




Distance between the hub centre and back-rest ㉓

- ☐ -2 cm (ext. wheel plate)
☐ 3,5 cm
☐ 5,5 cm
☐ 7,5 cm

(with HP=44 cm)

Adjustable swing away Lateral supports ㉔




- ☐ Right
☐ Left

Specify the height "h" from the center of the swing away to the seat
..... cm

Head-rest ㉕

Backrest extension



- ☐ Fixed
☐ Height adjustable
Default: h = 20 cm d = 15 c m

Orbital

- ☐ Anatomic



- ☐ Moulded



- ☐ With lateral supports




Table ㉖

- ☐ Transparent
☐ Rx OLx Hemiplegic transp.*
☐ Plastic grey
☐ Soft padded

Specify the fastening:

- ☐ single ☐ double



(Only with arm-rest)
(*with reinforced armrest) OXL if W>50cm

Lined cushion in foam polyurethane ㉗



- ☐ 3 cm high
☐ 5 cm high


Brakes ㉘

- ☐ Classic - short handle
☐ Classic - long handle
☐ Drum*
☐ Right single brake
☐ Left single brake
☐ Drum*+classic
(* only with drum rear wheels)

Options

- 29 ☐ "M" velcro strip (hook) on seat
30 ☐ "M" velcro strip for Vicair
31 ☐ "F" velcro strip (loop) on seat
32 ☐ Seat/back-rest connection flap
33 ☐ Chest harness
34 ☐ Chest belt
35 ☐ Abductor belt
36 ☐ Front protective covers
37 ☐ Stretched bar
38 ☐ Catheter holder pocket (note: only with entire seat)
39 ☐ R O L Crutch hold./single
40 ☐ R O L Crutch hold./double

Abductor ㉙



- ☐ 1 ☐ 3 ☐ 5 ☐ 7 cm

- ☐ cm

Choose or specify the height from seat cover

- 41 ☐ Pump for pneumatics
42 ☐ Suitcase holder
43 ☐ 45° pelvic belt
44 ☐ 45° padded pelvic belt
45 ☐ Black backpack
46 ☐ IV bar + poles
47 ☐ IV poles
48 ☐ Padded rigid cushion (L52max)
49 ☐ Padded rigid cushion (L56-60)
50 ☐ Q.rel.axles Ø12mm w/42mm bushes (def.)
51 ☐ Q.rel.a. Ø1/2" (12,7mm) w/42mm bushes
52 ☐ Q.rel.a. Ø1/2" (12,7mm) w/50mm bushes
53 ☐ Stickers for car fastenings (escl.: load>120Kg, R12", UltraL.RW, hinged bckr)