

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____ User's weight Kg _____

Key ☐ standard
☐ option

PAGE 1/2
 May 2025

Dealer _____

L ☐ RANGE 1 W ☐ 28 ☐ 30 ☐ 32 ☐ 34 ☐ 36 ☐ 38
☐ RANGE 2 W P > 36cm ☐ 32 ☐ 34 ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44
☐ RANGE 3 W P > 36cm ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48 ☐ 50 ☐ 52
(compulsory expandable stretched bar to be added)

P Seat depth ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46
 Adjustable up to 6cm (max.50cm) ☐ 48 ☐ 50

HS Back-rest height ☐ 29 ☐ 31 ☐ 33 ☐ 35 ☐ 37 ☐ 39 ☐ 41 ☐ 45
☐ 47 ☐ 49 ☐ 51 ☐ 55 ☐ 60 (included stretched bar that can be enlarged)

HA Front seat/floor height ☐ 38(fc80max) ☐ 39(fc100max) ☐ 40(fc125max) ☐ 42(fc150max) ☐ 45(fc175max)
 with MAX size of front castors (fc)
 (NOTE: 175 front castors do not allow 82° detachable footrests) ☐ 44(fc150max) ☐ 46(fc175max)


HP Rear seat/floor height ☐ 34 ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 45 ☐ 46 ☐ 48
 NOTE: 24" rear wheels (if HP min.38/max 48) 22" (if HP min.36/max 45) 20" (if HP min.34/max 43)

DP Foot-rest/seat distance ☐ 34 ☐ 35 ☐ 37 ☐ 39 ☐ 41 ☐ 43 ☐ 45
 (from 17 to 32 cm footrest directly assembled on the frame) cm

ALTHEA19PLUS

LOWERED SET-UP

Width & depth adjust.
foldable alum.frame



Frame colour ☐ Matte black ①
☐ Metallic blue ☐ Pearl red ☐ Metallic grey
☐ Pearl white ☐ Glossy black
☐ from: Neon colour
(see "Colours" catalogue)

Backrest upholstery ②
☐ Tension adj. nylon. back
☐ Tension adj. breath. back
☐ Without any backrest upholstery
☐ With divided seat (joint by velcro, def.)
 Seat upholstery can grow +6cm max

Upholstery colour ③
 ★ (Coloured uphol. only with nylon)
☐ Black
☐ Red★ ☐ Orange★
☐ Light blue★ ☐ Green★
☐ Cornflower blue★

☐ **Not Detachable foot-rest** ④

☐ **Detachable foot-rest 82° angle**
 (DP 41 cm)
 (STANDARD)

☐ **Detachable foot-rest 68° angle**
 (DP 41 cm)

☐ **Elevating foot-rest (detachable)**

☐ **Plastic separate** ⑤

☐ Cover for footplates
(recommended for seat width >52cm)
☐ Foot-holding belts (cover plates N.A.)
☐ Heel stops (cover plates N.A.)

Separate aluminium profile tilt adjustable footplates: ⑥

☐ forward ☐ backward

☐ Cover for footplates
(recommended for seat width >52cm)
☐ Foot-holding belts (cover plates N.A.)

Single aluminium profile tilt adjustable footplate: ⑦
(limited range of width adjustment)
 (not available with elevating foot-rest)

☐ forward ☐ backward

☐ Foot-holding belts
☐ Heel stops

Side-guards ⑧

☐ Plastic

☐ Aluminium
(same colour)

Detachable DESK armrest - tip-up ⑨

☐ OPEN H22 with side-guard
(only if HP ≥44 cm)

☐ OPEN H27 with side-guard
(if HP ≤42cm)

☐ CLOSED H22
(from HP40/R24"x1 only if without any table fastenings)

☐ CLOSED H27 cm from seat
(compulsory for HP=42cm&R24")

Elbow-rest

☐ Classic (23 cm long)
☐ Straight (26 cm long)
☐ Tubular (only if open armrest)
☐ Hemiplegic RX
☐ Hemiplegic LX

note: see also Hemi modular armrest ord.form

Detachable SPORT armrest - tip-up ⑩

☐ OPEN H22 with side-guard
(only if HP ≥44 cm)

☐ OPEN H27 with side-guard
(if HP ≤42cm)

☐ CLOSED H22
(from HP40/R24"x1 only if without any table fastenings)

☐ CLOSED H27 cm from seat
(compulsory for HP=42cm&R24")

Elbow-rest

☐ Classic (23 cm long)
☐ Straight (26 cm long)
☐ Tubular (only if open armrest)
☐ Hemiplegic RX
☐ Hemiplegic LX

note: see also Hemi mod.armrest ord.form

Detachable "U" armrest - tip-up ⑪

☐ OPEN H22 with side-guard
(only if HP ≥44 cm)

☐ OPEN H27 with side-guard
(if HP ≤42cm)

☐ CLOSED H22
(from HP40/R24"x1 only if without any table fastenings)

☐ CLOSED H27 cm from seat
(compulsory for HP=42cm&R24")

Elbow-rest

☐ Straight (26 cm long)
☐ Straight (36 cm long)
☐ Tubular (only if open armrest)
☐ Hemiplegic RX
☐ Hemiplegic LX

note: see also Hemi mod.armrest ord.form

Height adjustable tip-up arm-rest (only closed) ⑫

☐ Desk model

☐ Sport model

☐ "U" model

Elbow-rest

☐ Classic (30 cm long)
☐ Classic (36 cm long)
☐ Straight (26 cm long)
☐ Straight (36 cm long)
☐ Hemiplegic RX
☐ Hemiplegic LX

note: see also Hemi mod.armrest ord.form

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____

PAGE 2/2
 May 2025

ALTHEA19PLUS

ALTHEA LOWERED SET-UP

Dealer _____

Rear wheels ⑬

- ☐ 24" 36 spokes ☐ 12" solid*
☐ 22" 36 spokes ☐ 12" pneumatic*
☐ 20" 36 spokes (N.A. if HP=42 if HP44: only curved rear frame)
☐ 24" MTB citybike tyre
☐ 24" MTB offroad tyre
☐ 24" Light (24 spokes)
☐ Pneum. Plugged (x 1" 3/8) **Tyres**
 Pneum. high press. ☐ Plugged ☐ Smooth
☐ Solid grey x 1" 3/8
☐ Solid black KIK 24"
☐ 24" Pneum. Marathon Plus

Pushrims ⑭

- ☐ Aluminium silver
☐ Chromium-plated steel 24"
☐ Aluminium dark
☐ Rubber coated
☐ OVAAL al.silver (*)
☐ OVAALS al.dark (*)
☐ OVAALG al.dark w/sil. insert (*)
☐ Carbolife CURVE L 24"(*)
☐ Carbolife GEKKO 24" (*)
☐ Silicon covers 24"
 (*) not available with MTB wheels

Spoke guards ⑮

- ☐ Transparent
☐ Black
☐ Fantasy
 code pattern:
 (see catalogue of spokeguards)

Front castors ⑯

- ☐ 100 (4"x1.1/4) solid - plastic rim
☐ 100 (4"x1.1/4) solid - alum. rim
☐ 125 (5"x1.1/4) solid - plastic rim
☐ 125 (5"x1.1/4) solid - alum. rim
☐ 150 (6") pneumatic
☐ 150 (6") solid
☐ 175 (7") solid
☐ 175 (7") pneumatic

Rear frame ⑰

- ☐ STRAIGHT rear frame with small wheels for narrow passages (default)
☐ STRAIGHT rear frame WITHOUT any small wheels f/narrow passages
☐ STRAIGHT rear frame with LEVER-activated narrow passage wheels
☐ CURVED rear frame
 Max user weight 100 kg
 (Only if rear seat/floor height HP is between 44 and 48 cm)
 (It is not possible to assemble the small wheels for narrow passages)

Unbalancing system ⑱

- ☐ Right ☐ Right
☐ Left ☐ Left
 (Only with straight rear frame) (Only with curved rear frame)

Anti-tip device (Only with straight rear frame)

- ☐ Right
☐ Left


Revolving anti-tip device (Only with curved rear frame)

- ☐ Right
☐ Left


Single drive ⑲

- ☐ Right
☐ Left

(camber can be only 0°)
 (rear wheels can be only 22" or 24")
 (the axle can grow +2cm or +4cm)



Extended wheel plate ⑳



Push handles ㉑


- ☐ Standard ☐ Height adjustable
☐ Hinged backrest (no stretched bar) (max L=50)
 (with HP=44 cm)

Distance between the hub centre and back-rest ㉒

- ☐ -2 cm (ext. wheel plate)
☐ 3,5 cm
☐ 5,5 cm
☐ 7,5 cm
 (with HP=44 cm)

Adjustable swing away Lateral supports ㉓

(tension adj. upholstery is compulsory)



- ☐ Right
☐ Left

Specify the height "h" from the center of the swing away to the seat
 cm

Head-rest ㉔

- Padded** (cannot grow in width)
☐ Fixed
☐ Height adjustable
Orbital
☐ Anatomic
☐ Moulded
☐ With lateral supports
 Specify "h" & "d" where requested
 Standard: h = 20 cm d = 15 cm
 h = ... cm d = ... cm





Table ㉕

- ☐ Transparent
☐ RxOLx Hemiplegic transp.
☐ Plastic grey
☐ Soft padded
 Specify the fastening:
☐ single ☐ double



Lined cushion in foam polyurethane ㉖



- ☐ 3 cm high
☐ 5 cm high

Brakes ㉗

- ☐ Classic - short handle
☐ Classic - long handle
☐ Drum
☐ Right single brake
☐ Left single brake
☐ Drum+classic

Options

- 29 ☐ "M+F" velcro on seat
 30 ☐ "M" velcro strip (hook) on seat
 31 ☐ "F" velcro strip (loop) on seat
 32 ☐ Seat/back-rest connection flap
 33 ☐ Chest harness
 34 ☐ Chest belt
 35 ☐ Abductor belt
 36 ☐ Front protective covers
 37 ☐ Stretched bar (width adjust.)
 38 ☐ R OL Crutch hold./single
 39 ☐ R OL Crutch hold./double
 40 ☐ Pump for pneumatics

Abductor ㉘



- ☐ 1 ☐ 3 ☐ 5 ☐ 7 cm
☐ cm

Choose or specify the height from seat upholstery

- 41 ☐ 45° pelvic belt
 42 ☐ 45° padded pelvic belt
 43 ☐ Black backpack
 44 ☐ IV bar + poles
 45 ☐ IV poles

- ☐ Quick rel. axles Ø 12 mm with 42mm cambered bushes (default)
☐ Quick rel. axles Ø 1/2" (12,7 mm) with 42mm cambered bushes
☐ Quick rel. axles Ø 1/2" (12,7 mm) with 50mm cambered bushes