

Dealer

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____ User's weight Kg _____

PAGE 1/2
 Oct. 2026


Key

☐ standard☐ option

- L** Seat width ☐34 ☐36 ☐38 ☐40 ☐42 ☐44 ☐46 ☐48 ☐50
☐52 ☐56 ☐60 (including eventual stretched bar, tension adj. upholstery & cover for footplate according to setup)
- P** Seat depth ☐36 ☐38 ☒40 ☐42 ☐44 ☐46
☐48 ☐50 (Each depth is adjustable +2, +4, +6 cm up to max. 50 cm)
- HS** Back-rest height ☐29 ☐31 ☐33 ☐35 ☐37 ☐39 ☐41 ☐43 ☐45
☐47 ☐49 ☐51 ☐55 ☐60 (stretched bar included)
- HA** Front seat/floor height ☐42 ☐43 ☐44 ☐46 ☐48 ☒50 ☐52
☐54
- HP** Rear seat/floor height ☐38 ☐40 ☐42 ☐44 ☒46 ☐48 ☐50
- DP** Foot-rest/seat distance ☐34 ☐35 ☐37 ☐39 ☐41 ☐43 ☐45 ☐47 ☐49 ☐51
 (from 17 to 32 cm foot-rest mounted on the frame, not possible with plastic footrests) cm

ALTHEA

Alum. foldable frame
 with double crossbar
 crash tested



Frame colour

- ☐ Metallic blue ☐ Pearl red ☐ Matte black
- ☐ Basic colour.....
- ☐ Special colour..... from:
 (see "Colours" catalogue)

①

Upholstery

- ☐ Tension adj. NYLON backrest uphol.
☐ Tens.adj. BREATHABL. backrest uphol.
☐ Without any back-rest upholstery
- ☐ With divided seat (joint by velcro - std)
☐ With one entire seat (not depth adjustable)

②

Upholstery colour

- ★ (Coloured uphol. only with nylon)
- ☐ Black
☐ Red★ ☐ Orange★
☐ Light blue★ ☐ Green★
☐ Cornflower blue★


③

☐ Not Detachable foot-rest☐ Detachable foot-rest
82° angle

(DP 41 cm)
 (82° footrests do not allow any 175 front castors)

☐ Detachable foot-rest
68° angle

(DP 41cm)

☐ Elevating foot-rest
(detachable)

④

☐ Plastic separate

(only W= from 34 to 56 cm)

- ☐ Cover f/footplates *
☐ Foot-holding belts (W50max)
☐ Heel stops (W50max)
 (* eventual cover included for W 52 to 60)

⑤

Separate aluminium profile
tilt adjustable footplates:

- ☐ forward ☐ backward
- ☐ Foot-holding belts (W50max)
 (* eventual cover included for W 52 to 60)

⑥

Single aluminium profile tilt
adjustable footplate:

(not available w/detachable or elevating footrests)
 (Note: it might interfere with DP)


- ☐ forward ☐ backward
- ☐ Foot-holding belts
☐ Heel stops

⑦

☐ AUTOMATIC closure
foot-plate in carbon

(up to W46 for W 48cm)
 (not available w/detachable or elevating footrests)
 (Note: it might interfere with DP)
 (max W 48cm)

⑧




☐ Foot-holding belts

Side-guards

⑨

☐ Plastic

☐ Aluminium
 (same colour)

Detachable DESK
armrest - tip-up

⑩

- ☐ OPEN H22
with side-guard
(only if HP ≥44 cm)
- ☐ OPEN H27
with side-guard
(if HP ≤42cm)

- ☐ CLOSED H22
(from HP40/R24"x1
only if without
any table fastenings)
- ☐ CLOSED H27 cm
from seat
(compulsory for
HP=42cm&R24")

Elbow-rest

- ☐ Classic (23 cm long)
☐ Straight (26 cm long)
☐ Tubular (only if open armrest)
☐ Hemiplegic RX
☐ Hemiplegic LX

note: see also Hemi modular armrest ord.form

Detachable SPORT
armrest - tip-up

⑪

- ☐ OPEN H22
with side-guard
(only if HP ≥44 cm)
- ☐ OPEN H27
with side-guard
(if HP ≤42cm)

- ☐ CLOSED H22
(from HP40/R24"x1
only if without
any table fastenings)
- ☐ CLOSED H27 cm
from seat
(compulsory for
HP=42cm&R24")

Elbow-rest

- ☐ Classic (23 cm long)
☐ Straight (26 cm long)
☐ Tubular (only if open armrest)
☐ Hemiplegic RX
☐ Hemiplegic LX

note: see also Hemi mod.armrest ord.form

Detachable "U"
armrest - tip-up

⑫

- ☐ OPEN H22
with side-guard
(only if HP ≥44 cm)
- ☐ OPEN H27
with side-guard
(if HP ≤42cm)

- ☐ CLOSED H22
(from HP40/R24"x1
only if without
any table fastenings)
- ☐ CLOSED H27 cm
from seat
(compulsory for
HP=42cm&R24")

Elbow-rest

- ☐ Straight (26 cm long)
☐ Hemiplegic RX
☐ Hemiplegic LX

note: see also Hemi mod.armrest ord.form

Height adjustable
tip-up armrest (only closed)

⑬

- ☐ Desk model
☐ Sport model
☐ "U" model

Elbow-rest

- ☐ Classic (30 cm long)
☐ Classic (36 cm long)
☐ Straight (26 cm long)
☐ Straight (36 cm long)
☐ Hemiplegic RX
☐ Hemiplegic LX

note: see also Hemi modular armrest ord.form

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____

PAGE 2/2
 Oct. 2026

ALTHEA

Dealer

Rear wheels

14

- ☐ 24" 36 spokes (def.) ☐ 12" solid*
☐ 22" 36 spokes ☐ 12" pneumatic*
☐ 20" 36 spokes (N.A. if HP=42 or HP≥46cm if HP44: only curved rear frame)
☐ 24" MTB citybike tyre
☐ 24" MTB offroad tyre
☐ 24" Light (24 spokes)

- ☐ Pneum. Plugged (x 1.3/8") **Tyres**
☐ Plugged ☐ Smooth Pneum. high press.
☐ Solid grey 24x1.3/8"/22"x1.3/8"
☐ Solid black KIK 24"
☐ 24" Pneum. Marathon Plus

Pushrims

15

- ☐ Aluminium silver
☐ Shiny inox 24"
☐ Aluminium dark
☐ Rubber coated
☐ Carbolife CURVE L 24"(*)
☐ OVAAL silver alum. (*)
☐ OVAALS dark alum. (*)
☐ OVAALG dark w/sil. insert (*)
☐ Carbolife GEKKO 24" (*)
☐ Silicon covers 24"
 (*) not available with MTB wheels

Spokeguards

16

- ☐ Transparent
☐ Black
☐ Fantasy
 code pattern:
 (see catalogue of spokeguards)

Castor wheels


17

- ☐ 100 (4"x1.1/4) solid - plastic rim
☐ 100 (4"x1.1/4) sol. silver alum. rim
☐ 100 (4"x1.1/4) sol. black alum. rim
☐ 125 (5"x1.1/4) solid - plastic rim
☐ 125 (5"x1.1/4) sol. silver alum. rim
☐ 125 (5"x1.1/4) sol. black alum. rim
☐ 150 (6") pneumatic
☐ 150 (6") solid
☐ 175 (7") solid **
☐ 175 (7") pneumatic **
 (**) only if HA ≥ 48cm and with 68° detach. footrests


Rear frame

18


- ☐ STRAIGHT rear frame with small wheels for narrow passages (default)



- ☐ STRAIGHT rear frame WITHOUT any small wheels f/narrow passages




- ☐ STRAIGHT rear frame with LEVER-activated narrow passage wheels



Max
user weight
100 kg

(Only if rear seat/floor height
HP is between 44 and 48 cm)

- ☐ CURVED rear frame



Unbalancing system

19

- ☐ Right ☐ Right
☐ Left ☐ Left
 (Only with straight rear frame) (Only with curved rear frame)

Anti-tip device (Only with straight rear frame)



- ☐ Right
☐ Left

Revolving anti-tip device (Only with curved rear frame)




- ☐ Right
☐ Left

Single drive

20


(only w/24" or 22single dr. rear w.)

- ☐ Right
☐ Left



Extended wheel plate


21



Push handles

22

- ☐ Standard ☐ Height adjustable(*) ☐ Hinged back-rest (no stretched bar)



COG Distance between the hub centre and backrest

23


- ☐ -2 cm (ext. wheel plate)
☐ 3,5 cm
☐ 5,5 cm
☐ 7,5 cm



(with HP=44 cm)

Adjustable swing away Lateral supports

23




- ☐ Right
☐ Left

Specify the height "h"
from the center of the
swing away to the seat
.....
cm

Head-rest

24

Backrest extension




- ☐ Fixed
☐ Height adjustable
 Default: h = 20 cm d = 15 cm

Orbital


- ☐ Anatomic



- ☐ Moulded



- ☐ With lateral supports



Table

25

- ☐ Transparent
☐ Rx OLx Hemiplegic transp.*
☐ Plastic grey
☐ Soft padded

Specify the fastening:
☐ single ☐ double



(Only with arm-rest) OXL if W>50cm
 (* only with reinforced armrest)

Lined cushion in foam polyurethane

28



- ☐ 3 cm high
☐ 5 cm high

Brakes

26


- ☐ Classic - short handle
☐ Classic - long handle
☐ Drum*
☐ Right single brake
☐ Left single brake
☐ Drum*+classic
 (* needed drum wheels)

Options

- 29 ☐ "F" velcro strip (loop) on seat
 30 ☐ "M" velcro strip (hook) on seat
 31 ☐ "M" velcro for Vicair
 32 ☐ Seat/back-rest connection flap
 33 ☐ Chest harness
 34 ☐ Chest belt
 35 ☐ Abductor belt
 36 ☐ Front protective covers
 37 ☐ Stretched bar
 38 ☐ Catheter holder pocket
 (note: only with entire seat)
 39 ☐ R ☐ L Crutch hold./single
 40 ☐ R ☐ L Crutch hold./double
 41 ☐ Pump for pneumatics

Abductor

27



- ☐ 1 ☐ 3 ☐ 5 ☐ 7 cm

- ☐ cm

Choose or specify the height from seat upholstery

- 42 ☐ Suitcase holder
 43 ☐ 45° pelvic belt
 44 ☐ 45° padded pelvic belt
 45 ☐ Black backpack
 46 ☐ IV bar + poles
 47 ☐ IV poles
 48 ☐ Padded rigid cushion(W52max)
 49 ☐ Padded rigid cushion(W56-60)
 50 ☐ Q.rel.axles Ø12mm w/42mm bushes (def.)
 51 ☐ Q.rel.a. Ø1/2" (12,7mm) w/42mm bushes
 52 ☐ Q.rel.a. Ø1/2" (12,7mm) w/50mm bushes
 53 ☐ Stickers for car fastenings
 (escl.: load>120Kg, R12", UltraRW, hinged bekr)