

☐ Order N° \_\_\_\_\_ ☐ Quotation  
Date \_\_\_\_\_ Signature \_\_\_\_\_  
Customer Ref. \_\_\_\_\_ User's weight Kg \_\_\_\_\_

**PAGE 1/2**  
Oct.2026


Dealer

- L** ☐ **RANGE 1 W** P > 36cm ☐ 34 ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46  
☐ **RANGE 2 W** P > 36cm ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48 ☐ 50 ☐ 52  
☐ **RANGE 3 W** P > 36cm (compulsory expandable stretched bar to be added) ☐ 46 ☐ 48 ☐ 50 ☐ 52 ☐ 54 ☐ 56 ☐ 58
- P** Seat depth ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46  
Adjustable up to 6 cm(max. 50 cm)
- HS** Back-rest height ☐ 29 ☐ 31 ☐ 33 ☐ 35 ☐ 37 ☐ 39 ☐ 41 ☐ 43 ☐ 45  
☐ 47 ☐ 49 ☐ 51 ☐ 55 ☐ 60 (stretched bar - that can grow - included)
- HA** Front seat/floor height ☐ 42 ☐ 43 ☐ 44 ☐ 46 ☐ 48 ☐ 50 ☐ 52  
☐ 54
- HP** Rear seat/floor height ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48 ☐ 50
- DP** Foot-rest/seat distance ☐ 34 ☐ 35 ☐ 37 ☐ 39 ☐ 41 ☐ 43 ☐ 45 ☐ 47 ☐ 49 ☐ 51  
(from 17 to 32 cm foot-rest mounted on the frame, not possible with plastic footrests)

Seat width (adjustable within the selected range)

# ALTHEA PLUS

Width & depth adjust.  
fold.alum.frame



- Frame colour** ☐ Met.blue ☐ Pearl red ☐ Met.grey ☐ Matte black  
☐ Pearl white ☐ Shiny black  
☐ Basic colours .....  
☐ Special colours .....  
(see "Colours" catalogue)

- Backrest upholstery** ☐ Tension adj. nylon. back  
☐ Tension adj. breath. back  
☐ Without any back-rest upholstery  
☐ With divided seat (joint by velcro, def.)  
Seat upholstery can grow +6cm max

- Upholstery colour** ★ (Coloured uphol. only with nylon)  
☐ Black  
☐ Red★ ☐ Orange★  
☐ Light blue★ ☐ Green★  
☐ Cornflower blue★

- ☐ **Not Detachable foot-rest** **(4)**
- ☐ **Detachable foot-rest 82° angle**  
(DP 41 cm)  
(82° footrests do not allow any 175 front castors) **(5)**
- ☐ **Detachable foot-rest 68° angle**  
(DP 41cm) **(6)**
- ☐ **Elevating foot-rest (detachable)** **(7)**

- ☐ **Plastic separate** **(5)**
- ☐ Cover for footplates (recommended for seat width >52cm)  
☐ Foot-holding belts (cover plates N.A.)  
☐ Heel stops (cover plates N.A.)

- Separate aluminium profile tilt adjustable footplates:** **(6)**
- ☐ forward ☐ backward
- ☐ Cover for footplates (recommended for seat width >52cm)  
☐ Foot-holding belts (cover plates N.A.)

- Single aluminium profile tilt adjustable footplate:** **(7)**  
(Limited range of width adjustment)  
(Not available with elevating footrests)  
(Note: it might interfere with DP)
- ☐ forward ☐ backward
- ☐ Foot-holding belts  
☐ Heel stops

- Side-guards** **(8)**
- ☐ Plastic   
☐ Aluminium (same colour)
- Detachable DESK armrest - tip-up** **(9)**
- ☐ OPEN H22 with side-guard (only if HP ≥44 cm)  
☐ OPEN H27 with side-guard (if HP≤42cm)  
☐ CLOSED H22 (from HP40/R24"x1 only if without any table fastenings)  
☐ CLOSED H27 cm from seat (compulsory for HP=42cm&R24")
- Elbow-rest**  
☐ Classic (23 cm long)  
☐ Straight (26 cm long)  
☐ Tubular (only if open armrest)  
☐ Hemiplegic RX  
☐ Hemiplegic LX  
note: see also Hemi modular armrest ord.form

- Detachable SPORT armrest - tip-up** **(10)**
- ☐ OPEN H22 with side-guard (only if HP ≥44 cm)  
☐ OPEN H27 with side-guard (if HP≤42cm)  
☐ CLOSED H22 (from HP40/R24"x1 only if without any table fastenings)  
☐ CLOSED H27 cm from seat (compulsory for HP=42cm&R24")
- Elbow-rest**  
☐ Classic (23 cm long)  
☐ Straight (26 cm long)  
☐ Tubular (only if open armrest)  
☐ Hemiplegic RX  
☐ Hemiplegic LX  
note: see also Hemi mod.armrest ord.form

- Detachable "U" armrest - tip-up** **(11)**
- ☐ OPEN H22 with side-guard (only if HP ≥44 cm)  
☐ OPEN H27 with side-guard (if HP≤42cm)  
☐ CLOSED H22 (from HP40/R24"x1 only if without any table fastenings)  
☐ CLOSED H27 cm from seat (compulsory for HP=42cm&R24")
- Elbow-rest**  
☐ Straight (26 cm long)  
☐ Straight (36 cm long)  
☐ Tubular (only if open armrest)  
☐ Hemiplegic RX  
☐ Hemiplegic LX  
note: see also Hemi mod.armrest ord.form

- Height adjustable tip-up arm-rest (only closed)** **(12)**
- ☐ Desk model   
☐ Sport model   
☐ "U" model
- Elbow-rest**  
☐ Classic (30 cm long)  
☐ Classic (36 cm long)  
☐ Straight (26 cm long)  
☐ Straight (36 cm long)  
☐ Hemiplegic RX  
☐ Hemiplegic LX  
note: see also Hemi mod.armrest ord.form

☐ Order N° \_\_\_\_\_ ☐ Quotation  
 Date \_\_\_\_\_ Signature \_\_\_\_\_  
 Customer Ref. \_\_\_\_\_

PAGE 2/2  
 Oct.2026

# ALTHEA PLUS

Dealer

## Rear wheels

13

- ☐ 24" 36 spokes ☐ 12" solid\*  
☐ 22" 36 spokes ☐ 12" pneumatic\*  
☐ 20" 36 spokes (N.A. if HP=42 or HP≥ 46 if HP44: only curved rear frame)  
☐ 24" MTB citybike tyre  
☐ 24" MTB offroad tyre  
☐ 24" Light (24 spokes)  
☐ Pneum. Plugged (x 1"3/8) **Tyres**  
 Pneum. high press. ☐ Plugged ☐ Smooth  
☐ Solid grey 24x1.3/8"/22"x1or1.3/8"  
☐ Solid black KIK 24"  
☐ 24"/22" Pneum. Marathon Plus

## Pushrims

14

- ☐ Aluminium silver  
☐ Shiny inox 24"  
☐ Aluminium dark  
☐ Rubber coated  
☐ OVAAL al.silver (\*)  
☐ OVAALS al.dark (\*)  
☐ OVAALG al.dark w/sil. insert (\*)  
☐ Carbolife CURVE L 24"(\*)  
☐ Carbolife GEKKO 24" (\*)  
☐ Silicon covers 24"  
 (\*) not available with MTB wheels

## Spoke guards

15

- ☐ Transparent  
☐ Black  
☐ Fantasy  
 code pattern: .....  
 (see catalogue of spokeguards)

## Front castors


16

- ☐ 100 (4"x1.1/4) solid - plastic rim  
☐ 100 (4"x1.1/4) sol.silver alum. rim  
☐ 100 (4"x1.1/4) sol.back alum. rim  
☐ 125 (5"x1.1/4) solid - plastic rim  
☐ 125 (5"x1.1/4) sol.silver alum. rim  
☐ 125 (5"x1.1/4) sol.black alum. rim  
☐ 150 (6") pneumatic  
☐ 150 (6") solid  
☐ 175 (7") solid \*\*  
☐ 175 (7") pneumatic \*\*  
 (\*\*) only if HA 48cm as minimum and 68"detach.footr.


## Rear frame

17


- ☐ STRAIGHT rear frame with small wheels for narrow passages (default)




- ☐ STRAIGHT rear frame WITHOUT any small wheels f/narrow passages



- ☐ STRAIGHT rear frame with LEVER-activated narrow passage wheels



- ☐ CURVED rear frame



(Only if rear seat/floor height HP is between 44 and 48 cm)


## Unbalancing system

18

- ☐ Right  
☐ Left  
 (Only with straight rear frame)

- ☐ Right  
☐ Left  
 (Only with curved rear frame)

## Anti-tip device (Only with straight rear frame)



- ☐ Right  
☐ Left

## Revolving anti-tip device (Only with curved rear frame)




- ☐ Right  
☐ Left

## Single drive

19


- ☐ Right  
☐ Left

(camber can be only 0°)  
 (rear wheels can be only 22" or 24")  
 (the axle can grow +2cm or +4cm)



## Extended wheel plate


20



## Push handles

21

- ☐ Standard ☐ Height adjustable(\*) ☐ Hinged back-rest (no stretched bar)



## Distance between the hub centre and backrest

22


- ☐ -2 cm (ext. wheel plate)  
☐ 3,5 cm  
☐ 5,5 cm  
☐ 7,5 cm



(with HP=44 cm)

## Adjustable swing away Lateral supports

23




- ☐ Right  
☐ Left

Specify the height "h" from the center of the swing away to the seat  
 ..... cm

## Head-rest

24

### Backrest extension (cannot grow in width)



- ☐ Fixed  
☐ Height adjustable  
 Default: h = 20 cm d = 15 cm

### Orbital


- ☐ Anatomic



- ☐ Moulded



- ☐ With lateral supports



## Table

25

- ☐ Transparent  
☐ Rx OLx Hemiplegic transp.  
☐ Plastic grey  
☐ Soft padded

Specify the fastening:  
☐ single ☐ double



(Only with arm-rest) if L>50cm

## Lined cushion in foam polyurethan (cannot grow in width)

28



- ☐ 3 cm high  
☐ 5 cm high

## Brakes

26


- ☐ Classic - short handle  
☐ Classic - long handle  
☐ Drum  
☐ Right single brake  
☐ Left single brake  
☐ Drum+classic

## Options

- 29 ☐ "M+F" (hook&loop) velcro strip on seat  
 30 ☐ "M" on seat ☐ "M" f/VICAIR on seat  
 31 ☐ "F" velcro strip (loop) on seat  
 32 ☐ Seat/backrest connection flap  
 33 ☐ Chest harness  
 34 ☐ Chest belt  
 35 ☐ Abductor belt  
 36 ☐ Front protective covers  
 37 ☐ Stretched bar (width adj.)  
 38 ☐ R ☐ OL Crutch hold./single  
 39 ☐ R ☐ OL Crutch hold./double  
 40 ☐ Pump for pneumatics

## Abductor

27



- ☐ 1 ☐ 3 ☐ 5 ☐ 7 cm  
☐ ..... cm

Choose or specify the height from seat cover

- 41 ☐ 45° pelvic belt  
 42 ☐ 45° padded pelvic belt  
 43 ☐ Black backpack  
 44 ☐ IV bar + poles  
 45 ☐ IV poles  
 46 ☐ Q.rel.axles Ø12mm w/42mm bushes (def.)  
 47 ☐ Q.rel.a. Ø1/2" (12,7mm) w/42mm bushes  
 48 ☐ Q.rel.a. Ø1/2" (12,7mm) w/50mm bushes