

☐ Order N ° \_\_\_\_\_ ☐ Quotation  
 Date \_\_\_\_\_ Signature \_\_\_\_\_  
 Customer ref. \_\_\_\_\_

PAGE 1/2  
 Jan 2024

Dealer \_\_\_\_\_


**Key**

☐ standard

☐ option

- L** Seat width ☐24 ☐26 ☐28 ☐30 ☐32 ☐34 ☐36
- P** Seat depth ☐26 ☐28 ☐30 ☐32 ☐34 ☐36
- HS** Back-rest height ☐21 ☐23 ☐25 ☐27 ☐29 ☐31 ☐33 ☐35 ☐37 ☐39
- HA** Front seat/floor height with MAX size of front castors (fc) ☐38(fc80max) ☐39(fc100max) ☐40(fc125max) ☐42(fc150max) ☐45(fc175max) ☐48  
☐44(fc150max) ☐46(fc175max)
- HP** Rear seat/floor height ☐34 ☐36 ☐38 ☐40 ☐42 ☐44 ☐46 ☐48  
 NOTE: 24" rear wheels if min.HP 38-max HP 48 22" if min.HP 36- max HP 45 20" if min.HP 34- max HP 43
- DP** Foot-rest/seat distance ☐25 ☐27 ☐29 ☐31 ☐33 ☐35 ☐37 ☐39 ☐41 ☐43 cm

# ARIEL



**Frame colour**

①

- ☐ Met. blue ☐ Pearl red ☐ Pastel yellow  
☐ Deep sky light blue ☐ Met. pink ☐ Matte black  
☐ Neon orange ☐ Neon green  
☐ ..... (see "Colours" catalogue)

**Upholstery**

②

- ☐ Tension adj. nylon back-rest  
☐ Tension adj. breathable backr.  
☐ Without back-rest  
☐ Without seat

**Uphostery colour**

③

- (nylon only)  
☐ Black ☐ Light blue  
☐ Red ☐ Green  
☐ Cornflower blue ☐ Orange

**Footrests**


④

- ☐ Not detachable footrests, **STRAIGHT** front frame (Compulsory pull to brake if P<30cm)  
☐ Detachable footrests, **STRAIGHT** front frame (Compulsory pull to brake)  
☐ **ABDUCTED** front frame with NOT detachable foot-rests (Need technical approval)  
☐ **ABDUCTED** front frame with detachable foot-rests (Need technical approval)  
☐ **Elevating footrests** (Need technical approval)  
 (Std brake starting from P32cm)

☐ **Separate aluminium profile foot-plates**

⑤

(from DP 17 to DP 31 cm the footplate is assembled on the frame)  
 (for L24-26-28cm recommended single footplate)




☐ Foot-holding belts

☐ **Single aluminium profile tilt adjust. footplate**

⑥

(from DP 17 to DP 31 cm the footplate is assembled on the frame)



☐ Foot-holding belts

**Rear wheel camber**


⑦

- ☐ 0° ☐ 1° ☐ 3° ☐ 5°

- ☐ Quick rel. axles Ø 12 mm (default)  
☐ Quick rel. axles Ø 1/2" (12,7 mm)  
☐ Cambered bushes length 42 mm (def.)  
☐ Cambered bushes length 50 mm

**Side guards**

☐ Plastic



☐ Aluminium (same colour)

☐ **Height adjustable "T" armrest (with plastic side-guards)**




**Elbow-rest**

- ☐ Classic (23 cm long)  
☐ Classic (30 cm long)  
☐ Straight (26 cm long)

☐ **Height adjustable "T" armrest (with aluminium side-guards)**

⑨




**Elbow-rest**

- ☐ Classic (23 cm long)  
☐ Classic (30 cm long)  
☐ Straight (26 cm long)

☐ **Tip-up height adjustable from 22 to 28cm "L" armrest**

⑩



**Elbow-rest**

- ☐ Tubular  
☐ Straight (26 cm long)  
☐ Classic (23 cm long)  
☐ Classic (30 cm long)

**Rear wheels**

⑫

- ☐ 24" (36 spokes) (\*)  
☐ 22" (36 spokes)  
☐ 20" (36 spokes)  
☐ 24" Light (24 spokes) (\*)  
☐ 22" Light (24 spokes)  
☐ 20" Light (20 spokes)  
☐ 12" (300) pneumatic (\*\*)  
☐ 12" (300) solid (\*\*)  
 (\*\*) need technical approval
- Tyres**  
☐ Pneum. HP smooth 20"/22"/24"  
☐ Pneum. HP plugged 22"/24"  
☐ Solid 24x1.3/8"  
☐ Solid KIK 24" black  
☐ Marathon Plus 22" or 24"

(\*) Only with P=36 or P34cm

**Pushrims**

⑬

- ☐ Aluminium silver  
☐ Chromium-plated steel 24"  
☐ Aluminium dark  
☐ Rubber coated  
☐ OVAAL alum. silver 24" (\*)  
☐ OVAALS alum. dark 24" (\*)  
☐ OVAALG dark w/silic. insert 24" (\*)  
☐ Carbolife CURVE L 24" (\*)  
☐ Carbolife GEKKO 24" (\*)  
☐ Silicon covers

(\*) not available with MTB wheels

**Spoke guards**

⑭

- ☐ Transparent  
☐ Black  
☐ Fantasy  
 code pattern: .....  
 (see catalogue of spokeguards)  
☐ Fantasy with name  
 code pattern: .....  
 (see catalogue of spokeguards)  
 name to be printed:  
 "....."

**Front castors**

⑮


- ☐ 80 (3") solid with lights  
☐ 80 (3") solid  
☐ 100 (4") solid with lights  
☐ 100 (4") solid  
☐ 100 (4") solid with alum. rim  
☐ 125 (5") solid  
☐ 125 (5") solid with lights  
☐ 125 (5") solid with alum. rim  
☐ 150 (6") pneumatic  
☐ 150 (6") solid

Dealer


## Push handles

17


☐ Standard



☐ Height adjustable




☐ Abducted height adjustable




## ☐ Curved rear frame


(it is not possible to assemble the small wheels for narrow passages)



## ☐ With straight rear frame



## ☐ With small wheels for narrow passages



(need technical approval)

## Single drive

(camber can be only 0°)  
 (rear wheels can be only 22" or 24")


21



☐ Right  
☐ Left

## ☐ Extended wheel plate

22



## Unbalancing system

20



☐ Right

☐ Left




☐ Right

☐ Left

(Only with curved rear frame)

## Anti-tip device



☐ Right

☐ Left

☐ Right & left

## Revolving anti-tip device

(Only with curved rear frame)



☐ Right


☐ Left

☐ Right & left

## Head-rest

23

Padded



☐ Height adjustable

Orbital


☐ Anatomic



☐ Moulded



☐ With lateral supports




Specify "h" & "d" where requested  
 Standard: h = 20 cm d = 15 cm

h = ..... cm d = ..... cm

## ☐ Abductor

24



☐ 1 ☐ 3 ☐ 5 ☐ 7 cm

☐ ..... cm

Choose or specify the height from the seat cover

## Anodized set

(Rear wheel rims & hub if Light and forks)

27

☐ Black

☐ Silver

☐ Light blue

☐ Red

☐ Orange

☐ Cyclamen light purple

(Light wheels only)

## Table

26

☐ Transparent

☐ Soft padded


Specify the fastening:

☐ single ☐ double

(Only with arm-rest)

## Adjustable lateral supports

25



☐ Right

☐ Left

Specify the height "h" from seat cover to the center of the swing away

## Lined cushion

in foam polyurethane

28



☐ 3 cm high

☐ 5 cm high

## Brakes

29

☐ Classic - short handle (from P32cm)

☐ Classic - pull to brake

☐ Classic long handle

☐ Brakes inside the sideguard\*

☐ Right single brake

☐ Left single brake

☐ Drum

☐ Drum + classic

\* compulsory for some configurations

## Options

- 30 ☐ "M" velcro (hook) on seat  
 31 ☐ "M" velcro for VICAIR cushion  
 32 ☐ "F" (loop) on seat  
 33 ☐ Seat/backrest connection flap  
 34 ☐ Padded rigid cushion  
 35 ☐ Chest harness  
 36 ☐ Chest belt  
 37 ☐ Abductor belt  
 38 ☐ Double cross bar \*  
 39 ☐ Stretched bar  
 40 ☐ Pump for pneumatics  
 41 ☐ 45° pelvic belt  
 42 ☐ 45° padded pelvic belt  
 43 ☐ Black backpack  
 44 ☐ Front protective covers

\* need technical approval

## Notes