

☐ Order N ° _____ ☐ Quotation
 Date _____ Signature _____
 Customer ref. _____

PAGE 1/2
 Oct.2026


Dealer _____

- L** Seat width (adjustable within the selected range) ☐ RANGE 1 W ☐ 24* ☐ 26* ☐ 28* ☐ 30* ☐ 32 ☐ RANGE 2 W ☐ 28* ☐ 30* ☐ 32 ☐ 34 ☐ 36
- * Recommended abducted front frame
- P** Seat depth ☐ 26 ☐ 28 ☐ 30 ☐ 32 ☐ 34 ☐ 36
- HS** Back-rest height ☐ 21 ☐ 23 ☐ 25 ☐ 27 ☐ 29 ☐ 31 ☐ 33 ☐ 35 ☐ 37 ☐ 39
- HA** Front seat/floor height with MAX size of front castors (fc) ☐ 38(fc80max) ☐ 39(fc100max) ☐ 40(fc125max) ☐ 42(fc150max) ☐ 45 ☐ 48 ☐ 44(fc150max) ☐ 46
- HP** Rear seat/floor height ☐ 34 ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48
 NOTE: 24" rear wheels if min.HP 38-max HP 48 22" if min.HP 36- max HP 45 20" if min.HP 34- max HP 43
- DP** Foot-rest/seat distance ☐ 25 ☐ 27 ☐ 29 ☐ 31 ☐ 33 ☐ 35 ☐ 37 ☐ 39 ☐ 41 cm

Key
☐ standard
☐ option

ARIEL PLUS

Width & depth adjust.
 KID fold.alum.frame



Frame colour

- ☐ Met. blue ☐ Pearl red ☐ Pastel yellow
☐ Deep sky light blue ☐ Met. pink ☐ Matte black
☐ Basic colours
☐ Special colours
 (see "Colours" catalogue)

①

Backrest upholstery

- ☐ Tension adj. nylon back-rest
☐ Tension adj. breathable backr.
☐ Without any back-rest upholstery

②

Uphostery colour

(Nylon only)


- ☐ Black ☐ Light blue
☐ Red ☐ Green
☐ Cornflower blue ☐ Orange

③

Straight front frame

Recommended for seat width L > 32 cm


- ☐ Not detachable footrests
☐ Detachable footrests
 (Compulsory pull to brake if depth < 30 cm)
☐ Elevating
 (Need technical approval)
 (Push brake starting from P32cm)



④

Separate aluminium profile footplates

(from DP 17 to DP 31 cm the footplate is assembled on the frame)
 (for W24-26-28cm recommended single footplate)




- ☐ Foot-holding belts

⑤

Single aluminium profile tilt adjust. footplate

(from DP 17 to DP 31 cm the footplate is assembled on the frame)



- ☐ Foot-holding belts

⑥

Rear wheel camber

- ☐ 0° ☐ 1° ☐ 3° ☐ 5°


- ☐ Q.rel.axles Ø12mm w/42mm bushes (def.)
☐ Q.rel.a. Ø1/2" (12,7mm) w/42mm bushes
☐ Q.rel.a. Ø1/2" (12,7mm) w/50mm bushes

⑦

Abducted front frame

Recommended for seat width W < 32 cm

- ☐ Abducted front frame with NOT detach.foot-rests
☐ Abducted front frame with detach.foot-rests
☐ Elevating
 (To be added to abducted frame)
 (Push brake starting from P32cm)
 (Need technical approval)



⑧

Side-guards and ev. integrated armrests

- ☐ Plastic side-guard ☐ Aluminium side-guard (frame colour)
☐ Height adjustable "T" armrest (with plastic side-guards) ☐ Height adjustable "T" armrest (with aluminium side-guards)




Elbow-rest

- ☐ Classic (23 cm long)
☐ Classic (30 cm long)
☐ Straight (26 cm long)

⑨

Tip-up height adjustable from 22 to 28cm "L" armrest

(sideguard is needed)



⑩

Elbow-rest

- ☐ Tubular
☐ Straight (26 cm long)
☐ Classic (23 cm long)
☐ Classic (30 cm long)

Rear wheels

⑪

- ☐ 24" (36 spokes) (*)
☐ 22" (36 spokes)
☐ 20" (36 spokes)
☐ 24" Light (24 spokes) (*)
☐ 22" Light (24 spokes)
☐ 20" Light (20 spokes)
☐ 12" (300) pneumatic (**)
☐ 12" (300) solid (**)
 (**) need technical approval
Tyres
☐ Pneum. HP smooth 20"/22"/24"
☐ Pneum. HP plugged 22"/24"
☐ Solid Grey 20"/22"x1/24"x1.3/8"
☐ Solid Black KIK 24"
☐ Marathon Plus 22" or 24"

(*) Only with P=36 or P34cm

Pushrims

⑫

- ☐ Aluminium silver
☐ Shiny inox 24"
☐ Aluminium dark
☐ Rubber coated
☐ OVAAL al.silver 24"
☐ OVAAL al.dark 24"
☐ OVAALG al.dark w/sil. insert 24"
☐ Carbolife CURVE L 24"
☐ Carbolife GEKKO 24"
☐ Silicon covers 24"

Spoke guards

⑬

- ☐ Trasparent
☐ Black
☐ Fantasy
 code pattern:
 (see catalogue of spokeguards)
☐ Fantasy with name
 code pattern:
 (see catalogue of spokeguards)
 name to be printed:
 "....."

Front castors

⑭


- ☐ 80 (3") solid with lights
☐ 80 (3") solid
☐ 100 (4") solid with lights
☐ 100 (4") solid
☐ 100 (4") solid w/silver alum.rim
☐ 100 (4") solid w/black alum.rim
☐ 125 (5") solid
☐ 125 (5") solid w/silver alum.rim
☐ 125 (5") solid w/black alum.rim
☐ 150 (6") pneumatic
☐ 150 (6") solid

Dealer


Push handles

15


☐ Standard




☐ Height adjustable




☐ Abducted height adjustable




☐ **Curved rear frame**
(it is not possible to assemble the small wheels for narrow passages)



☐ **With straight rear frame**



☐ **With small wheels for narrow passages**



(need technical approval)

Single drive

18


(camber can be only 0°)
(rear wheels can be only 22" or 24")
(the axle can grow +2cm or +4cm)



☐ Right
☐ Left

☐ **Extended wheel plate**

19



Unbalancing system

17

☐ Right
☐ Left
(Only with curved rear frame)

☐ Right
☐ Left
(Only with straight rear frame)

Revolving anti-tip device

(Only with curved rear frame)

☐ Right
☐ Left
☐ Right & left

Anti-tip device


(Only with straight rear frame)

☐ Right
☐ Left
☐ Right & left

Head-rest

20

Backrest extension
(cannot grow in width)



☐ Height adjustable
Default: h = 20 cm d = 15 cm


☐ **Abductor**

21

onto the frame for P >28cm



onto the seat for P <28cm




☐ 1 ☐ 3 ☐ 5 ☐ 7 cm
☐ cm

Choose or specify the height from the seat cover

Adjustable lateral supports

22



☐ Right
☐ Left

Specify the height "h" from the center of the swing away to the seat

Anodized set

23

(Rear wheel rims & hub if Light and forks)

☐ Black
☐ Silver
☐ Light blue
☐ Red
☐ Orange
☐ Cyclamen light purple

(Light wheels only)

Lined cushion

24

in foam polyurethane
(cannot grow in width)



☐ 3 cm high
☐ 5 cm high

Brakes

26

☐ Classic push (short handle) from P32cm
☐ Classic pull (short handle)
☐ Classic push long handle
☐ Classic pull long handle
☐ Brakes inside the sideguard*
☐ Right single brake
☐ Left single brake
☐ Drum
☐ Drum + classic

* compulsory for some configurations

Options


- 27 ☐ "M+F" velcro on seat
- 28 ☐ "M" velcro (hook) on seat
- 29 ☐ "M" velcro for VICAIR cushion
- 30 ☐ "F" (loop) on seat
- 31 ☐ Seat/backrest connection flap
- 32 ☐ Chest harness
- 33 ☐ Chest belt
- 34 ☐ Abductor belt
- 35 ☐ Stretched bar (width adj.)
- 36 ☐ Pump for pneumatics
- 37 ☐ 45° pelvic belt
- 38 ☐ 45° padded pelvic belt
- 39 ☐ Black backpack
- 40 ☐ Front protective covers

Orbital


☐ Anatomic



☐ Moulded



☐ With lateral supports



Table

25

☐ Transparent
☐ Soft padded

(Only with armrests and always with double fastening)

Notes