

Key ☐ standard ☐ option

EASYTREKKING

- Aluminium frame
- Takes up less space: the device can be disassembled
- Height adjustable footplate, with foot housings and foot-holding belts
- Single 20" Mountainbike with shock absorber
- Triple parking stands
- Front removable handle bars
- Removable and adjustable rear handle bars
- Seat width: 34 cm
- Seat depth: 40 cm
- Backrest height: 35 cm
- Anatomic seat and backrest in microfiber
- Tilt adjustable (4 positions) backrest
- Fix H20cm armrests with elbowrest
- 45° pelvic belt
- Brake
- Seat height from the floor 65 cm
- Weight 26 Kg
- Max load: 50 kg



Options:

- ☐ Pulling belt (code 8TK004)
- ☐ Chest harness (code 802/1)
- ☐ Moulded and H adjust.headrest (code 8TK001)

Sizes of the disassembled device:

- Width 42 cm (seat width 34 cm)
- Height 43 cm (backrest height 35 cm)
- Depth 78 cm (seat depth 40 cm)



Notes