

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____ User's weight Kg _____

PAGE 1/2
 Oct.2026*

Dealer

Key

☐ standard

☐ option

☐ **EOS Fix axle**


☐ **EOS Adj axle**

L Seat width (*) ☐ 34 ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46* ☐ 48*
P Seat depth (*) ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48* ☐ 50*
HS Back-rest height ☐ 27 ☐ 29 ☐ 31 ☐ 33 ☐ 35 ☐ 37 ☐ 39 ☐ 41
HA Front seat/floor height ☐ 46 ☐ 48 ☐ 50
HP Rear seat/floor height ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48
DP Foot-rest/seat distance ☐ 37 ☐ 39 ☐ 41 ☐ 43 ☐ 45 ☐ 47 ☐ 49 cm

(*) compulsory reinforced frame if the user's weight is between 100 and 120Kg

EOS

Titanium rigid frame



Frame colour

☐ Shiny titanium ☐ Natural titanium
☐ Shiny anthracite ☐ Matte black
☐ Glossy black
☐ Basic colours
☐ Special colours
 (see "Colours" catalogue)

①

Backrest: Upholstery

☐ Tension adj nylon ☐ Fix (stand.)
☐ Breathable tension adj ☐ Without seat
☐ Without any backrest uphol. ☐ Rigid in carbon fibre
☐ **FAMEKH** (separate order form)


Seat:

Backrest upholstery colour


☐ Black
☐ Coloured (nylon).....
 colours: orange, red, green, cornflower, light blue
☐ Coloured (breath).....
 colours: red, grey

③


Footplate




☐ **Tubular OPEN**
 titanium footplate
 (for Y or U shaped front frame)



☐ Carbon fibre
CLOSED
 titanium footplate
 (for Y or U shaped front frame)



☐ **TILT Adjustable**
 Carbon fibre footplate
 and titanium support
 (for any type of front frame)



☐ **Retro-fit FLIP UP**
footplate
 (aluminium footplate
 and titanium support)
 (for any type of front frame)

④


Push handles

☐ Without any pushing handles ☐ Height adj. pushing handles
☐ Fixed pushing handles ☐ Foldable pushing handles

⑤

Distance between the hub centre and backrest


⑥



☐ **9 cm (standard)** ☐ 7 cm
☐ 0 cm (biamputated) ☐ 11 cm
☐ 3 cm ☐ 12 cm
☐ 5 cm

Backrest tilt


⑦




☐ 89°
☐ **91° (standard)**
☐ 93°

Front frame angle

⑧



☐ **83° A = 14 cm (default)**



☐ **75° A = 19 cm**

Rear wheel camber

⑨

☐ 1° ☐ 2° ☐ 3° ☐ 4°
☐ **Quick release axles Ø12mm with 42 cambered bushes (default)**
☐ Q.r.a. Ø 1/2" (12,7 mm) w/42 camb.bushes
☐ Q.r.a. Ø 1/2" (12,7 mm) w/50 camb.bushes

Revolving anti-tip device


⑩




☐ Right
☐ Left

Side-guards

⑪



☐ Aluminium with carbon fibre wing (fix sideguard)



☐ Detachable with carbon f.wing (nat.tit. or matte black support)




☐ Alumin. with straight carbon fibre



☐ Detachable w/straight crbn.sideg.

☐ **Height adjustable from 22 to 28cm "L" arm-rest, tip up, detachable**



Elbow-rest

☐ Tubular
☐ Straight (26 cm long)
☐ Classic (23 cm long)
☐ Classic (30 cm long)

Brakes

⑬



☐ Scissor brake



☐ Composite brake

☐ Standard

☐ Pull brake




Pockets

⑭



☐ Seat pocket



☐ Big under the seat




☐ Back-rest pocket

Backrest bar

☐ Standard: dia 19mm
 P (bar) 7,5cm;
 H (bar) 19cm.

☐ Custom:
 P (bar).....cm;
 H (bar).....cm



Notes

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____

PAGE 2/2
 Oct.2026

EOS

Dealer _____

Rear wheels

16

- ☐ 24" Ultralight (24 spokes)
☐ 24" Light (24 spokes)
☐ 24" Star 6x3 (18 spokes, silver or black)
☐ 24" (36 spokes)
☐ 24" Mountain-bike with Citybike tyre
☐ 24" Mountain-bike with Offroad tyre
☐ 24" Spinergy LX
 spokes colour.....(black def.)
☐ 25" Spinergy LX
 spokes colour.....(black def.)
☐ 25" Light (24 spokes) black

Tyres

- ☐ Pneum. Marathon Plus 24"/25"
☐ Pneum. high pressure Smooth
☐ Pneum. high pressure Plugged 24"
☐ Pn.Schwalbe Right Run 24"
☐ Solid black KIK 24"
 not available for Ultralight wheels

Pushrims

17

- ☐ Aluminium silver
☐ Aluminium dark
☐ OVAAL alum. silver 24"
☐ OVAAL alum. dark 24"
☐ OVAALG(dark w/silicon ins.)24"
☐ Shiny stainless steel 24"
☐ Rubber coated
☐ Titanium
☐ Carbolife CURVE L24"
☐ CL Curve L TETRAGRIP 24"
☐ CL QUADRO rubber ctd 24"
☐ Carbolife GEKKO 24"
☐ Silicon covers

OVAAL/S/G only in size 24"
 Only OVAAL/GEKKO available for Star6x3
 OVAAL & Carbolife not available for MTB

Spokeguards

18

- ☐ Transparent 24"
☐ Black 24"
☐ Fantasy
 code pattern
 (see catalogue of spokeguards)

Anodized set

(Rear wheel hub & rims if LIGHT 24",
 sideguard supports, forks and fittings)

- ☐ Black
☐ Silver
☐ Light blue
☐ Red
☐ Orange
☐ Cyclamen light purple

20

Castor wheels

19

- ☐ 80 (3") solid
☐ 80 (3") solid with lights

☐ 100 (4"x1.1/4) solid with plastic rim
☐ 100 (4"x1) solid with lights
☐ 100 (4"x1.1/4) solid with silver alum. rim
☐ 100 (4"x1.1/4) solid with black alum. rim


☐ 125 (5"x1.1/4) solid with plastic rim
☐ 125 (5"x1.1/4) solid with silver alum. rim
☐ 125 (5"x1.1/4) solid with black alum. rim

Lined cushion in foam polyurethan




- ☐ 3 cm high
☐ 5 cm high

21



(carbon fibre rigid seat only suitable for standard seat)



22

Options

- 26 ☐ "M+F" velcro strips (hook and loop) on seat
 27 ☐ "M" velcro strips (hook) on seat
 28 ☐ "M" velcro strips for Vicair
 29 ☐ "F" velcro strips (loop) on seat
 30 ☐ Chest harness
 31 ☐ Chest belt
 32 ☐ Abductor belt
 33 ☐ Catheter holder pocket
 34 ☐ Rx ☐ Lx Crutch holder
 35 ☐ Pump for pneumatics
 36 ☐ Suitcase holder
 37 ☐ 45° pelvic belt
 38 ☐ 45° padded pelvic belt
 39 ☐ Black backpack
 41 ☐ Front protective covers


7,5 kg configuration (FIX) 7,8 kg configuration (ADJ)

- Wheelchair 38x40x31 cm
 - Nylon upholstery
 - Titanium frame colour: Shiny Titanium
 - Open foot-rest
 - Without pushing handles
 - Aluminium side guards with carbon fibre wing
 - Scissor brakes
 - Ultralight wheels
 - Smooth high press. pneum. tyres
 - Ultralight pushrim
 - 80 mm front castors
 - Without pockets
 - Max load: 100 kg

25

TAPERED SEAT (L-12cm)


23



(Tapered seat is possible only if the seat width is not lower than 36cm)


Custom front opening (i.e. custom TAPERED seat) cm

Tapered seat requires straight front frame:




STRAIGHT SEAT (def.)

24




Y SHAPED FRONT FRAME (default)



Custom Y/V front opening cm


- ☐ Reinforced frame (max load 100-120 Kg)

V SHAPED FRONT FRAME



(The tilt adjustable footplate only, or the flip up one, are suitable with the V shaped front frame)

STRAIGHT FRONT FRAME



(Recommended if seat width is 34cm)
 (Compulsory if the seat is tapered)

Notes