

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____ User's weight Kg _____

PAGE 1/2
 Oct.2026

Dealer _____

Key


☐ standard

☐ option



☐ EOS³ Fix axle

☐ EOS³ Adj axle

L Seat width ☐ 34 ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44
 P Seat depth ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46
 HS Back-rest height ☐ 27 ☐ 29 ☐ 31 ☐ 33 ☐ 35 ☐ 37 ☐ 39 ☐ 41
 HA Front seat/floor height ☐ 46 ☐ 48 ☐ 50
 HP Rear seat/floor height ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48
 DP Foot-rest/seat distance ☐ 37 ☐ 39 ☐ 41 ☐ 43 ☐ 45 ☐ 47 ☐ 49 cm



EOS³
HINGED front frame
 Rigid frame in Titanium



Frame colour

①

☐ Shiny titanium ☐ Natural titanium

note: frame of Eos3 cannot be varnished

Backrest:

Upholstery

②

☐ Tension adj nylon

Seat:

☐ Breathable tension adj

☐ Fixed (stand.)

☐ Without any backrest uph. ☐ Rigid in Carbon fibre

☐ **AMEKH** (separate order form)

Backrest upholstery colour

③


☐ Black

☐ Coloured (nylon).....
 colours: orange, red, green, cornflower, light blue

☐ Coloured (breath.).....
 colours: red, grey)


Tubular single titanium footplate

④



☐ Open ☐ Carbon fibre closed ☐ Tilt adjustable Carbon fib.closed

Pushing handles



☐ Without any pushing handles


☐ Fixed pushing handles

☐ Height adjustable pushing handles

☐ Foldable pushing handles

Distance between the hub centre and back-rest (COG)


⑥



☐ 9 cm (standard) ☐ 7 cm
☐ 0 cm (biamputated) ☐ 11 cm
☐ 3 cm ☐ 12 cm
☐ 5 cm

Backrest tilt


⑦




☐ 89°
☐ 91° (standard)
☐ 93°

Front frame angle

⑧



☐ 83° A = 14 cm (default)



☐ 75° A = 19 cm

Rear wheel camber

⑨

☐ 1° ☐ 2° ☐ 3° ☐ 4°

☐ Quick release axles Ø12mm with 42 cambered bushes (default)
☐ Q.r.a. Ø 1/2" (12,7 mm) w/42 camb.bushes
☐ Q.r.a. Ø 1/2" (12,7 mm) w/50 camb.bushes

Revolving anti-tip device

⑩



☐ Right
☐ Left

Side-guards

⑪



☐ Aluminium with carbon fibre wing

☐ Detachable CRB wing



☐ Aluminium w/straight carbon fibre

☐ Detachable straight CRB sideguard

☐ Height adjustable from 22 to 28cm "L" arm-rest, tip up, detachable

⑫



Elbow-rest

☐ Tubular
☐ Straight (26 cm long)
☐ Classic (23 cm long)
☐ Classic (30 cm long)

Brakes

⑬

☐ Scissor brake

☐ Composite brake



Note: the model of brake might change due to technical conflicts folding the front frame

Pockets

⑭




☐ Seat pocket

☐ Back-rest pocket

Backrest bar

⑮



☐ Standard: dia 19mm P(bar) 7,5cm; H(bar) 19cm.
☐ Custom: P(bar).....cm; H(bar).....cm.

Notes

Dealer

Rear wheels

- ☐ 24" Ultralight (24 spokes)
- ☐ 24" Light (24 spokes)
- ☐ 24" Star 6x3 (18 spokes, silver or black)
- ☐ 24" (36 spokes)
- ☐ 24" Mountain-bike with Citybike tyre
- ☐ 24" Mountain-bike with Offroad tyre
- ☐ 24" Spinergy LX
spokes colour.....(black def.)
- ☐ 25" Spinergy LX
spokes colour.....(black def.)
- ☐ 25" Light (24 spokes) black

Tyres

- ☐ Pneum. Marathon Plus
 - ☐ Pneum. high pressure Smooth
 - ☐ Pneum. high pressure Plugged 24"
 - ☐ Pn.Schwalbe Right Run 24"
 - ☒ Solid black KIK 24"
- not available for Ultralight wheels

Pushrim

- ☐ Aluminium silver
 - ☐ Aluminium dark
 - ☐ OVAAL alum. silver 24"
 - ☐ OVAAL alum. dark 24"
 - ☐ OVAALG(dark w/silicon ins.)24"
 - ☐ Shiny stainless steel 24"
 - ☐ Rubber coated
 - ☐ Titanium
 - ☐ Carbolife CURVE L
 - ☐ CL Curve L TETRAGRIP 24"
 - ☐ CL QUADRO rubber ctd 24"
 - ☐ Carbolife GEKKO
 - ☐ Silicon covers
- OVAAL/S/G only in size 24"
- Only OVAAL/GEKKO available for Star6x3
- OVAAL & Carbolife not available for MTB

Spoke guards

- ☐ Transparent 24"
☐ Black 24"
☐ Fantasy
 code pattern

(see catalogue of spokequards)

Anodized set

(Rear wheel hub & rims if LIGHT 24",
sideguard supports, forks and fittings)



- ☐ Black
☐ Silver
☐ Light blue
☐ Red
☐ Orange
☐ Cyclamen light purple

Front castors

- ☐ 80 (3") solid
- ☐ 80 (3") solid with lights
- ☐ 100 (4"x1.1/4) solid with plastic rim
- ☐ 100 (4"x1) solid with lights
- ☐ 100 (4"x1.1/4) solid with silver alum. rim
- ☐ 100 (4"x1.1/4) solid with black alum. rim
- ☐ 125 (5"x1.1/4) solid with plastic rim
- ☐ 125 (5"x1.1/4) solid with silver alum. rim
- ☐ 125 (5"x1.1/4) solid with black alum. rim


Lined cushion
in foam polyurethan

- ☐ 3 cm high
- ☐ 5 cm high

STANDARD
SEAT (def.) ☐


ERGO SEAT

Note: Carbon fibre rigid seat is suitable only for standard seat



☐ **TAPERED SEAT**
(default)

(TAPERED seat is the default choice for any seat width starting from 36cm onwards)



☐ STRAIGHT SEAT (24)

(Upon request the seat can be straight.
For a seat width of 34cm
the STRAIGHT seat is
the default choice)

Options


- 26 ☐ "M+F" velcro strips (hook and loop) on seat
- 27a ☐ "M" velcro strips (hook) on seat
- 27b ☐ "M" velcro strips for Vicair on seat
- 28 ☐ "F" velcro strips (loop) on seat
- 29 ☐ Catheter holder pocket
- 30 ☐ Pump for pneumatics
- 31 ☐ 45° pelvic belt
- 32 ☐ 45° padded pelvic belt
- 33 ☐ Black backpack
- 34 ☐ Front protective covers
- 35 ☐ Transport bag for EOS3
- 36 ☐ Bag for 2 rear wheels 24"

8,5 kg configuration (FIX)
8,8 kg configuration (ADJ)

- Wheelchair 38x40x31 cm
- Nylon upholstery
- Titanium frame colour:
Shiny Titanium
- Open foot-rest
- Without pushing handles
- Aluminium side guards
with carbon fibre wing
- Scissor brakes
- Ultralight wheels
- Smooth high press.
pneum. tyres
- Ultralight pushrim
- 80 mm front castors
- Without pockets
- Max load: 100 kg

Notes

- ☐ STRAIGHT FRONT FRAME
(default)



- ☐ **Custom front opening**
(i.e. custom TAPERED seat)
.....cm