

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____ User's weight Kg _____

PAGE 1/2
 Oct.2026


Dealer _____

Key ☐ standard ☐ option

- L** Seat width ☐ 34 ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44
- P** Seat depth ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46
- HS** Backrest height ☐ 29 ☐ 31 ☐ 33 ☐ 35 ☐ 37 ☐ 39 ☐ 41
- HA** Front seat/floor height ☐ 46 ☐ 48 ☐ 50
- HP** Rear seat/floor height ☐ 40 ☐ 42 ☐ 44 ☐ 46
- DP** Footrest/seat distance ☐ 35 ☐ 37 ☐ 39 ☐ 41 ☐ 43 ☐ 45 ☐ 47 cm

(max load 100Kg)

HALLEY.TI TITANIUM fold.frame



Frame colour

- ☐ Shiny titanium ☐ Natural titanium
☐ Glossy black ☐ Matte black
☐ Basic colours
☐ Special colours
 (see "Colours" catalogue)

①

Upholstery

- ☐ Tension adjustable nylon back-rest
☐ Tens.adjustable breathable backrest
☐ Without any backrest upholstery
☐ **FAMEKH** custom postural backrest
 (fill up the proper order form)

②



Upholstery colour

- ☐ Black
☐ Coloured.....
 Available colours:
 NYLON: orange, red, green, cornflower blue, light blue;
 BREATHABLE: red, grey.

③

Front frame shape and angle


- ☐ Y shaped (default) ☐ Straight



④

AUTOMATIC closure foot-plate in carbon fibre

(Note: it might interfere with DP)





- ☐ Foot-holding belts

⑤

Separate aluminium profile tilt adjustable foot-plates



- ☐ forward ☐ backward



- ☐ Foot-holding belts

⑥


- ☐ Front frame angle (with DP=43 cm and P=40 cm) (default) ☐ Extended front frame angle (with DP=43 cm and P=40 cm)



(approx.measurements)

Closed carbon fibre single foot-plate

(Note: it might interfere with DP)



- ☐ Foot-holding belts

⑦

Single aluminium profile tilt adjustable foot-plate

(Note: it might interfere with DP)

- ☐ forward ☐ backward




- ☐ up & outward

- ☐ Foot-holding belts


⑧

Side guards

- ☒ Carbon fibre



- ☐ Tip-up, detachable sport arm-rest h. 22 cm in aluminium




- ☐ in titanium

Elbow-rest

- ☐ Tubular
☐ Classic (23 cm long)

⑩

- ☐ Tip-up height adjustable from 22 to 28cm "L" arm-rest (detachable, open with sideguards)



Elbow-rest

- ☐ Tubular
☐ Straight (26 cm long)
☐ Classic (23 cm long)
☐ Classic (30 cm long)

⑪

Rear wheels

- ☐ 24" Ultralight (24 spokes) black
☐ 24" Light (24 spokes)
☐ 24" Star 6x3 (18 sps) silver or black
☐ 24" (36 spokes)
☐ 24" Mountain-bike CITYBIKE tyre
☐ 24" Mountain-bike OFFROAD tyre
☐ 24" Spinergy LX
 spokes colour.....(black def.)

Tyres

- ☐ Pneumatic Marathon Plus
☐ Pneum. high press.plugged
☐ Pneum. high press. smooth
☐ Solid black KIK 24"
 (not avail.with Ultralight)

⑫

Pushrims

- ☐ Aluminium silver
☐ Aluminium dark
☐ OVAAL 24" aluminium silver **
☐ OVAALS 24" aluminium dark **
☐ OVAALG 24" alum.dark w/silicon insert **
☐ Titanium
☐ Rubber coated
☐ Carbolife CURVE L 24"***
☐ CL CurveL TETRAGRIP 24"***
☐ CL QUADRO rubber ctd 24"***
☐ Carbolife GEKKO 24"***
☐ Silicon covers

OVAAL/S/G only in size 24"
 Only OVAAL/GEKKO available for Star6x3
 ** OVAAL & Carbolife not available for MTB

⑬

Order N° _____ ☐ Quotation
 Date _____ Signature _____
☐ Customer _____

PAGE 2/2
 Oct.2026

HALLEY.TI

Dealer

Spoke guards

14

- ☐ Transparent 24"
- ☐ Black 24"
- ☐ Fantasy
- code pattern
- (see catalogue of spokeguards)

Front castors

15

- ☐ 100 (4"x1.1/4) black solid plastic rim
- ☐ 100 (4"x1.1/4) solid with silver alum. rim
- ☐ 100 (4"x1.1/4) solid with black alum. rim
- ☐ 125 (5"x1.1/4) black solid plastic rim
- ☐ 125 (5"x1.1/4) solid with silver alum. rim
- ☐ 125 (5"x1.1/4) solid with black alum. rim
- ☐ 150 (6") pneumatic
- ☐ 150 (6") black solid

Revolving anti-tip device

16

- ☐ Right
- ☐ Left




Unbalancing Pedal

- ☐ Right
- ☐ Left



Back-rest pushing handles

17

- ☐ Standard ☐ Without
- ☐ Height adjustable ☐ Foldable pushing handles
- 
- 
- 
- 

Lined cushion in foam polyurethane

18



- ☐ 3 cm high
- ☐ 5 cm high

Anodized set

(Rear wheel rims & hubs only if 24" Light rear wheels, forks and fittings)

19

- ☐ Black
- ☐ Silver
- ☐ Light blue
- ☐ Red
- ☐ Orange
- ☐ Cyclamen light purple

(24" Light wheels only)

Brakes

20

- ☐ Classic - short handle
- ☐ Classic - pull to brake short handle
- ☐ Classic - long handle
- ☐ Scissor
- ☐ Composite

Opzioni

- 23 ☐ "M+F" velcro strip on seat
- 24 ☐ "M" velcro strip on seat
- 25 ☐ "M" velcro for Vicair
- 26 ☐ "F" velcro strip (loop) on seat
- 27 ☐ Seat/backrest connection flap
- 28 ☐ Abductor belt
- 29 ☐ Catheter holder pocket
- 30 ☐ Back-rest stretched bar

Opzioni

- 31 ☐ Rx ☐ Lx Crutch hold./single
- 32 ☐ Rx ☐ Lx Crutch hold./double
- 33 ☐ Pump for pneumatics
- 34 ☐ Suitcase holder
- 35 ☐ 45° padded pelvic belt
- 36 ☐ Black backpack
- 37 ☐ Front protective covers


Rear wheel camber

21

- ☐ 1° ☐ 2° ☐ 3° ☐ 4°
- ☐ Q.rel.axles Ø12mm w/42mm bushes (def.)
- ☐ Q.rel.a. Ø1/2" (12,7mm) w/42mm bushes
- ☐ Q.rel.a. Ø1/2" (12,7mm) w/50mm bushes

Distance between the hub and the backrest

22



- ☐ cautious setup (*)
- ☐ medium setup (*)
- ☐ very active setup (*)

(*) these values might change according to the selected HP rear seat height

Notes