

☐ Order N° _____ ☐ Quotation
Date _____ Signature _____
Customer Ref. _____ User's weight Kg _____

PAGE 1/2
Oct.2026

Dealer

Key

☐ standard



☐ option

L Seat width	☐ 34	☐ 35	☐ 36	☐ 37	☐ 38	☐ 39	☐ 40	☐ 41
	☐ 42	☐ 43	☐ 44	☐ 45	☐ 46	☐ 47	☐ 48	
P Seat depth	☐ 34	☐ 36	☐ 38	☐ 40	☐ 42	☐ 44	☐ 46	☐ 48
	☐ 50							
HS Back-rest height	☐ 27	☐ 29	☐ 31	☐ 33	☐ 35	☐ 37	☐ 39	
HA Front seat/floor height	☐ 40	☐ 41	☐ 42	☐ 43	☐ 44	☐ 45	☐ 46	
	☐ 47	☐ 48	☐ 49	☐ 50	☐ 51	☐ 52		
HP Rear seat/floor height	☐ 37	☐ 38	☐ 39	☐ 40	☐ 41	☐ 42	☐ 43	
	☐ 44	☐ 45	☐ 46	☐ 47	☐ 48	☐ 49	☐ 50	
DP Foot-rest/seat distance	☐ 35	☐ 36	☐ 37	☐ 38	☐ 39	☐ 40	☐ 41	
	☐ 42	☐ 43	☐ 44	☐ 45	☐ 46	☐ 47		

cm
(max load 100kg)

HEKA

Ergonomic Titanium custom rig.frame



Frame colour (including forks if varnished)

- ☐ **Natural titanium** ☐ **Shiny titanium**
☐ Matte black ☐ Glossy black
☐ Matte purple ☐ Metallic blue
☐ Matte London smoke ☐ Metallic grey
☐ Tiffany ☐ Champagne
☐ Matte green ☐ Pearl white
☐ with front stripes: ☐ black ☐ white ☐ red ☐ blue
(Note: not all colour matches are suitable / any selection needs approval)
☐ Colour (see "Colour catalogue").....
☐ RAL/Chameleon col./Camouflage.....shiny or matte, from:

①

Backrest: Upholstery

- ☐ Tension adjustable nylon
☐ Coloured (nylon).....
colours: orange, red, green, cornflower, light blue
☐ Breathable tension adjustable
☐ Coloured (breath.colours: red, grey).....
☐ Without any backrest upholstery

☐ **AMEKH** (separate order form)

Seat:

- ☐ Fix (stand.)
☐ Rigid in carbon fibre

②

Footplate

- ☐ Tubular single in titanium **OPEN** (def. backwards) (only for Y or U shaped front frame)
☐ Tubular single in titanium **CLOSED** in carbon (def. backwards) (only for Y or U shaped front frame)
☐ Closed **tilt adjustable** for any front frame (V, Y or U shaped) (carbon plate, titanium support)
☐ **Retro-fit liftable footplate** (in aluminium and titanium support) for any front frame (V, Y or U shaped)

☐ event. custom feature (to be quoted)

③


Pushing handles

- ☐ Without any pushing handles
☐ Fixed pushing handles
☐ Height adjustable pushing handles
☐ Foldable pushing handles

④

Distance between the hub centre and back-rest COG
(range from 0cm to 15cm max)


- ☐ cm 0
☐ cm.....
☐ cm 15



⑤


Backrest type & tilt

- ☐ 89°
☐ 91° (default)
☐ 93°
☐ **FOLDABLE** and tilt adjustable backrest (side-guards can be fix or removable)
☐ **WELDED**, fix backrest (side-guards can be only fix)



⑥

Front frame angle



- ☐ **83°** (default)
A = 16 cm, with DP 41 cm
☐ 75°
A = 19 cm, with DP 41 cm

⑦

Rear wheel camber

- ☐ 1° ☐ 2° ☐ 3° ☐ 4°
☐ **Quick release axles Ø12mm with 42 cambered bushes** (default)
☐ Quick release axles Ø 1/2" (12,7 mm) with 42 mm cambered bushes
☐ Quick release axles Ø 1/2" (12,7 mm) with 50mm cambered bushes

⑧

Side-guards (and event. armrests)

Carbon fiber side-guards with FOLDABLE backrest

- ☐ fix with CRB wing & matte black support (def.)
☐ fix with straight CRB
☐ detachable with CRB wing
☐ detachable with straight CRB



Carbon fiber side-guards with WELDED backrest

- ☐ fix with CRB wing (def.)
☐ fix with straight CRB



- ☐ Height adj., removable "L" armrest



- ☐ Sport, open, removable, custom armrest



Elbowrest

- ☐ Tubular
☐ Straight (26 cm)
☐ Classic (23 cm)

⑨

Revolving anti-tip device




- ☐ Right
☐ Left

⑩

Brakes

- ☐ composite
☐ Scissor
☐ standard



⑪

Notes

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____

PAGE 2/2
 Oct.2026

HEKA

Dealer

Rear wheels

12

- ☐ 24" Ultralight (24 spokes)
☐ 24" Light (24 spks) ** see section 15 f/colours
☐ 24" Star 6x3 (18 spokes, silver or black)
☐ 24" (36 spokes)
☐ 24" Mountain-bike with Citybike tyre
☐ 24" Mountain-bike with Offroad tyre
☐ 24" Spinergy LX
 spokes colour.....(black def.)
☐ 25" Spinergy LX
 spokes colour.....(black def.)
☐ 25" Light (24 spokes) black

Tyres

- ☐ Pneum. Marathon Plus
☐ Pneum. high pressure Smooth
☐ Pneum. high pressure Plugged 24"
☐ Pn.Schwalbe Right Run 24"
☐ Solid black KIK 24"
 not available for Ultralight wheels

Pushrims

13

- ☐ Aluminium silver
☐ Aluminium dark
☐ OVAAL alum. silver 24"
☐ OVAAL alum. dark 24"
☐ OVAALG(dark w/silicon ins.)24"
☐ Rubber coated
☐ Titanium
☐ Carbolife CURVE L24"
☐ CL Curve L TETRAGRIP 24"
☐ CL QUADRO rubber ctd 24"
☐ Carbolife GEKKO 24"
☐ Silicon covers

OVAAL/S/G only in size 24"
 Only OVAAL/GEKKO available for Star6x3
 OVAAL & Carbolife not available for MTB

Spokeguards

14

- ☐ Transparent 24"
☐ Black 24"
☐ Fantasy
 code pattern
 (see catalogue of spokeguards)

Anodized set (backrest plate forks if anodized & Light wheels 24")

15

- ☐ Black
☐ Silver
☐ Light blue
☐ Red
☐ Orange
☐ Cyclamen light purple

Front castors

16

- ☐ 80 (3") solid
☐ 80 (3") solid with lights
☐ 100 (4"x1.1/4) solid with plastic rim
☐ 100 (4"x1) solid with lights
☐ 100 (4"x1.1/4) solid silver alum. rim
☐ 100 (4"x1.1/4) solid black alum. rim
☐ 125 (5"x1.1/4) solid with plastic rim
☐ 125 (5"x1.1/4) solid silver alum. rim
☐ 125 (5"x1.1/4) solid black alum. rim

Lined cushion in foam polyurethan

17

- ☐ 3 cm high
☐ 5 cm high



(carbon fibre rigid seat only suitable for standard seat)



18

Inclined backrest bar

21

- ☐ Standard: dia 22mm
 P (bar) 7,5cm;
 H (bar) 19cm.

- ☐ Custom:
 H (bar).....cm


Pockets

19



- ☐ Small under the seat


or



- ☐ Hanging from seat



- ☐ Big under the seat




- ☐ Backrest pocket


Options

- 22 ☐ "M+F" velcro strips (hook and loop) on seat
 23 ☐ "M" velcro strips (hook) on seat
 24 ☐ "M" velcro for Vicair
 25 ☐ "F" velcro strips (loop) on seat
 26 ☐ Catheter holder pocket
 27 ☐ Pump for pneumatics
 28 ☐ 45° pelvic belt
 29 ☐ 45° padded pelvic belt
 30 ☐ Black backpack
 31 ☐ Front protective covers

- ☐ Straight seat (def.)



- ☐ Tapered seat
 (front frame can only be straight)




- ☐ Custom tapered seat
 cm

20

Front frame shape


21

- ☐ V shaped front frame




The CRB tilt adjustable footplate
 only is suitable with the V shaped
 front frame
 or eventually the retro-fit liftable one

- ☐ STRAIGHT
 front frame



Straight front frame or Y shaped front frame are
 suitable with any type of footplate

- ☐ Y shaped
 front frame



- ☐ Custom Y front opening cm
☐ Custom V front opening cm

Notes: