

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____ User's weight Kg _____

PAGE 1/2
 May 2025

Dealer

Key

☐ standard

☐ option

☐ IDRA2.0 Fix axle



☐ IDRA 2.0 Adj axle

L Seat width ☐ 34 ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44
P Seat depth ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46
HS Back-rest height ☐ 27 ☐ 29 ☐ 31 ☐ 33 ☐ 35 ☐ 37 ☐ 39 ☐ 41
HA Front seat/floor height ☐ 46 ☐ 48 ☐ 50
HP Rear seat/floor height ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48
DP Foot-rest/seat distance ☐ 37 ☐ 39 ☐ 41 ☐ 43 ☐ 45 ☐ 47 ☐ 49

cm

IDRA2.0

Aluminium rigid frame



Frame colour

☐ Glossy black ☐ Matte black

☐ Pearl white ☐ Pearl red

☐ (see "Colours" catalogue) from:

☐ Neon green ☐ Neon orange

Upholstery

Back-rest:

☐ Tension adj nylon

☐ Breathable tension adj

☐ Without any backrest upholstery

☐ **FAMEKH** (separate order form)

Seat:

☐ Fix (standard, screws)

☐ Tension adj. (velcro)




Backrest upholstery colour

☐ Black

☐ Coloured (nylon).....
 colours: orange, red, green, cornflower, light blue





☐ Coloured (breathable).....
 colours: red, grey)

Single tubular footplate




☐ Open ☐ Carbon fibre closed ☐ Tilt adjustable Carbon fibre closed

Push handles




☐ Without pushing handles ☐ Fixed pushing handles ☐ Height adjustable pushing handles ☐ Foldable pushing handles

Distance between the hub centre and backrest




☐ 9 cm (standard) ☐ 7 cm
☐ 0 cm (biamputated) ☐ 11 cm
☐ 3 cm ☐ 13 cm
☐ 5 cm

Backrest tilt



☐ 89°
☐ 91° (standard)
☐ 93°

Front angle



☐ 83° (standard) (A = 16 cm)
☐ 75° (A = 21 cm)

Rear wheel camber

☐ 1° ☐ 2° ☐ 3° ☐ 4°

☐ Quick rel. axles Ø 12mm with cambered bushes length 42 mm (default)
☐ Quick rel. axles Ø 1/2" (12,7mm) w/42mm bushes
☐ Quick rel. axles Ø 1/2" (12,7mm) w/50mm bushes

Side guards



☐ Aluminium with carbon fibre wing ☐ Aluminium with straight carbon fibre


☐ Height adjustable from 22 to 28cm "L" arm-rest, detachable and tip-up




Elbow-rest

☐ Tubular
☐ Straight (26 cm long)
☐ Classic (23 cm long)
☐ Classic (30 cm long)

Brakes




☐ Scissor brake ☐ Composite



☐ Standard ☐ Pull brake

Anti-tip device



☐ Right ☐ Left

Notes

