

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____ User's weight Kg _____

PAGE 1/2
 Oct.2026

Dealer _____

Key

☐ standard

☐ option

☐ **QUASAR Fix axle**

☐ **QUASAR Adj axle**

L Seat width ☐ 34 ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46*☐ 48*

P Seat depth ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48*☐ 50*

HS Back-rest height ☐ 27 ☐ 29 ☐ 31 ☐ 33 ☐ 35 ☐ 37 ☐ 39 ☐ 41

HA Front seat/floor height ☐ 46 ☐ 48 ☐ 50

HP Rear seat/floor height ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48


DP Foot-rest/seat distance ☐ 37 ☐ 39 ☐ 41 ☐ 43 ☐ 45 ☐ 47 ☐ 49

(*) compulsory reinforced frame if the user's weight is between 100 and 120Kg

cm

QUASAR

Titanium rigid frame



Frame colour

①

☐ Shiny titanium ☐ Natural titanium

☐ Shiny anthracite ☐ Matte black

☐ Glossy black

☐ Basic colours

☐ Special colours
 (see "Colours" catalogue)

Back-rest: Upholstery

☐ Tension adj nylon

☐ Breathable tension adj

☐ Without any back-rest uph.

☐ **FAMEKH**

custom postural backrest

(fill up the proper order form)

Seat:

☐ Fix (standard, screws)

☐ Tension adj. (velcro)

☐ Carbon fibre

Backrest upholstery colour


③

☐ Black

☐ Coloured (nylon).....
 colours: orange, red, green, cornflower, light blue

☐ Coloured (breathable).....
 colours: red, grey)

Footplate



☐ **Tub.OPEN**
 titanium
 (for Y or U
 shaped
 front frame)

☐ Carbon fibre
CLOSED tit.
 (for Y or U shaped
 front frame)

☐ **TILT Adjustable**
 Carbon fib.tit.
 (for any type of
 front frame)

☐ **Retro-fit FLIP UP**
 (aluminium footplate
 and titanium support)
 (for any type of front frame)

☐ Separate aluminium
 footplates
 with welded rear arch
 (for Y or U shaped
 front frame)

☐ Single aluminium
 footplates
 with welded rear arch
 (for Y or U shaped
 front frame)

Push handles

☐ Without
 any pushing
 handles


☐ Fixed
 pushing
 handles

☐ Height adj.
 pushing handles

☐ Foldable
 pushing handles

**Distance between
 the hub centre and back-rest**


⑥



☐ **9 cm (standard)** ☐ 7 cm
☐ 0 cm (biamputated) ☐ 11 cm
☐ 3 cm ☐ 13 cm
☐ 5 cm

Backrest tilt


⑦



☐ 89°
☐ **91° (standard)**
☐ 93°

Front angle

⑧



☐ **83° (standard) (A = 16 cm)**
☐ 75° (A = 21 cm)

**Rear
 wheel camber**


⑨

☐ 1° ☐ 2°
☐ 3° ☐ 4°

☐ Quick rel. axles Ø 12mm
 with cambered bushes
 length 42 mm (default)
☐ Quick rel. axles Ø 1/2"
 (12,7mm) w/42mm bushes
☐ Quick rel. axles Ø 1/2"
 (12,7mm) w/50mm bushes

Side-guards

⑩



☐ Aluminium with
 carbon fibre wing

☐ Detachable
 with carbon f.wing
 (nat.tit. or matte black support)

☐ Aluminium with
 straight CRB

☐ Detachable
 w/straight crbn.sideg.

Arm-rests

⑪

☐ **"L" arm-rest**
 detach.H22-28cm adj.
 tip-up

☐ **SPORT arm-rest**
 open,
 tip-up

Elbow-rest

☐ Tubular
☐ Straight (26 cm long)
☐ Classic (23 cm long)
☐ Classic (30 cm long) only f/L armrest

Brakes


⑫

☐ Scissor brake
☐ Composite

☐ Standard
☐ Pull brake

Anti-tip device

⑬



☐ Right
☐ Left

Notes

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____

PAGE 2/2
 Oct.2026

QUASAR

Dealer

Rear wheels

14

- ☐ 24" Ultralight (24 spokes)
☐ 24" Light (24 spokes)
☐ 24" Star 6x3 (18 spokes, silver or black)
☐ 24" (36 spokes)
☐ 24" Mountain-bike with Citybike tyre
☐ 24" Mountain-bike with Offroad tyre
☐ 24" Spinergy LX
 spokes colour.....(black def.)
☐ 25" Spinergy LX
 spokes colour.....(black def.)
☐ 25" Light (24 spokes) black

Tyres

- ☐ Pneum. Marathon Plus 24"/25"
☐ Pneum. high pressure Smooth
☐ Pneum. high pressure Plugged 24"
☐ Pn.Schwalbe Right Run 24"
☐ Solid black KIK 24"
 not available for Ultralight wheels

Pushrims

15

- ☐ Aluminium silver
☐ Aluminium dark
☐ Shiny stainless steel 24"
☐ Rubber coated
☐ Titanium
☐ OVAAL (silver aluminium)
☐ OVAALS (dark aluminium)
☐ OVAALG (w/sil.insert)
☐ CarboLife CURVE L 24" *
☐ CL CurveL TETRAGRIP rubb.c.24*
☐ CL QUADRO rubber ctd 24*
☐ CarboLife GEKKO 24" *
☐ Silicon covers 24"
 OVAAL/S/G only in size 24"
 Only OVAAL/GEKKO available for Star63
 OVAAL & Carbolife not available for MTB


Spokes guards

16

- ☐ Transparent 24"
☐ Black 24"
☐ Fantasy
 code pattern
 (see catalogue of spokeguards)

Lined cushion in foam polyurethan

18



- ☐ 3 cm high
☐ 5 cm high


Front castors

17


- ☐ 80 (3") solid
☐ 80 (3") solid with lights
☐ 100 (4"x1.1/4) solid with plastic rim
☐ 100 (4"x1) solid with lights
☐ 100 (4"x1.1/4) solid with silver alum. rim
☐ 100 (4"x1.1/4) solid with black alum. rim
☐ 125 (5"x1.1/4) solid with plastic rim
☐ 125 (5"x1.1/4) solid with silver alum. rim
☐ 125 (5"x1.1/4) solid with black alum. rim
☐ 150 (6") pneumatic
☐ 150 (6") solid

Pockets

19



- ☐ Small under the fix seat
 or
☐ Small hanging from the seat



- ☐ Big under the seat



- ☐ On the back-rest

Anodized set

(Rear wheel rims & hubs if LIGHT 24",
 sideguard supports, forks and fittings)


20

- ☐ Black
☐ Silver
☐ Light blue
☐ Red
☐ Orange
☐ Cyclamen light purple

(only 24" Light wheels)

Backrest bar

21




- ☐ Standard:
 P(bar) 7,5cm;
 H(bar) 19cm;
 dia 19mm.
☐ Custom:
 P.....cm;
 H.....cm;

Options

- 23 ☐ "M+F" velcro strips (hook and loop) on seat
 24 ☐ "M" velcro strips (hook) on seat
 25 ☐ "M" velcro strips for VICAIR
 26 ☐ "F" velcro strips (loop) on seat
 27 ☐ Chest harness
 28 ☐ Chest belt
 29 ☐ Abductor belt
 30 ☐ Catheter holder pocket
 31 ☐ Rx ☐ Lx Crutch holder
 32 ☐ Pump for pneumatics
 33 ☐ Suitcase holder
 34 ☐ 45° pelvic belt
 35 ☐ 45° padded pelvic belt
 36 ☐ Black backpack
 38 ☐ Front protective covers

STANDARD SEAT


* Carbon fibre rigid seat
 (only suitable for standard seat)



ERGO SEAT




"Y" SHAPED FRONT FRAME (def.)




(suitable for any type of footplate)

"V" SHAPED FRONT FRAME



(suitable only for TILT adj. or FLIP UP footplate)

STRAIGHT FRONT FRAME



(recommended for seat width 34cm)
 (suitable for any type of footplate)

- ☐ Custom Y or V shaped front opening.....cm
☐ Reinforced frame (max load 100-120 Kg)

7,0 kg configuration (FIXED) 7,4 kg configuration (ADJ.)

- Wheelchair 38x40x31 cm
 - Nylon upholstery
 - Tit. frame colour: Shiny titanium
 - Open foot-rest
 - No pushing handles
 - Aluminium side guards with carbon fibre wing
 - Scissor brakes
 - Ultralight wheels
 - Pneum. high press. smooth tyres
 - Ultralight handrims
 - 80 mm front castors
 - No pockets
 Max load: 100 kg