

Key

☐ standard

☐ option

☐ **QUASAR GT Fix axle**


☐ **QUASAR GT Adj axle**

L Seat width ☐ 34 ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44
P Seat depth ☐ 36 ☐ 38 ☒ 40 ☐ 42 ☐ 44 ☐ 46
HS Back-rest height ☐ 27 ☐ 29 ☐ 31 ☐ 33 ☐ 35 ☐ 37 ☐ 39 ☐ 41
HA Front seat/floor height ☐ 46 ☐ 48 ☐ 50
HP Rear seat/floor height ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48
DP Foot-rest/seat distance ☐ 37 ☐ 39 ☐ 41 ☐ 43 ☐ 45 ☐ 47 ☐ 49 cm

crash tested

QUASAR GT

Titanium rigid frame



Frame colour

- ☐ Shiny titanium ☐ Natural titanium
☐ Shiny anthracite ☐ Matte black
☐ Glossy black
☐ Basic colours
☐ Special colours
(see "Colours" catalogue)

①

Back-rest: Upholstery

- ☐ Tension adj nylon
☐ Breathable tension adj
☐ Without any backrest uph.
☐ **FAMEKH**
Seat:
☐ Fix (standard, screws)
☐ Tension adj. (velcro)
☐ Carbon fibre




②

Backrest upholstery colour

- ☐ Black
☐ Coloured (nylon).....
colours: orange, red, green, cornflower, light blue
☐ Coloured (breathable).....
colours: red, grey

③

Footplate




④

Push handles

- ☐ Without any pushing handles
☐ Fixed pushing handles
☐ Height adj. pushing handles
☐ Foldable pushing handles

⑤


Distance between the hub centre and back-rest



- ☐ **9 cm (standard)** ☐ 7 cm
☐ 0 cm (biamputated) ☐ 11 cm
☐ 3 cm ☐ 13 cm
☐ 5 cm

⑥


Backrest tilt



- ☐ 89°
☐ **91° (standard)**
☐ 93°

⑦

Front angle



- ☐ **83° (standard) (A = 16 cm)**
☐ 75° (A = 21 cm)



⑧

Rear wheel camber

- ☐ 1° ☐ 2°
☐ 3° ☐ 4°
☐ Quick rel. axles Ø 12mm with cambered bushes length 42 mm (default)
☐ Quick rel. axles Ø 1/2" (12,7mm) w/42mm bushes
☐ Quick rel. axles Ø 1/2" (12,7mm) w/50mm bushes

⑨

Side-guards




⑩

☒ **Security belt (included)**

☒ **Car fastening stickers (included)**


⑪

Brakes



⑫

Anti-tip device



- ☐ Right
☐ Left

⑬

Notes

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____

PAGE 2/2
 Oct.2026

QUASAR GT

Dealer _____

Rear wheels

14

- ☐ 24" Ultralight (24 spokes)
☐ 24" Light (24 spokes)
☐ 24" Star 6x3 (18 spokes, silver or black rim)
☐ 24" (36 spokes)
☐ 24" Mountain-bike with Citybike tyre
☐ 24" Mountain-bike with Offroad tyre
☐ 24" Spinerly LX
 spokes colour.....(black def.)
☐ 25" Spinerly LX
 spokes colour.....(black def.)
☐ 25" Light (24 spokes) black

Tyres

- ☐ Pneum. Marathon Plus
☐ Pneum. high pressure Smooth
☐ Pneum. high pressure Plugged 24"
☐ Pn.Schwalbe Right Run 24"
☐ Solid black KIK 24"
 not available for Ultralight wheels

Pushrims

15

- ☐ Aluminium silver
☐ Aluminium dark
☐ Shiny stainless steel 24"
☐ Rubber coated
☐ Titanium
☐ OVAAL (silver aluminium)
☐ OVAALS (dark aluminium)
☐ OVAALG (w/sil.insert)
☐ CarboLife CURVE L 24" *
☐ CL Curvel TETRAGRIP rubb.c.24"
☐ CL QUADRO rubber ctd 24"*
☐ CarboLife GEKKO 24" *
☐ Silicon covers 24"

* All OVAAL & CL pushrims: not suitable for MTB wheels
 * Only OVAAL & GEKKO suitable for Star 6x3 wheels


Spokes guards

16

- ☐ Transparent 24"
☐ Black 24"
☐ Fantasy
 code pattern
 (see catalogue of spokeguards)

Lined cushion in foam polyurethan

18



- ☐ 3 cm high
☐ 5 cm high


Front castors

17

- ☐ 80 (3") solid
☐ 80 (3") solid with lights
☐ 100 (4"x1.1/4) solid with plastic rim
☐ 100 (4"x1) solid with lights
☐ 100 (4"x1.1/4) solid with silver alum. rim
☐ 100 (4"x1.1/4) solid with black alum. rim
☐ 125 (5"x1.1/4) solid with plastic rim
☐ 125 (5"x1.1/4) solid with silver alum. rim
☐ 125 (5"x1.1/4) solid with black alum. rim
☐ 150 (6") pneumatic
☐ 150 (6") solid

Pockets

19



- ☐ Small under the fix seat
 or
☐ Small hanging from the seat
☐ Big under the seat
☐ On the back-rest

Anodized set

(Rear wheel rims & hubs if LIGHT 24",
 sideguard supports, forks and fittings)


20

- ☐ Black
☐ Silver
☐ Light blue
☐ Red
☐ Orange
☐ Cyclamen light purple

(only 24" Light wheels)

Backrest bar

21




- ☐ Standard:
 P(bar) 7,5cm;
 H(bar) 19cm;
 dia 19mm.
☐ Custom:
 P.....cm;
 H.....cm;

Options

- 23 ☐ "M+F" velcro strips (hook and loop) on seat
 24 ☐ "M" velcro strips (hook) on seat
 25 ☐ "F" velcro strips (loop) on seat
 26 ☐ Padded rigid cushion h 5
 27 ☐ Chest harness
 28 ☐ Chest belt
 29 ☐ Abductor belt
 30 ☐ Catheter holder pocket
 31 ☐ Rx ☐ Lx Crutch holder
 32 ☐ Pump for pneumatics
 33 ☐ Black backpack
 34 ☐ Front protective covers

STANDARD SEAT


* Carbon fibre rigid seat
 (only suitable for standard seat)



ERGO SEAT




"Y" SHAPED FRONT FRAME (def.)




(suitable for any type of footplate)

"V" SHAPED FRONT FRAME



(suitable only for TILT adj. or FLIP UP footplate)

STRAIGHT FRONT FRAME



(recommended for seat width 34cm)

22

7,0 kg configuration (FIXED) 7,4 kg configuration (ADJ.)

- Wheelchair 38x40x31 cm
 - Nylon upholstery
 - Tit. frame colour: Shiny titanium
 - Open foot-rest
 - No pushing handles
 - Aluminium side guards with carbon fibre wing
 - Scissor brakes
 - Ultralight wheels
 - Pneum. high press. smooth tyres
 - Ultralight handrims
 - 80 mm front castors
 - No pockets
 Max load: 100 kg

☐ Custom Y or V shaped front opening.....cm