

Dealer

Date:
☐ Quotation
☐ Order N°:
Customer Ref.:
2026 2027 · Signature

SAMEKH by



Alumium shell colour:

only for ALUMINIUM shell

- ☐ Matte black (def.)
☐ Colour

see "Colours" catalogue

POSTURAL BACK SUPPORT:

with **CARBON FIBRE** shell

- ☐ ACTIVE..... 6 cm

with **ALUMINIUM** shell

- ☐ ACTIVE contour..... 6 cm
☐ DEEP contour..... 9 cm
☐ EXTRA DEEP contour..... 18 cm

ACTIVE contour allows more movement

DEEP contour more ergonomic to better answer the accomodation demands on either the frontal and sagittal plane


- ☐ Shell length: ☐ 16 cm ☐ 24 cm ☐ 30 cm ☐ 37 cm

SIZE (* sizes available with carbon fibre shell):


Samekh SIZE	Wheelchair WIDTH
26	26 - 28
28	28 - 30
30	30 - 32
32	32 - 34 - 36
34*	34 - 36 - 38
36*	36 - 38 - 40
38*	38 - 40 - 42
40*	40 - 42 - 44
42*	42 - 44 - 46
44*	44 - 46 - 48
46	46 - 48 - 50

Code & Serial No.....
(to be filled up by OFFCARR staff)

- ☐ **Single central fastening** (H adjustable ± 3.5 cm)



- ☐ **Lateral fastenings** (H and width adjustable)



Lateral extensions: (only for DEEP contour)
(included in the EXTRA DEEP contour)

- ☐ Right Extension
☐ Left Extension

Swing Away (only for aluminium shell)

- ☐ Right lat.support / size XS S M L
☐ Left lat.support / size XS S M L

☐ **Privacy flap** 8SK012

Assembling notes

Usually (def.) assembled at 7 cm over the seat upholstery or 6 cm over the cushion. If ordered in combination with a wheelchair, pls eventually specify different assembling requests: