



☐ Order N° _____ ☐ Quotation
Date _____ Signature _____
Customer ref. _____ User's weight Kg _____

PAGE 1/2
May 2025

Dealer


Key ☐ standard ☐ option

L Seat width ☐ 34 ☐ 35 ☐ 36 ☐ 37 ☐ 38 ☐ 39 ☐ 40 ☐ 41 ☐ 42 ☐ 43 ☐ 44 ☐ 45 ☐ 46 ☐ 47 ☐ 48
P Seat depth ☐ 34 ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48 ☐ 50
HS Backrest height ☐ 27 ☐ 29 ☐ 31 ☐ 35 ☐ 37 ☐ 39
HA Front seat/floor H ☐ 40 ☐ 41 ☐ 42 ☐ 43 ☐ 44 ☐ 45 ☐ 46 ☐ 47 ☐ 48 ☐ 49 ☐ 50 ☐ 51 ☐ 52
HP Rear seta/floor H ☐ 37 ☐ 38 ☐ 39 ☐ 40 ☐ 41 ☐ 42 ☐ 43 ☐ 44 ☐ 45 ☐ 46 ☐ 47 ☐ 48 ☐ 49 ☐ 50
DP Footrest/seat dist. ☐ 35 ☐ 36 ☐ 37 ☐ 38 ☐ 39 ☐ 40 ☐ 41 ☐ 42 ☐ 43 ☐ 44 ☐ 45 ☐ 46 ☐ 47
cm
(max load 100kg)

Colour:
- **frame** ☐ Glossy black ☐ Matte black
- **front forks** ☐ Met. blue ☐ Pearl white
- **footplate** ☐ Colour (see "Colour catalogue")..... starting from:
- **sideguard supp.** ☐ RAL/Changing col./Camouflage.....
.....shiny or matte, from:

Upholstery
☐ Tension adj. nylon backrest
☐ Tension adj. breathable backrest
☐ Upholst.colour: Black (default)
Seat
☐ Nylon black (default)
☐ Carbon fibre

THEMIS
PRESTIGE
Ergonomic Aluminium
custom rigid frame



Front frame

☐ **V shaped**
(allows only tilt adj. footplate)
☐ **Straight**

Footplate

☐ Tilt adjustable
Carbon fib.closed
(available with V shaped or straight front frame)
☐ Open
(available only with straight front frame)
☐ Carbon fibre closed
(available only with straight front frame)

Push handles

☐ Without any pushing handles
☐ Fixed pushing handles
☐ Foldable pushing handles
☐ Height adjustable pushing handles

COG Distance between the hub centre and back-rest

☐ cm 0
☐ cm.....
☐ cm 15

Backrest tilt

☐ 89°
☐ 91° (def.)
☐ 93°
☐ Foldable and tilt adjustable backrest (default)
☐ Fixed, welded backrest

Front frame angle

☐ 83° (default) A = 16 cm
☐ 75° A = 19 cm, with DP 41 cm

Carbon fiber side-guards with FOLDABLE backrest:


☐ Fix with CRB wing & matte black support (def.)
☐ Fix w/STRAIGHT CRB&matte black supp.(def.)
☐ Detachable w/CRB wing & m.black supp.(def.)
☐ Det.w/STRAIGHT CRB&m.black supp.(def.)

Carbon fiber side-guards with WELDED backrest:

☐ fixed with CRB wing (def.)
☐ Fix w/STRAIGHT CRB

Backrest bar

☐ Default: P(bar depth) 7,5cm; H(bar height) 19cm; dia 21mm.
☐ Custom: P.....cm; H.....cm.



Brakes

☐ Composite brake
☐ Scissor brake
☐ Classic push brake

Revolving anti-tip

☐ Right
☐ Left

Rear wheel camber ☐ 1° ☐ 2° ☐ 3° ☐ 4°

☐ Quick release axles Ø 12mm with cambered bushes 42mm long (def.)
☐ Quick release axles Ø 1/2" (12,7mm) with cambered bushes 42mm long
☐ Quick release axles Ø 1/2" (12,7mm) with cambered bushes 50mm long

Dealer

Rear wheels

12

- ☐ 24" Ultralight (24 spokes)
☐ 24" Light (24 spks) ** see section 18 f/colours
☐ 24" Star 6x3 (18 spokes, silver or black)
☐ 24" (36 spokes)
☐ 24" Mountain-bike with Citybike tyre
☐ 24" Mountain-bike with Offroad tyre
☐ 24" Spinergy LX
spokes colour.....(black def.)
☐ 25" Spinergy LX
spokes colour.....(black def.)
☐ 25" Light (24 spokes) black

Tyres

- ☐ Pneum. Marathon Plus
☐ Pneum. high pressure Smooth
☐ Pneum. high pressure Plugged 24"
☐ Pn.Schwalbe Right Run 24"
☐ Solid KIK 24" not available for
black Ultralight wheels

Pushrims

13

- ☐ Aluminium silver 24"
☐ Aluminium dark 24"/25"
☐ OVAAL aluminium silver 24" *
☐ OVAALS aluminium dark 24" *
☐ OVAALG black w/sil.insert 24" *
☐ Rubber coated
☐ Titanium
☐ CarboLife CURVE L 24" *
☐ CL CurveL TetraGrip rubb.c.24" *
☐ CL Quadro rubber ctd 24" *
☐ CarboLife GEKKO 24" *
☐ Silicon covers 24"

* not available for Star 6x3 or MTB wheels
*CL pushrims: only GEKKO suitable for Star 6x3 wheels

Spokes guards

14

- ☐ Transparent 24"
☐ Black 24"
☐ Fantasy
code pattern
(see catalogue of spokeguards)


Front castors

15

- ☐ 80 (3") solid black W 25 mm
☐ 80 (3") solid with lights W 25 mm
☐ 100 (4"x1.1/4) black w/plastic rim
☐ 100 (4"x1) solid with lights
☐ 100 (4"x1.1/4) solid w/alum.rim W 32mm
☐ 125 (5"x1.1/4) black w/plastic rim
☐ 125 (5"x1.1/4) solid w/alum.rim W 32mm
☐ 125 (5"x1) solid with lights

Pockets

16



☐ Seat pocket

☐ Big under
the seat

☐ Backrest pocket

Lined cushion
in foam polyurethan

17



- ☐ 3 cm high
☐ 5 cm high

Eventual anodized set:
forks if anodized and
24" LIGHT rear wheels

18

- ☐ Black
☐ Silver
☐ Light blue
☐ Red
☐ Orange
☐ Cyclamen light purple

Options

- 19 ☐ "M+F" velcro strips (hook and loop) on seat
20 ☐ "M" velcro strips for VICAIR on seat
21 ☐ "M" velcro strips (hook) on seat
22 ☐ "F" velcro strips (loop) on seat
23 ☐ Abductor belt
24 ☐ Catheter holder pocket
25 ☐ Pump for pneumatics
26 ☐ 45° padded pelvic belt
27 ☐ Black backpack
28 ☐ Front protective covers

Notes