

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____

PAGE 1/2
 Oct.2026

Dealer _____

Key ☐ standard ☐ option

CHILDREN 3000

L Seat width

☐ 24* ☐ 26* ☐ 28* ☐ 30 ☐ 32 ☐ 34 ☐ 36

* Not available with detachable foot-rests

CHILDREN 3000 PLUS

L Seat width (adjustable within the selected range)

☐ RANGE 1 ☐ 24* ☐ 26* ☐ 28* ☐ 30 ☐ 32 ☐ 34

☐ RANGE 2 ☐ 26* ☐ 28* ☐ 30 ☐ 32 ☐ 34 ☐ 36 ☐ 38

* Not available with detachable foot-rests

P Seat depth

☐ 28 ☐ 30 (adjustable from 28 to 34 cm)

☐ 32 ☐ 34 ☐ 36 ☐ 38 (adjustable from 32 to 38 cm)

HS Back-rest height

☐ 25 ☐ 27 ☐ 29 ☐ 31 ☐ 33 ☐ 35 ☐ 37 ☐ 39 ☐ 41

HA Front seat/floor height

☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48

HP Rear seat/floor height

☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46

DP Foot-rest/seat distance


☐ 17* ☐ 19* ☐ 21* ☐ 23 ☐ 25 ☐ 27 ☐ 29 ☐ 31 ☐ 33

* Not available with detachable foot-rests (unless the front frame is abducted)

cm

CHILDREN 3000

Fold.alum.frame for Kid crash tested



Frame colour

☐ Deep sky light blue ☐ Orange ☐ Pearl red ☐ Metallic pink

☐ Basic colours

☐ Special colours
 (see "Colours" catalogue)

①

Backrest upholstery

☐ Tension adjustable nylon back-rest

☐ Tension adj. breathable b.

☐ Without back-rest

☐ Without seat

②

Upholstery colour

(only nylon)

☐ Black

☐ Light blue

☐ Red

☐ Green


☐ Cornflower blue ☐ Orange

③

Front frame

☐ Standard
 (Not suggested for L < 32cm)

☐ Abducted



④

Foot-rest


☐ Detachable

☐ Not detachable

☐ Elevating (detachable)
 Subject to technical approval

⑤


Separate aluminium profile footplates



☐ Foot-holding belts

⑥

Single aluminium profile footplate





☐ Foot-holding belts

⑦

Side-guards

☐ Plastic (standard)


(Not available with reclining backrest)



☐ Aluminium (curved, same colour)

⑧

Height adjustable "T" armrest



Elbow-rest


☐ Classic (23 cm long)

☐ Classic (30 cm long)


☐ Straight (26 cm long)

⑨


Standard push handles



Straight height adjustable push handles



Abducted height adjustable push handles



⑩

Rear wheels

☐ 24" Light (24 spokes) (*)

☐ 22" Light (24 spokes)

☐ 20" Light (20 spokes)

Tyres

☐ Pneum. high pressure smooth

☐ Pneum. high pressure plugged

☐ Solid Grey 20"/22"x1/24x1.3/8

☐ Solid black KIK 24" x1

☐ Marathon Plus 22" or 24"

(*) Only with P = 34 or 36 cm and HP from 40 to 46 cm

⑪

Handrims

☐ Aluminium silver

☐ Shiny inox

☐ Dark aluminium

☐ Rubber coated

☐ Carbolife CURVE L 24"

☐ Carbolife GEKKO 24"

☐ Silicon covers

⑫

Spoke guards

☐ Transparent

☐ Black

☐ Fantasy

code pattern:

(see catalogue of spokeguards)

⑬

Front castors

☐ 80 (3") solid with lights

☐ 80 (3") solid

☐ 100 (4") solid with lights

☐ 100 (4") solid

☐ 100 (4") solid w/silver alum.rim

☐ 100 (4") solid w/black alum.rim

☐ 125 (5") solid

☐ 125 (5") solid w/silver alum.rim

☐ 125 (5") solid w/black alum.rim

⑭

Rear wheel camber

☐ 1°

☐ 3°

☐ 5°

☐ Quick release axles Ø12mm with 42mm bushes (def.)

☐ Quick release axles Ø1/2" (12,7mm) with 42mm bushes

☐ Quick release axles Ø1/2" (12,7mm) with 50mm bushes

⑮


☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____

PAGE 2/2
 Oct.2026

CHILDREN 3000


Dealer

Liftable anti-tip device Children 3000
 for Children 3000 rear frame




- ☐ Right
☐ Left
☐ Right and left

Revolving anti-tip device
 for CURVED rear frame




- ☐ Right
☐ Left
☐ Right and left

☐ **Extended wheel plate**




Rear frame of Children 3000 (requires liftable anti-tip for Children 3000 and is compulsory with reclining back-rest)



Curved rear frame (requires revolving anti-tip device)

Single drive


(Only with 0° rear wheel camber)
 (Only with 22" or 24" rear wheels)



- ☐ Right
☐ Left

Head-rest

Backrest extension
 (cannot grow in width)




- ☐ Height adjustable
 Default: h = 20 cm d = 15 cm

Orbital


- ☐ Anatomic



- ☐ Moulded




- ☐ With lateral supports



Specify "h" & "d" where requested
 Standard: h = 15 cm d = 10 cm

h = ... cm d = ... cm


Abductor



- ☐ 1 ☐ 3 ☐ 5 ☐ 7 cm
☐ cm

Choose or specify the height from the seat cover

Lateral supports



- ☐ Right
☐ Left

Specify the height "h"

Table

- ☐ Transparent
☐ Black soft padded

specify the fastening:
☐ single ☐ double



(Only with arm-rest)

Anodized set

(Rear wheel rim&hub and front fork)

- ☐ Black
☐ Silver
☐ Light blue
☐ Red
☐ Orange
☐ Cyclamen light purple

Brakes

- ☐ Classic - short handle
☐ Classic - pulling brake
☐ Classic - long handle

- ☐ Right single brake
☐ Left single brake

- ☐ Drum
☐ Drum + classic
 (needed drum wheels)

Lined cushion

in foam polyurethane



- ☐ 3 cm high
☐ 5 cm high


Options

- 28 ☐ "M" & "F" velcro on seat
 29 ☐ "M" velcro strip (hook) on seat
 30 ☐ "F" velcro strip (loop) on seat
 31 ☐ Label of connection seat/back-rest
 32 ☐ Padded rigid cushion
 33 ☐ Chest harness
 34 ☐ Chest belt
 35 ☐ Abductor belt
 36 ☐ Pump for pneumatics
 37 ☐ 45° pelvic belt
 38 ☐ 45° padded pelvic belt
 39 ☐ Black backpack
 40 ☐ **Stickers for car fastenings** **crash tested**

☐ **Reclining backrest**

from 90° to 150°
 (compulsory:

- rear frame for Children 3000
- straight or abducted height adjustable push handles
- stretched bar
- extended wheel plate



Notes