

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____ User's weight Kg _____

Dealer

Key


☐ standard

☐ option

FUNKY

- L** Seat width ☐34 ☐36 ☐38 ☐40 ☐42 ☐44 ☐46
P Seat depth ☐36 ☐38 ☒40 ☐42 ☐44
HS Back-rest height ☐27 ☐29 ☐31 ☐33 ☐35 ☐37 ☐39
HA Front seat/floor height ☐46 ☒48 ☐50
HP Rear seat/floor height ☐40 ☐42 ☒44 ☐46 ☐48
DP Foot-rest/seat distance ☐35 ☐37 ☐39 ☐41 ☐43 ☐45 ☐47 ☐49

cm



Frame colour

- ☐ Metallic grey ☐ Pearl red ☐ Matte black ☐ Metallic blue
☐ (see "Colours" catalogue) from:
☐ Neon green ☐ Neon orange

①


Upholstery

- ☐ Tension adjustable nylon backrest
☐ Tension adj. breathable backrest
☐ Without any back-rest upholstery
☐ ~~AMEKH~~ (separate order form)
☐ other seating system:


②

Single tubular foot-plate


③



☐ Open




☐ Closed in aluminium




☐ Carbon closed footplate

Pushing handles


④



☐ Without pushing handles




☐ Standard pushing handles



☐ Height adjustable pushing handles

Distance between the hub centre and back-rest


⑤



- ☐ 7 cm
☒ 9 cm (standard)
☐ 11 cm
☐ 13 cm
☐ 15 cm

Back-rest tilt


⑥



- ☐ 89° ☒ 91° (standard) ☐ 93°

Front frame angle

⑦




- ☐ 90° (standard) (A = 16 cm)
☐ 95° (A = 22 cm)

(approximate standard measures)

Clothes guards


⑧

- ☐ Plastic
☐ Aluminium straight
☐ Carbon fibre



Detachable "L" tip-up arm-rest h 22 cm

⑨




Elbow-rest

- ☐ Tubular
☐ Straight (26 cm long)
☐ Classic (23 cm long)
☐ Classic (30 cm long)

Brakes


⑩

- ☐ Push brake (std)
☐ Pull brake
☐ Scissor brake



Revolving anti-tip device

⑪



- ☐ Right
☐ Left

Notes

☐ Order N° _____ ☐ Quotation
Date _____ Signature _____
Customer Ref. _____

PAGE 2/2
May 2024

FUNKY

Dealer

Rear wheels

12

- ☐ 24" Light(24 spokes)
☐ 24" Star 6x3 (18 spokes)
☐ 24" (36 spokes)
☐ 24" MTB with Citybike tyres
☐ 24" MTB with Offroad tyres
☐ 24" Ultralight(24 spokes)
☐ 24" Spinergy LX

Tyres

- ☐ Pneum. Marathon Plus
☐ Pneum. high press. smooth
☐ Pneum. high press. plugged
☐ Solid KIK 24" not available for Ultralight wheels

Handrims

13

- ☐ Aluminium silver
☐ Aluminium dark
☐ OVAAL alum. silver 24"
☐ OVAAL alum. dark 24"
☐ OVAALG(dark w/silicon ins.)24"
☐ Chromium-plated steel 24"
☐ Rubber coated
☐ Titanium
☐ Carbolife CURVE L24"
☐ CL Curve L TetraGrip 24"
☐ CL Quadro rubber ctd 24"
☐ Carbolife GEKKO 24"
☐ Silicon covers

OVAAL/S/G only in size 24"
Only OVAAL/GEKKO available for Star6x3
OVAAL & Carbolife not available for MTB

Spoke guards

14

- ☐ Transparent 24"

- ☐ Black 24"

- ☐ Fantasy

code pattern
(see catalogue of spokeguards)

Rear wheel camber

17

- ☐ 1° ☐ 2°
☐ 3° ☐ 4°

- ☐ Quick release axles Ø12mm with 42 cambered bushes (default)
☐ Q.r.a. Ø 1/2" (12,7 mm) w/42 camb.bushes
☐ Q.r.a. Ø 1/2" (12,7 mm) w/50 camb.bushes


Front castors

16


- ☐ 80 (3") solid
☐ 80 (3") solid with lights
☐ 100 (4"x1.1/4) solid with plastic rim
☐ 100 (4"x1) solid with lights
☐ 100 (4"x1.1/4) solid alum. rim
☐ 125 (5"x1.1/4) solid with plastic rim
☐ 125 (5"x1.1/4) solid alum. rim

Pockets

18



- ☐ Small under the seat




- ☐ Big under the seat



- ☐ On the back-rest

Lined cushion
in foam polyurethane

19



- ☐ 3 cm high
☐ 5 cm high

Anodized set

(Rear wheel rims, rear wheel hubs, forks)

20

- ☐ Black
☐ Silver
☐ Light blue
☐ Red
☐ Orange
☐ Cyclamen light purple

(only 24" Light wheels)

Options

- 21 ☐ "M+F" velcro strips (hook and loop) on seat
22 ☐ "M" velcro strips (hook) on seat
23 ☐ "F" velcro strips (loop) on seat
24 ☐ Padded rigid cushion h 5
25 ☐ Chest harness
26 ☐ Chest belt
27 ☐ Abductor belt
28 ☐ Catheter holder pocket
29 ☐ Rx ☐ Lx Crutch holder
30 ☐ Pump for pneumatics
31 ☐ Suitcase holder
32 ☐ 45° pelvic belt
33 ☐ 45° padded pelvic belt
34 ☐ Black backpack
35 ☐ Front protective covers

Notes