


Key standard ☐ option

ZODIAC

Aluminium rigid frame



L Seat width	☐ 34 ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44
P Seat depth	☐ 36 ☐ 38 ☒ 40 ☐ 42 ☐ 44
HS Back-rest height	☐ 27 ☐ 29 ☐ 31 ☐ 33 ☐ 35 ☐ 37 ☐ 39
HA Front seat/floor height	☐ 46 ☐ 48 ☐ 50
HP Rear seat/floor height	☐ 40 ☐ 42 ☐ 44 ☐ 46
DP Foot-rest/seat distance	☐ 35 ☐ 37 ☐ 39 ☐ 41 ☐ 43 ☐ 45 ☐ 47 ☐ 49

cm

Frame colour

- ☐ Metallic grey ☐ Pearl white
☐ Glossy black ☐ Matte black
☐ (see "Colours" catalogue) from:
☐ Neon green ☐ Neon orange

①


Upholstery

- ☐ Tension adjustable nylon back-rest
☐ Tension adj. breathable back-rest
☐ Without any back-rest upholstery
☒ **SAMEKH**
☐ other seating system:.....

②


Single tubular foot-plate

③



(not height adjustable)


☐ Open




☐ Closed

Push handles


④



☐ Without pushing handles




☐ Standard pushing handles



☐ Height adjustable pushing handles

Distance between the hub centre and back-rest


⑤



- ☐ 7 cm ☐ 13 cm
☒ 9 cm (standard) ☐ 15 cm
☐ 11 cm

Back-rest tilt


⑥



- ☐ 89°
☒ 91° (standard)
☐ 93°

Front angle

⑦



- ☐ 83° (standard) (A=16cm)
☐ 75° (A=23cm)

Rear wheel camber

⑧

- ☐ 1° ☐ 2° ☐ 3° ☐ 4°


- ☐ Quick release axles Ø12mm with 42 cambered bushes (default)
☐ Q.r.a. Ø 1/2" (12,7 mm) w/42 camb.bushes
☐ Q.r.a. Ø 1/2" (12,7 mm) w/50 camb.bushes

Side-guards

⑨



☐ Aluminium straight



☐ Aluminium with carbon fibre wing

☐ Aluminium with carbon fibre wing - detachable

"L" arm-rest, detachable, h 22cm

⑩

☐ Tip-up



Elbow-rest

- ☐ Tubular
☐ Straight (26 cm long)
☐ Classic (23 cm long)
☐ Classic (30 cm long)

Brakes

⑪


☐ Standard



☐ Scissor




☐ Composite



Revolving anti-tip device

⑫



☐ Right

☐ Left

Notes

☐ Order N° _____ ☐ Quotation
Date _____ Signature _____
Customer Ref. _____

PAGE 2/2
May 2024

ZODIAC

Rear wheels**12**

- ☐ 24" Light(24 spokes)
☐ 24" Star 6x3 (18 spokes, black or silver)
☐ 24" (36 spokes)
☐ 24" Mountain-bike with Citybike tyre
☐ 24" Mountain-bike with Offroad tyre
☐ 24" Ultralight (24 spokes) black
☐ 24" Spinergy LX
spokes colour.....(black def.)
☐ 25" Spinergy LX
spokes colour.....(black def.)
☐ 25" Light (24 spokes) black

Tyres

- ☐ Pneum. Marathon Plus
☐ Pneum. high pressure Smooth
☐ Pneum. high pressure Plugged 24"
☐ Pn.Schwalbe Right Run 24"
☐ Solid black KIK 24"
not available for Ultralight wheels

Pushrims**13**

- ☐ Aluminium silver
☐ Aluminium dark
☐ OVAAL alum. silver 24"
☐ OVAAL alum. dark 24"
☐ OVAALG(dark w/silicon ins.)24"
☐ Chromium-plated steel 24"
☐ Rubber coated
☐ Titanium
☐ Carbolife CURVE L24"
☐ CL Curve L TetraGrip 24"
☐ CL Quadro rubber ctd 24"
☐ Carbolife GEKKO 24"
☐ Silicon covers

OVAAL/S/G only in size 24"
Only OVAAL/GEKKO available for Star6x3
OVAAL & Carbolife not available for MTB

Spokeguards**14**


- ☐ Transparent 24"
☐ Black 24"
☐ Fantasy
code pattern
(see catalogue of spokeguards)

**Lined cushion
in foam polyurethan****16**


- ☐ 3 cm high
☐ 5 cm high

Front castors**15**

- ☐ 80 (3") solid
☐ 80 (3") solid with lights
☐ 100 (4"x1.1/4) solid with plastic rim
☐ 100 (4"x1) solid with lights
☐ 100 (4"x1.1/4) solid with alum. rim
☐ 125 (5"x1.1/4) solid with plastic rim
☐ 125 (5"x1.1/4) solid with alum. rim
☐ 125 (5"x1) solid with lights

Pockets**17**


- ☐ Small under the seat



- ☐ Big under the seat



- ☐ On the back-rest

Backrest bar

- ☐ Standard:
P(bar) 7,5cm;
H(bar) 19cm;
dia 21mm.
☐ Custom:
P.....cm;
H.....cm.

Anodized set

(Rear wheel rims & hubs if LIGHT 24",
forks and fittings)

- ☐ Black
☐ Silver
☐ Light blue
☐ Red
☐ Orange
☐ Cyclamen light purple

(only 24" Light wheels)

Options

- 20 ☐ "M+F" velcro strips (hook and loop) on seat
21 ☐ "M" velcro strips (hook) on seat
22 ☐ "F" velcro strips (loop) on seat
23 ☐ Padded rigid seat cushion h 5
24 ☐ Chest harness
25 ☐ Chest belt
26 ☐ Abductor belt
27 ☐ Catheter holder pocket
28 ☐ Rx ☐ Lx Crutch holder
29 ☐ Pump for pneumatics
30 ☐ Suitcase holder
31 ☐ 45° pelvic belt
32 ☐ 45° padded pelvic belt
33 ☐ Black backpack
34 ☐ Front protective covers

Notes