

☐ Order N° _____ ☐ Quotation
Date _____ Signature _____
Customer Ref. _____ User's weight Kg _____


PAGE 1/2
May 2024

Dealer _____

Key ☐ standard ☐ option

- L** Seat width ☐34 ☐36 ☐38 ☐40 ☐42 ☐44
- P** Seat depth ☐36 ☐38 ☐40 ☐42 ☐44 ☐46
- HS** Backrest height ☐29 ☐31 ☐33 ☐35 ☐37 ☐39 ☐41
- HA** Front seat/floor height ☐46 ☐48 ☐50
- HP** Rear seat/floor height ☐40 ☐42 ☐44 ☐46
- DP** Footrest/seat distance ☐35 ☐37 ☐39 ☐41 ☐43 ☐45 ☐47 cm

HALLEY.TI TITANIUM fold.frame



Frame colour

- ☐ Shiny titanium ☐ Natural titanium
☐ Glossy black ☐ Matte black
☐ Colour
(see "Colours" catalogue)

①

Upholstery

- ☐ Tension adjustable nylon back-rest
☐ Tens.adjustable breathable backrest
☐ Without any backrest upholstery
☐ ~~AMEKH~~ custom postural backrest
(fill up the proper order form)

②



Upholstery colour

- ☐ Black
☐ Coloured.....
Available colours:
NYLON: orange, red, green, cornflower blue, light blue;
BREATHABLE: red, grey.



③

Front frame

- ☐ Y shaped (default) ☐ Straight




- ☐ Front frame angle (with DP=43 cm and P=40 cm) (default) ☐ Extended front frame angle (with DP=43 cm and P=40 cm)



(approx.measurements)

④

Automatic closure foot-plate in carbon fibre





- ☐ Foot-holding belts

⑤

Separate aluminium profile tilt adjustable foot-plates


- ☐ forward ☐ backward



- ☐ Foot-holding belts

⑥

Closed carbon fibre single foot-plate





- ☐ Foot-holding belts

⑦

Single aluminium profile tilt adjustable foot-plate

- ☐ forward ☐ backward




- ☐ up & outward

- ☐ Foot-holding belts


⑧

Side guards

- ☒ Carbon fibre



Tip-up, detachable sport arm-rest h. 22 cm in aluminium




- ☐ in titanium

Elbow-rest

- ☐ Tubular
☐ Classic (23 cm long)

⑩

Tip-up height adjustable from 22 to 28cm "L" arm-rest (detachable, open with sideguards)



Elbow-rest

- ☐ Tubular
☐ Straight (26 cm long)
☐ Classic (23 cm long)
☐ Classic (30 cm long)

⑪

Rear wheels

- ☐ 24" Ultralight (24 spokes) black
☐ 24" Light (24 spokes)
☐ 24" Star 6x3 (18 sps) silver or black
☐ 24" (36 spokes)
☐ 24" Mountain-bike CITYBIKE tyre
☐ 24" Mountain-bike OFFROAD tyre
☐ 24" Spinergy LX
spokes colour.....(black def.)

Tyres

- ☐ Pneumatic Marathon Plus
☐ Pneum. high press.plugged
☐ Pneum. high press. smooth
☐ Solid black KIK 24"
(not avail.with Ultralight)

⑫

Pushrims

- ☐ Aluminium silver
☐ Aluminium dark
☐ OVAAL 24" aluminium silver **
☐ OVAALS 24" aluminium dark **
☐ OVAALG 24" alum.dark w/silicon insert **
☐ Titanium
☐ Rubber coated
☐ Carbolife CURVE L 24"***
☐ CL Curve L TetraGrip 24"***
☐ CL Quadro rubber ctd 24"***
☐ Carbolife GEKKO 24"***
☐ Silicon covers
OVAAL/S/G only in size 24"
Only OVAAL/GEKKO available for Star6x3
** OVAAL & Carbolife not available for MTB

⑬

Order N° _____ ☐ Quotation
Date _____ Signature _____
☐ Customer _____

PAGE 2/2
May 2024

HALLEY.TI

Dealer

Spoke guards

14

- ☐ Transparent 24"
☐ Black 24"
☐ Fantasy
code pattern
(see catalogue of spokeguards)

Front castors

15

- ☐ 100 (4"x1.1/4) black solid plastic rim
☐ 100 (4"x1.1/4) solid with alum. rim
☐ 125 (5"x1.1/4) black solid plastic rim
☐ 125 (5"x1.1/4) solid with alum. rim
☐ 150 (6") pneumatic
☐ 150 (6") black solid

Revolving anti-tip device

16

- ☐ Right
☐ Left




Unbalancing Pedal

- ☐ Right
☐ Left



Back-rest pushing handles


17

- ☐ Standard ☐ Without
 
☐ Height adjustable ☐ Foldable pushing handles
 

Lined cushion in foam polyurethane

18

- ☐ 3 cm high
☐ 5 cm high



Anodized set

(Rear wheel rims & hubs only if 24" Light rear wheels, forks and fittings)

19

- ☐ Black
☐ Silver
☐ Light blue
☐ Red
☐ Orange
☐ Cyclamen light purple
(24" Light wheels only)

Brakes

20

- ☐ Classic - short handle
☐ Classic - pull to brake short handle
☐ Classic - long handle
☐ Scissor
☐ Composite

Opzioni

- 22 ☐ "M+F" velcro strip on seat
23 ☐ "M" velcro strip on seat
24 ☐ "M" velcro for Vicair
25 ☐ "F" velcro strip (loop) on seat
26 ☐ Seat/backrest connection flap
27 ☐ Abductor belt
28 ☐ Catheter holder pocket
29 ☐ Back-rest stretched bar

Opzioni

- 30 ☐ Rx ☐ Lx Crutch hold./single
31 ☐ Rx ☐ Lx Crutch hold./double
32 ☐ Pump for pneumatics
33 ☐ Suitcase holder
34 ☐ 45° padded pelvic belt
35 ☐ Black backpack
36 ☐ Front protective covers

Rear wheel camber

21

- ☐ 1° ☐ 2° ☐ 3° ☐ 4°
☐ Q.rel.axles Ø12mm w/42mm bushes (def.)
☐ Q.rel.a. Ø1/2" (12,7mm) w/42mm bushes
☐ Q.rel.a. Ø1/2" (12,7mm) w/50mm bushes

Notes