

☐ Order N° _____ ☐ Quotation
Date _____ Signature _____
Customer ref. _____ User's weight Kg _____

PAGE 1/2
May 2024

Dealer


Key ☐ standard ☐ option

L Seat width ☐ 34 ☐ 35 ☐ 36 ☐ 37 ☐ 38 ☐ 39 ☐ 40 ☐ 41 ☐ 42 ☐ 43 ☐ 44 ☐ 45 ☐ 46 ☐ 47 ☐ 48
P Seat depth ☐ 34 ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48 ☐ 50
HS Backrest height ☐ 27 ☐ 29 ☐ 31 ☐ 35 ☐ 37 ☐ 39
HA Front seat/floor H ☐ 40 ☐ 41 ☐ 42 ☐ 43 ☐ 44 ☐ 45 ☐ 46 ☐ 47 ☐ 48 ☐ 49 ☐ 50 ☐ 51 ☐ 52
HP Rear seta/floor H ☐ 37 ☐ 38 ☐ 39 ☐ 40 ☐ 41 ☐ 42 ☐ 43 ☐ 44 ☐ 45 ☐ 46 ☐ 47 ☐ 48 ☐ 49 ☐ 50
DP Footrest/seat dist. ☐ 35 ☐ 36 ☐ 37 ☐ 38 ☐ 39 ☐ 40 ☐ 41 ☐ 42 ☐ 43 ☐ 44 ☐ 45 ☐ 46 ☐ 47 cm

Colour:
- **frame** ☐ Glossy black ☐ Matte black
- **front forks** ☐ Met. blue ☐ Matte Blue Klein
- **footplate** ☐ Pearl white ☐ Matte purple
- **sideguard supp.** ☐ Tiffany ☐ Matte London smoke
☐ Champagne ☐ Colour.....
(excluded camouflage, neon & special colours)
☐ RAL/Changing col.....shiny or matte, from:
(see "Colours" catalogue)

Upholstery
☐ Tension adj. nylon backrest
☐ Tension adj. breathable backrest
☐ Upholst.colour: Black (default)
Seat
☐ Nylon black (default)
☐ Carbon fibre

THEMIS
PRESTIGE
Ergonomic Aluminium
custom rigid frame




Front frame

☐ **V shaped**
(allows only tilt adj. footplate)
☐ **Straight**

Footplate


☐ Tilt adjustable
Carbon fib.closed
(available with V shaped or straight front frame)
☐ Open
(available only with straight front frame)
☐ Carbon fibre closed
(available only with straight front frame)



Push handles


☐ Without any pushing handles
☐ Fixed pushing handles
☐ Foldable pushing handles
☐ Height adjustable pushing handles

COG Distance between the hub centre and back-rest




☐ 0 cm ☐ 8 cm
☐ 1 cm ☐ 9 cm
☐ 2 cm ☐ 10 cm
☐ 3 cm ☐ 11 cm
☐ 4 cm ☐ 12 cm
☐ 5 cm ☐ 13 cm
☐ 6 cm ☐ 14 cm
☐ 7 cm ☐ 15 cm

Backrest tilt



☐ 89°
☐ 91° (def.)
☐ 93°
☐ Foldable and tilt adjustable backrest (default)
☐ Fixed, welded backrest

Front frame angle



☐ 83° (default) A = 16 cm
☐ Custom A=cm
(approx measurements)

Carbon fiber side-guards with foldable backrest:

☐ fixed with CRB wing (def.)
☐ fixed with straight CRB
☐ detachable with CRB wing
☐ detachable with straight CRB




Carbon fiber side-guards with welded backrest:

☐ fixed with CRB wing (def.)
☐ fixed with straight CRB



Backrest bar


☐ Default:
P(bar depth) 7,5cm;
H(bar height) 19cm;
dia 21mm.
☐ Custom:
P.....cm;
H.....cm.



Brakes

☐ Composite brake
☐ Scissor brake
☐ Classic push brake

Revolving anti-tip



☐ Right
☐ Left

Rear wheel camber ☐ 1° ☐ 2° ☐ 3° ☐ 4°

☐ Quick release axles Ø 12mm with cambered bushes 42mm long (def.)
☐ Quick release axles Ø 1/2" (12,7mm) with cambered bushes 42mm long
☐ Quick release axles Ø 1/2" (12,7mm) with cambered bushes 50mm long

Dealer

Rear wheels

⑫

- ☐ 24" Ultralight (24 spokes)
- ☐ 24" Light (24 spks) ** see section 18 f/colours
- ☐ 24" Star 6x3 (18 spokes, silver or black)
- ☐ 24" (36 spokes)
- ☐ 24" Mountain-bike with Citybike tyre
- ☐ 24" Mountain-bike with Offroad tyre
- ☐ 24" Spinergy LX
spokes colour.....(black stand.)
- ☐ 25" Spinergy LX
spokes colour.....(black stand.)
- ☐ 25" Light (24 spokes) black

Tyres

- ☐ Pnum. Marathon Plus
- ☐ Pnum. high pressure Smooth
- ☐ Pnum. high pressure Plugged 24"
- ☐ Pn.Schwalbe Right Run 24"
- ☒ Solid KIK 24" not available for
black Ultralight wheels

Pushrims

⑬

- ☐ Aluminium silver 24"
- ☐ Aluminium dark 24"/25"
- ☐ OVAAL aluminium silver 24" *
- ☐ OVAALS aluminium dark 24" *
- ☐ OVAALG black w/sil.insert 24" *
- ☐ Rubber coated
- ☐ Titanium
- ☐ CarboLife CURVE L 24" *
- ☐ CL CurveL TetraGrip[®] rubb.c.24"*
- ☐ CL Quadro rubber ctd 24" *
- ☐ CarboLife GEKKO 24" *
- ☐ Silicon covers 24"

* not available for Star 6x3 or MTB wheels
*CL pushrims: only GEKKO suitable for Star 6x3 wheels

Spokes guards

14

- ☐ Transparent 24"
☐ Black 24"
☐ Fantasy
 code pattern
 (see catalogue of spokequards)


Front castors

15

- ☐ 80 (3'') solid black W 25 mm
- ☐ 80 (3'') solid with lights W 25 mm
- ☐ 100 (4"x1.1/4) black w/plastic rim
- ☐ 100 (4"x1) solid with lights
- ☐ 100 (4"x1.1/4) solid w/alum.rim W 32mm
- ☐ 125 (5"x1.1/4) black w/plastic rim
- ☐ 125 (5"x1.1/4) solid w/alum.rim W 32mm
- ☐ 125 (5"x1) solid with lights

Pockets

16



- ☐ Seat pocket
- ☐ Big under the seat
- ☐ Backrest pocket



Lined cushion
in foam polyurethan

17



- ☐ 3 cm high
- ☐ 5 cm high

Eventual anodized set:
forks and 24" LIGHT rear wheels

18

- ☐ Black
- ☐ Silver
- ☐ Light blue
- ☐ Red
- ☐ Orange
- ☐ Cyclamen light purple

Options

- 19 ☐ "M+F" velcro strips (hook and loop) on seat
- 20 ☐ "M" velcro strips for VICAIR on seat
- 21 ☐ "M" velcro strips (hook) on seat
- 22 ☐ "F" velcro strips (loop) on seat
- 23 ☐ Abductor belt
- 24 ☐ Catheter holder pocket
- 25 ☐ Pump for pneumatics
- 26 ☐ 45° padded pelvic belt
- 27 ☐ Black backpack
- 28 ☐ Front protective covers

Technical drawing for approval sent on: _____ approved on: _____

Note