


Dealer \_\_\_\_\_

**Key**       standard       option

- L** Seat width      34 36 38 40 42 44 46
- P** Seat depth      36 38 40 42 44
- HS** Back-rest height      27 29 31 33 35 37 39
- HA** Front seat/floor height      46 48 50
- HP** Rear seat/floor height      43 45 47
- DP** Foot-rest/seat distance      31 33 35 37 39 41 43 45 47 49 cm



- Frame colour**      ①
- Shiny titanium       Natural Titanium
  - Camouflage       Matte black
  - ..... (see "Colours" catalogue)
  - Neon green       Neon orange

- Upholstery**      ②
- Tension adjustable nylon backrest
  - Tension adj. breathable backrest
  - Without back-rest for seating system:
    - AMEKH** (separate order form)
    - other seating system: .....

- Upholstery colour**
- Black
  - Coloured (nylon).....  
colours: orange, red, green, cornflower, light blue
  - Coloured (breathable).....  
colours: red, dry green, grey

- Single tubular foot-plate**      ③
- 
- Open
  - Carbon fibre closed
  - Tilt adjustable Carbon fibre closed

- Pushing handles**      ④
- 
- Without pushing handles
  - Standard pushing handles
  - Height adjustable pushing handles
  - Foldable pushing handles

- Distance between the hub centre and back-rest**      ⑤
- 
- 4 cm
  - 5 cm
  - 7 cm (standard)
  - 9 cm
  - 11 cm

- Back-rest tilt**      ⑥
- 
- 89°
  - 91° (standard)
  - 93°

- Frog Legs (included)**      ⑦
- 
- User weight: ..... kg

- Side guards**      ⑧
- 
- Straight aluminium, same colour
  - Straight aluminum, matte black

- Brakes**      ⑨
- Push brake (stand.)
  - Pull brake
  - Scissor brake
- 

- Revolving anti-tip device**      ⑩
- 
- Right
  - Left

Notes \_\_\_\_\_

Dealer

<p><b>11 Rear wheels</b></p> <p><input type="checkbox"/> 24" (36 spokes)</p> <p><input type="checkbox"/> 24" Mountain-bike with Citybike tyre</p> <p><input type="checkbox"/> 24" Mountain-bike with Offroad tyre</p> <p><input type="radio"/> 24" Spinergy spokes colour.....(black stand.)</p> <p><input type="radio"/> 25" Spinergy spokes colour.....(black stand.)</p> <p><b>Tyres</b></p> <p><input type="checkbox"/> Pneum. Marathon Plus</p> <p><input type="checkbox"/> Pneum. high pressure Smooth</p> <p><input type="checkbox"/> Pneum. high pressure Plugged 24"</p> <p><input type="radio"/> Solid KIK 24"</p>	<p><b>12 Handrims</b></p> <p><input type="checkbox"/> Aluminium silver</p> <p><input type="checkbox"/> Aluminium dark</p> <p><input type="radio"/> Chromium-plated steel 24"</p> <p><input type="radio"/> Rubber coated</p> <p><input type="radio"/> Titanium</p> <p><input type="radio"/> CarboLife CURVE L 24" *</p> <p><input type="radio"/> CL Curve L TetraGrip 24" *</p> <p><input type="radio"/> CL Quadro rubber ctd 24" *</p> <p><input type="radio"/> CarboLife GEKKO 24" *</p> <p><input type="radio"/> Silicon covers 24"</p> <p>* not available with Mountain bike wheels</p>	<p><b>13 Rear wheel camber</b></p> <p><input type="checkbox"/> 1° <input type="checkbox"/> 2° <input type="checkbox"/> 3° <input type="checkbox"/> 4°</p> <p><input type="checkbox"/> Quick rel.axles <math>\varnothing</math> 12 mm (default)</p> <p><input type="checkbox"/> Quick rel. axles <math>\varnothing</math> 1/2" (12,7 mm)</p> <p><input type="checkbox"/> Cambered bushes length 42 mm (def.)</p> <p><input type="checkbox"/> Cambered bushes length 50 mm</p>	<p><b>14 Front castors</b></p> <p><input type="checkbox"/> 80 (3") solid</p> <p><input type="radio"/> 80 (3") solid with lights</p> <p><input type="checkbox"/> 100 (4"x1.1/4) solid with plastic rim</p> <p><input type="radio"/> 100 (4"x1) solid with lights</p> <p><input type="checkbox"/> 100 (4"x1.1/4) solid with alum. rim</p> <p><input type="checkbox"/> 125 (5"x1.1/4) solid with plastic rim</p> <p><input type="checkbox"/> 125 (5"x1.1/4) solid with alum. rim</p> <p><input type="radio"/> 125 (5"x1) solid with lights</p>	
		<p><b>15 Backrest bar</b></p> <p><input type="checkbox"/> Standard: P(bar) 7,5cm; H(bar) 19cm; dia 21mm.</p> <p><input type="checkbox"/> Custom: P.....cm; H.....cm.</p>		

<p><b>16 Pockets</b></p> <p><input type="checkbox"/> Small under the seat</p> <p><input type="radio"/> Big under the seat</p> <p><input type="checkbox"/> On the back-rest</p>	<p><b>17 Lined cushion in foam polyurethane</b></p> <p><input type="radio"/> 3 cm high</p> <p><input type="radio"/> 5 cm high</p>	<p><b>Options</b></p> <p>18 <input type="checkbox"/> "M+F" velcro strips (hook and loop) on seat</p> <p>19 <input type="checkbox"/> "M" velcro strips (hook) on seat</p> <p>20 <input type="checkbox"/> "F" velcro strips (loop) on seat</p> <p>21 <input type="radio"/> Catheter holder pocket</p> <p>22 <input type="radio"/> Pump for pneumatics</p> <p>23 <input type="radio"/> 45° pelvic belt</p> <p>24 <input type="radio"/> 45° padded pelvic belt</p> <p>25 <input type="radio"/> Black backpack</p> <p>26 <input type="radio"/> Umbrella holder</p> <p>27 <input type="radio"/> Front protective covers</p>
--	---	---

**Notes**

-----

-----

-----

-----

-----

-----

-----

-----

-----

-----

-----

-----

-----

-----