

Dealer

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____ User's weight Kg _____


PAG. 1/2
 May 2024

Key☐ standard☐ option

- L** Seat width ☐ 34 ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46
- P** Seat depth ☐ 34 ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46
- HS** Back-rest height ☐ 29 ☐ 31 ☐ 33 ☐ 35 ☐ 37 ☐ 39 ☐ 41
- HA** Front seat/floor height ☐ 42 ☐ 43 ☐ 44 ☐ 46 ☐ 48 ☐ 50
- HP** Rear seat/floor height ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46
- DP** Foot-rest/seat distance ☐ 33 ☐ 35 ☐ 37 ☐ 39 ☐ 41 ☐ 43 ☐ 45 ☐ 47 cm

ALHENA

Titanium folding frame

**Frame colour**

- ☐ Shiny titanium ☐ Matte black
☐ Pearl white ☐ Natural titanium
☐ Neon green ☐ Neon orange
☐from:
 (see "Colours" catalogue)


Backrest upholstery

- ☐ Tension adjustable nylon back-rest
☐ Tens.adjustable breathable backrest
☐ Without any backrest upholstery
☐ **FAMEKH** custom postural backrest
 (fill up the proper order form)


Upholstery colour

- ☐ Black
☐ Coloured.....
 Available colours:
 NYLON: orange, red, green, cornflower blue, light blue;
 BREATHABLE: red, grey.


☐ **Not detachable foot-rest** ④
 (with DP=43 cm and P=40 cm)
 (standard)



☐ **Not detachable foot-rest**
 (with DP=43 cm and P=40 cm)

☐ **Detachable foot-rest**

- ☐ 84° (std)
☐ 76°




(approx.measurements)




Front Frame ⑤☐ STANDARD

☐ Y
 (only in case of
 not detachable foot-rest)


☐ V
 (only in case of not detachable foot-rest
 & without any automatic closure footplate)

**Separate aluminium profile tilt adjustable footplates** ⑧

☐ forward ☐ backward

☐ Foot-holding belts☐ **Closed carbon single footplate** ⑥☐ Foot-holding belts☐ **Automatic closure footplate in carbon** ⑦☐ Foot-holding belts**Single aluminium profile tilt adjustable footplate** ⑨

☐ forward ☐ backward

☐ up & outward☐ Foot-holding belts**Titanium tube structure** ⑩

- ☐ open
☐ closed by carbon plate
☐ Forward ☐ Backward
 (default)

Side-guards ⑪

☐ Aluminium with carbon fibre wing



☐ Plastic




☐ Curved aluminium
 (same colour or black)



☐ Straight aluminium
 (same colour or black)



☐ **Tip-up, detachable sport armrest h. 22 cm** ⑪
 (open with sideguards)

**Elbow-rest**

- ☐ Tubular
☐ Classic (23 cm long)

☐ **Tip-up height adjustable from 22 to 28 cm "L" armrest** ⑫
 (detachable, open with sideguards)

☐ Not detachable**Elbow-rest**

- ☐ Tubular
☐ Straight (26 cm long)
☐ Classic (23 cm long)
☐ Classic (30 cm long)

Rear wheel ⑬

- ☐ 24" Ultralight (24 spokes) black
☐ 24" Light (24 spokes)
☐ 24" Star 6x3 (18 sps) silver or black
☐ 24" (36 spokes)
☐ 24" Mountain-bike CITYBIKE tyre
☐ 24" Mountain-bike OFFROAD tyre
☐ 24" Spinergy LX
 Spokes colour.....(black def.)
☐ 25" Spinergy LX
 Spokes colour.....(black def.)
☐ 25" Light (24 spokes) black

Tyre

- ☐ Pneumatic Marathon Plus
☐ Pneum. high press.plugged 24"
☐ Pneum. high press. smooth
☐ Solid black KIK 24"
 (not available with Ultralight)

Pushrim ⑭

- ☐ Aluminium silver
☐ Aluminium dark
☐ OVAAL alum. silver 24"
☐ OVAAL alum. dark 24"
☐ OVAALG(dark w/silicon ins.)24"
☐ Chromium-plated steel 24"
☐ Rubber coated
☐ Titanium
☐ Carbolife CURVE L24"
☐ CL Curve L TetraGrip 24"
☐ CL Quadro rubber ctd 24"
☐ Carbolife GEKKO 24"
☐ Silicon covers
 OVAAL/S/G only in size 24"
 Only OVAAL/GEKKO available for Star6x3
 OVAAL & Carbolife not available for MTB

Dealer

Spokeguards

15

- ☐ Transparent 24"
- ☐ Black 24"
- ☐ Fantasy
- code pattern
- (see catalogue of spokeguards)

Castor wheel

16

- ☐ 80 (3") solid
- ☐ 80 (3") solid with lights
- ☐ 100 (4"x1.1/4) black solid plastic rim
- ☐ 100 (4"x1) solid with lights
- ☐ 100 (4"x1.1/4) solid with alum. rim
- ☐ 125 (5"x1.1/4) black solid plastic rim
- ☐ 125 (5"x1.1/4) solid with alum. rim
- ☐ 150 (6") pneumatic
- ☐ 150 (6") black solid

Revolving anti-tip device

17

- ☐ Right
- ☐ Left




Unbalancing Pedal

- ☐ Right
- ☐ Left



Extended wheel plate

18



(Not available with scissor brakes)

Push handles


19

- ☐ Standard
- ☐ Without
- ☐ Height adjustable
- ☐ Foldable push handles



Distance between hub centre and backrest

20



- ☐ 3,5 cm
- ☐ 6 cm (standard)
- ☐ 8,5 cm
- (with HP=44 cm)

Anodized set

(Rear wheel rims, rear wheel hubs, forks and fittings)

21


- ☐ Black
- ☐ Silver
- ☐ Light blue
- ☐ Red
- ☐ Orange
- ☐ Cyclamen light purple

(Light 24" wheels only)


Lined cushion in foam polyurethane

22


- ☐ 3 cm high
- ☐ 5 cm high



Pockets



☐ Small under the seat



☐ On the back-rest

Camber

25

- ☐ 1° ☐ 3°
- ☐ 2° ☐ 4°

- ☐ Quick rel. axles Ø 12 mm with 42mm cambered bushes (default)
- ☐ Quick rel. axles Ø 1/2" (12,7 mm) with 42mm cambered bushes
- ☐ Quick rel. axles Ø 1/2" (12,7 mm) with 50mm cambered bushes

Brakes

23

- ☐ Classic - short handle
- ☐ Classic - pulling brake short handle
- ☐ Classic - long handle
- ☐ Scissor
- ☐ Composite

Notes

Options

- 26 ☐ "M" velcro strip (hook) on seat
- 27 ☐ "F" velcro strip (loop) on seat
- 28 ☐ Seat/backrest connection flap
- 29 ☐ Padded rigid cushion h 5
- 30 ☐ Chest harness
- 31 ☐ Chest belt
- 32 ☐ Abductor belt
- 33 ☐ Catheter holder pocket
- 34 ☐ Backrest stretched bar
- 35 ☐ Rx ☐ Lx Crutch hold./single
- 36 ☐ Rx ☐ Lx Crutch hold./double
- 37 ☐ Pump for pneumatics
- 38 ☐ Suitcase holder
- 39 ☐ 45° pelvic belt
- 40 ☐ 45° padded pelvic belt
- 41 ☐ Front protective covers
- 42 ☐ Reinforced frame (max load 100-120Kg)

9,7 kg configuration

- Wheelchair 38x40x31 cm
- Nylon upholstery
- Titanium frame col.: shiny titanium
- Automatic closure footrest
- Fixed backrest pushing handles
- Aluminium sideguards with carbon fibre wings
- Scissor brakes
- Ultralight wheels
- Pneumatic high pressure smooth tyres
- Ultralight handrims
- 100 mm front castors
- Without pockets
- Max load: 100 kg