

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____ User's weight Kg _____

PAGE 1/2
 May 2024

Dealer _____

Key


☐ standard

☐ option


☐ **EOS Fix axle**

☐ **EOS Adj axle**

L Seat width ☐ 34 ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44
P Seat depth ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46
HS Back-rest height ☐ 27 ☐ 29 ☐ 31 ☐ 33 ☐ 35 ☐ 37 ☐ 39 ☐ 41
HA Front seat/floor height ☐ 46 ☐ 48 ☐ 50
HP Rear seat/floor height ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48
DP Foot-rest/seat distance ☐ 37 ☐ 39 ☐ 41 ☐ 43 ☐ 45 ☐ 47 ☐ 49 cm



EOS
 Titanium rigid frame



Frame colour

☐ Shiny titanium ☐ Glossy black
☐ Pearl white ☐ Matte black
☐ Natural titanium
☐ (see "Colours" catalogue) from:
☐ Neon green ☐ Neon orange

①

Backrest: Upholstery

☐ Tension adj nylon ☐ Fix (stand.)
☐ Breathable tension adj ☐ Without seat
☐ Without any backrest uphol. ☐ Rigid in carbon fibre
☐ **FAMEKH** (separate order form)
☐ other seating system

②


Backrest upholstery colour

☐ Black
☐ Coloured (nylon).....
 colours: orange, red, green, cornflower, light blue
☐ Coloured (breath).....
 colours: red, grey)

③


Tubular single titanium footplate

④



☐ Open ☐ Carbon fibre closed ☐ Tilt adjustable Carbon fib.closed

Push handles




☐ Without pushing handles ☐ Fixed pushing handles ☐ Height adjustable pushing handles ☐ Foldable pushing handles

⑤

Distance between the hub centre and backrest


⑥



☐ 9 cm (standard) ☐ 7 cm
☐ 0 cm (biamputated) ☐ 11 cm
☐ 3 cm ☐ 12 cm
☐ 5 cm

Backrest tilt


⑦




☐ 89°
☐ 91° (standard)
☐ 93°

Front frame angle

⑧



☐ 83° A = 14 cm (default)



☐ 75° A = 19 cm

Rear wheel camber

⑨

☐ 1° ☐ 2° ☐ 3° ☐ 4°
☐ Quick release axles Ø12mm with 42 cambered bushes (default)
☐ Q.r.a. Ø 1/2" (12,7 mm) w/42 camb.bushes
☐ Q.r.a. Ø 1/2" (12,7 mm) w/50 camb.bushes

Revolving anti-tip device

⑩



☐ Right
☐ Left

Side-guards

⑪

☐ Aluminium with carbon fibre wing ☐ Detachable with carbon f.wing

☐ Alumin. with straight carbon fibre ☐ Detachable w/straight crbn.f.

☐ Height adjustable from 22 to 28cm "L" arm-rest, detachable

☐ Tip-up

Elbow-rest

☐ Tubular
☐ Straight (26 cm long)
☐ Classic (23 cm long)
☐ Classic (30 cm long)

Brakes

⑬

☐ Scissor brake
☐ Composite brake

☐ Standard
☐ Pull brake


Pockets

⑭

☐ Seat pocket ☐ Big under the seat ☐ Back-rest pocket

Backrest bar

☐ Standard: dia 19mm
 P (bar) 7,5cm;
 H (bar) 19cm.
☐ Custom:
 P (bar).....cm;
 H (bar).....cm



Notes









☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____

PAGE 2/2
 May 2024

EOS

Dealer

Rear wheels (16) ☐ 24" Ultralight (24 spokes) ☐ 24" Light (24 spokes) ☐ 24" Star 6x3 (18 spokes, silver or black) ☐ 24" (36 spokes) ☐ 24" Mountain-bike with Citybike tyre ☐ 24" Mountain-bike with Offroad tyre ☐ 24" Spinergy LX spokes colour.....(black def.) ☐ 25" Spinergy LX spokes colour.....(black def.) ☐ 25" Light (24 spokes) black Tyres ☐ Pneum. Marathon Plus ☐ Pneum. high pressure Smooth ☐ Pneum. high pressure Plugged 24" ☐ Pn.Schwalbe Right Run 24" ☐ Solid black KIK 24" not available for Ultralight wheels	Pushrims (17) ☐ Aluminium silver ☐ Aluminium dark ☐ OVAAL alum. silver 24" ☐ OVAAL alum. dark 24" ☐ OVAALG(dark w/silicon ins.)24" ☐ Chromium-plated steel 24" ☐ Rubber coated ☐ Titanium ☐ Carbolife CURVE L24" ☐ CL Curve L TetraGrip 24" ☐ CL Quadro rubber ctd 24" ☐ Carbolife GEKKO 24" ☐ Silicon covers OVAAL/S/G only in size 24" Only OVAAL/GEKKO available for Star6x3 OVAAL & Carbolife not available for MTB	Spokeguards (18) ☐ Transparent 24" ☐ Black 24" ☐ Fantasy code pattern (see catalogue of spokeguards) Anodized set (20) (Rear wheel hub & rims if LIGHT 24", sideguard supports, forks and fittings) (only 24" Light wheels) ☐ Black ☐ Silver ☐ Light blue ☐ Red ☐ Orange ☐ Cyclamen light purple	Castor wheels (19) ☐ 80 (3") solid ☐ 80 (3") solid with lights ☐ 100 (4"x1.1/4) solid with plastic rim ☐ 100 (4"x1) solid with lights ☐ 100 (4"x1.1/4) solid with alum. rim ☐ 125 (5"x1.1/4) solid with plastic rim ☐ 125 (5"x1.1/4) solid with alum. rim ☐ 125 (5"x1) solid with lights Lined cushion (21) in foam polyurethan  ☐ 3 cm high ☐ 5 cm high
---	--	---	--

 STANDARD SEAT (def.) ☐  ERGO SEAT ☐ (carbon fibre rigid seat only suitable for standard seat)	Options 26 ☐ "M+F" velcro strips (hook and loop) on seat 27 ☐ "M" velcro strips (hook) on seat 28 ☐ "F" velcro strips (loop) on seat 29 ☐ Padded rigid cushion h 5 30 ☐ Chest harness 31 ☐ Chest belt 32 ☐ Abductor belt 33 ☐ Catheter holder pocket 34 ☐ Rx ☐ Lx Crutch holder 35 ☐ Pump for pneumatics 36 ☐ Suitcase holder 37 ☐ 45° pelvic belt 38 ☐ 45° padded pelvic belt 39 ☐ Black backpack 41 ☐ Front protective covers	7,5 kg configuration (FIX) 7,8 kg configuration (ADJ) - Wheelchair 38x40x31 cm - Nylon upholstery - Titanium frame colour: Shiny Titanium - Open foot-rest - Without pushing handles - Aluminium side guards with carbon fibre wing - Scissor brakes - Ultralight wheels - Smooth high press. pneum. tyres - Ultralight pushrim - 80 mm front castors - Without pockets - Max load: 100 kg
☐ TAPERED SEAT (L-12cm) (23)  (Tapered seat is possible only if the seat width is not lower than 36cm) ☐ Custom front opening (i.e. custom TAPERED seat)cm Tapered seat requires straight front frame: 	☐ STRAIGHT SEAT (def.) (24)  ☐ Y SHAPED FRONT FRAME (default)  (Tapered seat not possible) default L -12cm ☐ Custom Y/V front openingcm ☐ Reinforced frame (max load 100-120 Kg)	☐ V SHAPED FRONT FRAME  (Tapered seat not possible) default L -12cm (The tilt adjustable footplate only is suitable with the V shaped front frame) ☐ STRAIGHT FRONT FRAME (25)  (Compulsory if seat width is 34cm or if the seat is tapered)

Notes