

☐ Order N° _____ ☐ Quotation
Date _____ Signature _____

Customer Ref. _____ User's weight Kg _____

PAGE 1/2
May 2024

Dealer _____

Key

☐ standard

☐ option


☐ **IDRA2.0 FIX** DETACH. FOOT.
COD. IDR2FE

☐ **IDRA 2.0 ADJ** DETACH. FOOT.
COD. IDR2RE

L Seat width ☐ 34 ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44
P Seat depth ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46
HS Back-rest height ☐ 27 ☐ 29 ☐ 31 ☐ 33 ☐ 35 ☐ 37 ☐ 39 ☐ 41
HA Front seat/floor height ☐ 46 ☐ 48 ☐ 50
HP Rear seat/floor height ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48
DP Foot-rest/seat distance ☐ 37 ☐ 39 ☐ 41 ☐ 43 ☐ 45 ☐ 47 ☐ 49 cm

IDRA2.0

with DETACHABLE footrests
Aluminium rigid frame



Frame colour

- ☐ Glossy black ☐ Matte black
☐ Pearl white ☐ Pearl red
☐ Glossy anthracite grey
☐ (see "Colours" catalogue) from:
☐ Neon green ☐ Neon orange

Upholstery

Backrest:

- ☐ Tension adj nylon
☐ Breathable tension adj
☐ Without any backrest upholstery
☐ **FAMEKH** (separate order form)
☐ other seating system


Seat:

- ☐ Fix (standard, screws)
☐ Tension adj. (velcro)


Upholstery colour

- ☐ Black
☐ Coloured (nylon).....
colours: orange, red, green, cornflower, light blue
☐ Coloured (breathable).....
colours: red, grey)


Detachable separate footrests



Push handles




Distance between the hub centre and backrest




- ☐ **9 cm (standard)** ☐ 7 cm
☐ 0 cm (biamputated) ☐ 11 cm
☐ 3 cm ☐ 13 cm
☐ 5 cm

Backrest tilt



- ☐ 89°
☐ **91° (standard)**
☐ 93°

Separate aluminium profile tilt adjustable foot-plates:




☐ Foot-holding belts

Rear wheel camber

- ☐ 1° ☐ 2° ☐ 3° ☐ 4°

- ☐ **Quick rel. axles Ø 12mm with cambered bushes length 42 mm (default)**
☐ Quick rel. axles Ø 1/2" (12,7mm) w/42mm bushes
☐ Quick rel. axles Ø 1/2" (12,7mm) w/50mm bushes

Side-guards



☐ Aluminium with carbon fibre wing

☐ Aluminium with straight carbon fibre

Height adjustable from 22 to 28cm "L" arm-rest, detachable and tip-up



Elbow-rest


- ☐ Tubular
☐ Straight (26 cm long)
☐ Classic (23 cm long)
☐ Classic (30 cm long)

Brakes



- ☐ Standard
☐ Pull brake

Anti-tip device



- ☐ Right
☐ Left

Notes _____

☐ Order N° _____ ☐ Quotation
Date _____ Signature _____
Customer Ref. _____

PAGE 2/2
May 2024

IDRA2.0
with DETACHABLE footrests

Dealer

Rear wheels

14

- ☐ 24" Ultralight (24 spokes)
☐ 24" Light (24 spokes)
☐ 24" Star 6x3 (18 spokes, silver or black)
☐ 24" (36 spokes)
☐ 24" Mountain-bike with Citybike tyre
☐ 24" Mountain-bike with Offroad tyre
☐ 24" Spinergy LX
spokes colour.....(black def.)
☐ 25" Spinergy LX
spokes colour.....(black def.)
☐ 25" Light (24 spokes) black

Tyres

- ☐ Pneum. Marathon Plus
☐ Pneum. high pressure Smooth
☐ Pneum. high pressure Plugged 24"
☐ Pn.Schwalbe Right Run 24"
☐ Solid black KIK 24"
not available for Ultralight wheels

Pushrims

15

- ☐ Aluminium silver
☐ Aluminium dark
☐ Chromium-plated steel 24"
☐ Rubber coated
☐ Titanium
☐ OVAAL (silver aluminium)
☐ OVAALS (dark aluminium)
☐ OVAALG (w/sil.insert)
☐ CarboLife CURVE L 24" *
☐ CL CurveL TetraGrip rubb.c.24" *
☐ CL Quadro rubber ctd 24" *
☐ CarboLife GEKKO 24" *
☐ Silicon covers 24"

* CL pushrims: not suitable for MTB wheels
* CL pushrims: only GEKKO suitable for Star 6x3 wheels


Spokes guards

16

- ☐ Transparent 24"
☐ Black 24"
☐ Fantasy
code pattern
(see catalogue of spokeguards)

Lined cushion in foam polyurethan

18



- ☐ 3 cm high
☐ 5 cm high


Front castors

17

- ☐ 80 (3") solid
☐ 80 (3") solid with lights
☐ 100 (4"x1.1/4) solid with plastic rim
☐ 100 (4"x1) solid with lights
☐ 100 (4"x1.1/4) solid with alum. rim
☐ 125 (5"x1.1/4) solid with plastic rim
☐ 125 (5"x1.1/4) solid with alum. rim
☐ 125 (5"x1) solid with lights
☐ 150 (6") pneumatic
☐ 150 (6") solid

Pockets

19



- ☐ Small under the fix seat
or
☐ Small hanging from the seat

☐ Big under the seat

☐ On the back-rest

Anodized set

(Rear wheel rims & hubs if LIGHT 24",
sideguard supports, forks and fittings)


20

- ☐ Black
☐ Silver
☐ Light blue
☐ Red
☐ Orange
☐ Cyclamen light purple

(only 24" Light wheels)

Backrest bar

21



- ☐ Standard:
P(bar) 7,5cm;
H(bar) 19cm;
dia 21mm.
☐ Custom:
P.....cm;
H.....cm.

Options

- 23 ☐ "M+F" velcro strips (hook and loop) on seat
24 ☐ "M" velcro strips (hook) on seat
25 ☐ "F" velcro strips (loop) on seat
26 ☐ Padded rigid seat cushion h 5
27 ☐ Chest harness
28 ☐ Chest belt
29 ☐ Abductor belt
30 ☐ Catheter holder pocket
31 ☐ Rx ☐ Lx Crutch holder
32 ☐ Pump for pneumatics
33 ☐ Suitcase holder
34 ☐ 45° pelvic belt
35 ☐ 45° padded pelvic belt
36 ☐ Black backpack
37 ☐ Front protective covers

Notes