

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____ User weight Kg _____

PAGE 1/2
May 2024

Dealer _____

Key


☐ standard

☐ option

VENUS ELITE

Titanium rigid frame

L Seat width ☐34 ☐36 ☐38 ☐40 ☐42 ☐44 ☐46
P Seat depth ☐36 ☐38 ☒40 ☐42 ☐44
HS Back-rest height ☐27 ☐29 ☐31 ☐33 ☐35 ☐37 ☐39
HA Front seat/floor height ☐46 ☒48 ☐50
HP Rear seat/floor height ☐39 ☐41 ☒43 ☐45 ☐47
DP Foot-rest/seat distance ☐31 ☐33 ☐35 ☐37 ☐39 ☐41 ☐43 ☐45 ☐47 ☐49
 cm



Frame colour

- ☐ Shiny titanium ☐ Glossy black
☐ Pearl white ☐ Matte black
☐ Natural titanium
☐ (see "Colours" catalogue)
☐ Neon green ☐ Neon orange

①

Upholstery

- ☐ Tension adjustable nylon backrest
☐ Tension adj. breathable backrest
☐ Without any back-rest upholstery
☐ **FAMEKH** (separate order form)
☐ other seating system:


②

Upholstery colour

- ☐ Black
☐ Coloured (nylon).....
 colours: orange, red, green, cornflower, light blue
☐ Coloured (breathable).....
 colours: red, grey)


Single tubular titanium footplate

③




Push handles

④




Distance between the hub centre and back-rest

⑤




Backrest tilt

⑥




Frame lenght

⑦



Seat with shock absorbers


⑧



User's weight : kg
 It is not possible to have two anti-tip devices

Side-guards

⑨



"L" detachable arm-rest - h 22 cm

⑩

☐ Tip-up




Elbow rest

- ☐ Tubular
☐ Straight (26 cm long)
☐ Classic (23 cm long)
☐ Classic (30 cm long)


Brakes

⑪



Anti-tip device

⑫



- ☐ Right
☐ Left

Notes

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____ User weight Kg _____

PAGE 2/2
 May 2024

VENUS
 ELITE

Dealer _____

Rear wheels

⑫

- ☐ 24" Ultralight (24 spokes)
☐ 24" Light (24 spokes)
☐ 24" Star 6x3 (18 spokes, silver or black)
☐ 24" (36 spokes)
☐ 24" Mountain-bike with Citybike tyre
☐ 24" Mountain-bike with Offroad tyre
☐ 24" Spinergy LX
 spokes colour.....(black def.)
☐ 25" Spinergy LX
 spokes colour.....(black def.)
☐ 25" Light (24 spokes) black

Tyres

- ☐ Pneum. Marathon Plus
☐ Pneum. high pressure Smooth
☐ Pneum. high pressure Plugged 24"
☐ Pn.Schwalbe Right Run 24"
☐ Solid KIK 24"
 not available for Ultralight wheels

Pushrims

⑬

- ☐ Aluminium silver
☐ Aluminium dark
☐ Chromium-plated steel 24"
☐ Rubber coated
☐ Titanium
☐ OVAAL (silver aluminium)
☐ OVAALS (dark aluminium)
☐ OVAALG (w/sil.insert)
☐ CarboLife CURVE L 24" *
☐ CL Curvel TetraGrip rubb.c.24" *
☐ CL Quadro rubber ctd 24" *
☐ CarboLife GEKKO 24" *
☐ Silicon covers 24"
 OVAAL/S/G only in size 24"
 Only OVAAL/GEKKO available for Star6x3
 OVAAL & Carbolife not available for MTB

Spokeguards

⑭

- ☐ Transparent 24"
☐ Black 24"
☐ Fantasy
 code pattern
 (see catalogue of spokeguards)

Front castors

⑮

- ☐ 80 (3") solid
☐ 80 (3") solid with lights
☐ 100 (4"x1.1/4) solid with plastic rim
☐ 100 (4"x1) solid with lights
☐ 100 (4"x1.1/4) solid with alum. rim
☐ 125 (5"x1.1/4) solid with plastic rim
☐ 125 (5"x1.1/4) solid with alum. rim
☐ 125 (5"x1) solid with lights


Rear wheel camber

⑯

- ☐ 1° ☐ 2° ☐ 3° ☐ 4°
☐ Quick rel.axles Ø 12mm with
 cambered bushes length 42 mm (default)
☐ Quick rel. axles Ø 1/2" (12,7mm)
 w/42mm bushes
☐ Quick rel. axles Ø 1/2"
 (12,7mm) w/50mm bushes

Backrest bar


⑰




- ☐ Standard:
 P(bar) 7,5cm;
 H(bar) 19cm;
 dia 21mm.
☐ Custom:
 P.....cm;
 H.....cm.

Pockets


⑱



- ☐ Small under the seat



- ☐ Big under the seat



- ☐ On the back-rest

Lined cushion in foam polyurethane

⑲



- ☐ 3 cm high
☐ 5 cm high

Anodized set

(Rear wheel rims & hubs if LIGHT 24",
 sideguard supports, forks and fittings)

⑳

- ☐ Silver
☐ Black
☐ Light blue
☐ Red
☐ Orange
☐ Cyclamen light purple
 (only 24" Light wheels)

Options

- 21 ☐ "M+F" velcro strips (hook and loop) on seat
 22 ☐ "M" velcro strips (hook) on seat
 23 ☐ "F" velcro strips (loop) on seat
 24 ☐ Padded rigid seat cushion h 5
 25 ☐ Chest harness
 26 ☐ Chest belt
 27 ☐ Abductor belt
 28 ☐ Catheter holder pocket
 29 ☐ Rx ☐ Lx Crutch holder
 30 ☐ Pump for pneumatics
 31 ☐ Suitcase holder
 32 ☐ 45° pelvic belt
 33 ☐ 45° padded pelvic belt
 34 ☐ Black backpack
 35 ☐ Front protective covers

Notes