

☐ Order N° _____ ☐ Quotation
Date _____ Signature _____
Customer Ref. _____ User's weight Kg _____

PAGE 1/2
May 2024

Dealer

Key


☐ standard

☐ option

L Seat width ☐34 ☐36 ☐38 ☐40 ☐42 ☐44 ☐46
P Seat depth ☐36 ☐38 ☐40 ☐42 ☐44
HS Back-rest height ☐27 ☐29 ☐31 ☐33 ☐35 ☐37 ☐39
HA Front seat/floor height ☐46 ☐48 ☐50
HP Rear seat/floor height ☐39 ☐41 ☐43 ☐45 ☐47
DP Foot-rest/seat distance ☐31 ☐33 ☐35 ☐37 ☐39 ☐41 ☐43 ☐45 ☐47 ☐49 cm

VENUS

Titanium & alum.rigid frame



Frame colour

- ☐ Shiny titanium ☐ Glossy black
☐ Pearl white ☐ Matte black
☐ Natural titanium
☐ (see "Colours" catalogue)
☐ Neon green ☐ Neon orange

①

Upholstery

②

- ☐ Tension adjustable nylon back-rest
☐ Tension adj. breathable back-rest
☐ Without back-rest

Upholstery colour

- ☐ Black
☐ Coloured (nylon).....
colours: orange, red, green, cornflower, light blue
☐ Coloured (breathable).....
colours: red, grey)

Single tubular footplate

③

- ☐ Open ☐ Carbon fibre closed
☐ Aluminium closed
☐ Tilt adjustable Carbon fibre closed
☐ Footrest colour: black (default)
☐ Foot-rest colour.....


Push handles

④

- ☐ Without pushing handles
☐ Standard pushing handles
☐ Height adjustable pushing handles
☐ Foldable pushing handles

Distance between the hub centre and backrest


⑤



- ☐ 4 cm
☐ 5 cm
☐ 7 cm (standard)
☐ 9 cm
☐ 11 cm

Backrest tilt


⑥



- ☐ 89° ☐ 91° (standard) ☐ 93°

Seat with shock absorbers

⑦



User weight : kg

It is not possible to have two anti-tip devices


Side guards

⑧

- ☐ Aluminium with fixed carbon fibre wing ☐ Straight aluminium
☐ Aluminium with detachable carbon fibre wing

Detachable "L" arm-rest - h 22 cm

⑨



- ☐ Tip-up
Elbow-rest
☐ Tubular
☐ Straight (26 cm long)
☐ Classic (23 cm long)
☐ Classic (30 cm long)

Brakes

⑩


- ☐ Push brake (stand.)
☐ Pull brake

- ☐ Composite brake
☐ Scissor brake



anti-tip device

⑪



- ☐ Right
☐ Left

Notes

☐ Order N° _____ ☐ Quotation
Date _____ Signature _____
Customer Ref. _____

PAGE 2/2
May 2024

VENUS

Dealer

12

Rear wheels

- ☐ 24" Ultralight (24 spokes)
☐ 24" Light (24 spokes)
☐ 24" Star 6x3 (18 spokes, silver or black)
☐ 24" (36 spokes)
☐ 24" Mountain-bike with Citybike tyre
☐ 24" Mountain-bike with Offroad tyre
☐ 24" Spinergy LX
spokes colour.....(black def.)
☐ 25" Spinergy LX
spokes colour.....(black def.)
☐ 25" Light (24 spokes) black

Tyres

- ☐ Pneum. Marathon Plus
☐ Pneum. high pressure Smooth
☐ Pneum. high pressure Plugged 24"
☐ Pn.Schwalbe Right Run 24"
☐ Solid black KIK 24"
not available for Ultralight wheels

Pushrims

- ☐ Aluminium silver
☐ Aluminium dark
☐ Chromium-plated steel 24"
☐ Rubber coated
☐ Titanium
☐ OVAAL (silver aluminium)
☐ OVAALS (dark aluminium)
☐ OVAALG (w/sil.insert)
☐ CarboLife CURVE L 24" *
☐ CL CurveL TetraGrip rubb.c.24" *
☐ CL Quadro rubber ctd 24" *
☐ CarboLife GEKKO 24" *
☐ Silicon covers 24"
OVAAL/S/G only in size 24"
Only OVAAL/GEKKO available for Star6x3
OVAAL & Carbolife not available for MTB

13

Spokes guards

- ☐ Transparent 24"
☐ Black 24"
☐ Fantasy
code pattern
(see catalogue of spokeguards)

14

Front castors

- ☐ 80 (3") solid
☐ 80 (3") solid with lights
☐ 100 (4"x1.1/4) solid with plastic rim
☐ 100 (4"x1) solid with lights
☐ 100 (4"x1.1/4) solid with alum. rim
☐ 125 (5"x1.1/4) solid with plastic rim
☐ 125 (5"x1.1/4) solid with alum. rim
☐ 125 (5"x1) solid with lights

15


Rear wheel camber

- ☐ 1° ☐ 2°
☐ 3° ☐ 4°

16

- ☐ Quick rel. axles Ø 12mm with
cambered bushes length 42 mm (default)
☐ Quick rel. axles Ø 1/2" (12,7mm)
w/42mm bushes
☐ Quick rel. axles Ø 1/2"
(12,7mm) w/50mm bushes


Backrest bar




- ☐ Standard:
P(bar) 7,5cm;
H(bar) 19cm;
dia 19mm.
☐ Custom:
P.....cm;
H.....cm.

17


Pockets



- ☐ Small under the seat



- ☐ Big under the seat



- ☐ On the back-rest

18

Lined cushion

in foam polyurethane



- ☐ 3 cm high
☐ 5 cm high

19

Anodized set

(Rear wheel rims & hubs if LIGHT 24",
sideguard supports, forks and fittings)

- ☐ Black
☐ Silver
☐ Light blue
☐ Red
☐ Orange
☐ Cyclamen light purple
(only 24" Light wheels)

20

Options

- 21 ☐ "M" velcro strip (hook) on seat
22 ☐ "F" velcro strip (loop) on seat
23 ☐ Padded rigid cushion h 5
24 ☐ Chest harness
25 ☐ Chest belt
26 ☐ Abductor belt
27 ☐ Catheter holder pocket
28 ☐ Rx ☐ Lx Crutch holder
29 ☐ Pump for pneumatics
30 ☐ Suitcase holder
31 ☐ 45° pelvic belt
32 ☐ 45° padded pelvic belt
33 ☐ Black holdall
34 ☐ Front protective covers

Notes