

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____ User weight Kg _____

PAGE 1/2
 May 2024

Dealer

Key

☐ standard

☐ option

☐ **XLR8 FIX**


☐ **XLR8 ADJ**

- L** Seat width ☐ 34* ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44
P Seat depth ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46
HS Back-rest height ☐ 27 ☐ 29 ☐ 31 ☐ 33 ☐ 35 ☐ 37 ☐ 39 ☐ 41
HA Front seat/floor height ☐ 46 ☐ 48 ☐ 50
HP Rear seat/floor height ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48
DP Foot-rest/seat distance ☐ 37 ☐ 39 ☐ 41 ☐ 43 ☐ 45 ☐ 47 ☐ 49 cm

* With a seat width of 34, it is recommended a straight seat and a straight front frame (see section 24)

XLR8

Aluminium rigid frame



Frame colour

①

- ☐ Shiny alluminium ☐ Matte black
☐ (see "Colours" catalogue) from:
☐ Neon green ☐ Neon orange

Back-rest: Upholstery

- ☐ Tension adj nylon
☐ Breathable tension adj
☐ Without back-rest:
☐ **AMEKH** (separate order form)
☐ other seating system

Seat:

- ☐ Fix (stand.)
☐ Without seat
☐ Carbon fibre

Upholstery colour

③

- ☐ Black
☐ Coloured (nylon).....
 colours: orange, red, green, cornflower, light blue
☐ Coloured (breathable).....
 colours: red, grey)

Tubular single foot-rest (matte black)

④


- ☐ Open
☐ Carbon fibre closed
☐ Tilt adjustable Carbon fib.closed
 (OPEN or CLOSED only if the front frame is STRAIGHT)
 (suitable for any type of front frame)

Push handles

- ☐ Without pushing handles
☐ Fixed pushing handles
☐ Height adjustable pushing handles
☐ Foldable pushing handles

Distance between the hub centre and backrest


⑥



- ☐ 9 cm (standard)
☐ 0 cm (biamputated)
☐ 3 cm
☐ 5 cm
☐ 7 cm
☐ 11 cm
☐ 12 cm

Backrest tilt


⑦



- ☐ 89°
☐ 91° (standard)
☐ 93°

Front frame opening

⑧



- ☐ 83° = standard


Rear wheel camber

⑨

- ☐ 1° ☐ 2° ☐ 3° ☐ 4°
☐ Quick release axles Ø12mm with 42 cambered bushes (default)
☐ Q.r.a. Ø 1/2" (12,7 mm) w/42 camb.bushes
☐ Q.r.a. Ø 1/2" (12,7 mm) w/50 camb.bushes

Revolving anti-tip device

⑩




Side-guards

(if the backrest is foldable)

⑪




- ☐ Fix sideguard with carbon fibre wing
☐ Detachable with carbon f.wing



- ☐ Fix sideg. with straight carbon fibre
☐ Detachable w/straight crbn.f.

Height adjustable from 22 to 28cm "L" arm-rest, detachable

⑫



- ☐ Tip-up
☐ Fix (def.)

Elbow-rest

- ☐ Tubular
☐ Straight (26 cm long)
☐ Classic (23 cm long)
☐ Classic (30 cm long)

Brakes

⑬



- ☐ Scissor brake
☐ Composite brake



- ☐ Standard
☐ Pull brake

Backrest


⑭

- ☐ Foldable backrest (def.)
 (Note: detachable sideguards are needed)
☐ Welded backrest
 (Note: fix sideguards are needed)



Backrest bar

⑮



- ☐ Standard: dia 25mm
 P (bar) 7,5cm;
 H (bar) 19cm.
☐ Custom:
 P (bar).....cm;
 H (bar).....cm.

Notes

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____

PAGE 2/2
 May 2024

XLR8

Dealer

Rear wheels

16

- ☐ 24" Ultralight (24 spokes)
☐ 24" Light (24 spokes)
☐ 24" Star 6x3 (18 spokes, silver or black)
☐ 24" (36 spokes)
☐ 24" Mountain-bike with Citybike tyre
☐ 24" Mountain-bike with Offroad tyre
☐ 24" Spinerger LX
 spokes colour.....(black stand.)
☐ 25" Spinerger LX
 spokes colour.....(black stand.)
☐ 25" Light (24 spokes) black

Tyres

- ☐ Pneum. Marathon Plus
☐ Pneum. high pressure Smooth
☐ Pneum. high pressure Plugged 24"
☐ Pn.Schwalbe Right Run 24"
☐ Solid black KIK 24"
 not available for Ultralight wheels

Pushrims

17

- ☐ Aluminium silver
☐ Aluminium dark
☐ OVAAL alum. silver 24"
☐ OVAAL alum. dark 24"
☐ OVAALG(dark w/silicon ins.)24"
☐ Chromium-plated steel 24"
☐ Rubber coated
☐ Titanium
☐ CarboLife CURVE L 24" *
☐ CL CurveL TetraGrip rubb.c. 24"*
☐ CL Quadro rubber ctd 24" *
☐ CarboLife GEKKO 24" *
☐ Silicon covers 24"
 OVAAL/S/G only in size 24"
 Only OVAAL/GEKKO available for Star6x3
 OVAAL & Carbolife not available for MTB

Spokes guards

18

- ☐ Transparent 24"
☐ Black 24"
☐ Fantasy
 code pattern
 (see catalogue of spokeguards)

Anodized set

(Light 24" hubs & rims, forks and fittings)

20

- ☐ Black
☐ Silver
☐ Light blue
☐ Red
☐ Orange
☐ Cyclamen light purple

Front castors

19

- ☐ 80 (3") solid
☐ 100 (4"x1.1/4) solid black w/plastic rim
☐ 100 (4"x1.1/4) solid with alum. rim
☐ 125 (5"x1.1/4) solid black w/plastic rim
☐ 125 (5"x1.1/4) solid with alum. rim


- ☐ **Welded Batec fastening**

21

**Lined cushion
 in foam polyurethan**


22

- ☐ 3 cm high
☐ 5 cm high



Pockets


23



- ☐ Seat pocket



- ☐ Big under the seat





- ☐ Back-rest pocket

Options



- 25 ☐ "M+F" velcro strips (hook and loop) on seat
 26 ☐ "M" velcro strips (hook) on seat
 27 ☐ "F" velcro strips (loop) on seat
 28 ☐ Chest harness
 29 ☐ Chest belt
 30 ☐ Abductor belt
 31 ☐ Catheter holder pocket
 32 ☐ Pump for pneumatics
 33 ☐ 45° pelvic belt
 34 ☐ 45° padded pelvic belt
 35 ☐ Black backpack
 36 ☐ Front protective covers

☐ **STRAIGHT SEAT & "V" shaped FRONT FRAME (def.)**



(with a "V" shaped front frame the footrest has to be **tilt adjustable**)



☐ **STRAIGHT SEAT & STRAIGHT FRONT FRAME**



(with a straight front frame all proposed footplate models in see section 4 are suitable)



☐ **CLOSING SEAT & STRAIGHT FRONT FRAME**

24



(The closing seat is suitable only if the seat width is not lower than L 36cm)
 (with a straight front frame all proposed footplate models in see section 4 are suitable)

☐ **CLOSING SEAT & "V" shaped FRONT FRAME**



(The closing seat is suitable only if the seat width is not lower than L 36cm)
 (with a "V" shaped front frame the footrest has to be **tilt adjustable**)

- ☐ Custom front opening (** from inside frame tubes)cm

Notes