




Dealer _____

Key

standard option



EOS³
HINGED front frame
 Rigid frame in Titanium



- EOS³ Fix axle EOS³ Adj axle
- L** Seat width 34 36 38 40 42 44
P Seat depth 36 38 40 42 44 46
HS Back-rest height 27 29 31 33 35 37 39 41
HA Front seat/floor height 46 48 50
HP Rear seat/floor height 40 42 44 46 48
DP Foot-rest/seat distance 37 39 41 43 45 47 49 cm

Frame colour

- Shiny titanium Natural titanium

note: frame of Eos3 cannot be varnished

Backrest:

- Tension adj nylon Breathable tension adj
 Without any backrest uph. Rigid in Carbon fibre
- AMEKH** (separate order form)
 other seating system

Upholstery

Seat:

- Fixed (stand.)
 Rigid in Carbon fibre

Backrest upholstery colour

- Black
 Coloured (nylon).....
colours: orange, red, green, cornflower, light blue
 Coloured (breath.).....
colours: red, grey

Tubular single titanium footplate

- Open Carbon fibre closed Tilt adjustable Carbon fib.closed

Pushing handles

- Without any pushing handles Fixed pushing handles Height adjustable pushing handles Foldable pushing handles

Distance between the hub centre and back-rest (COG)

- 9 cm (standard) 7 cm
 0 cm (biamputated) 11 cm
 3 cm 12 cm
 5 cm

Backrest tilt

- 89° 91° (standard) 93°

Front frame angle

- 83° A = 14 cm (default)
 75° A = 19 cm

Rear wheel camber

- 1° 2° 3° 4°
- Quick release axles Ø12mm with 42 cambered bushes (default)
 Q.r.a. Ø 1/2" (12,7 mm) w/42 camb.bushes
 Q.r.a. Ø 1/2" (12,7 mm) w/50 camb.bushes

Revolving anti-tip device

- Right Left

Side-guards

- Aluminium with carbon fibre wing Detachable CRB wing
 Aluminium w/straight carbon fibre Detachable straight CRB sideguard

Height adjustable from 22 to 28cm "L" arm-rest, detachable

- Tip-up
- Elbow-rest**
- Tubular
 Straight (26 cm long)
 Classic (23 cm long)
 Classic (30 cm long)

Brakes

- Scissor brake
 Composite brake
- Note: the model of brake might change due to technical conflicts folding the front frame

Pockets

- Seat pocket Back-rest pocket

Backrest bar

- Standard: dia 19mm Custom:
 P(bar) 7,5cm; P(bar).....cm;
 H(bar) 19cm. H(bar).....cm.

Notes

