

Dealer


Key standard option

IDRA2.0 FIX DETACH. FOOT. COD. IDR2FE **IDRA 2.0 ADJ** DETACH. FOOT. COD. IDR2RE

- L** Seat width 34 36 38 40 42 44
P Seat depth 36 38 40 42 44 46
HS Back-rest height 27 29 31 33 35 37 39 41
HA Front seat/floor height 46 48 50
HP Rear seat/floor height 40 42 44 46 48
DP Foot-rest/seat distance 37 39 41 43 45 47 49 cm

IDRA2.0

with **DETACHABLE** footrests
Aluminium rigid frame



Frame colour ①

- Glossy black Matte black
 Pearl white Pearl red
 Glossy anthracite grey
 (see "Colours" catalogue) from:
 Neon green Neon orange


Upholstery ②

- Backrest:**
 Tension adj nylon Breathable tension adj
 Without any backrest upholstery
 SAMEKH (separate order form)
 other seating system
- Seat:**
 FIX (standard, screws)
 Tension adj. (velcro)


Upholstery colour ③

- Black
 Coloured (nylon).....
colours: orange, red, green, cornflower, light blue
 Coloured (breathable).....
colours: red, grey

Detachable separate footrests ④




Push handles ⑤




- Without pushing handles Fixed pushing handles
 Height adjustable pushing handles Foldable pushing handles

Distance between the hub centre and backrest ⑥




- 9 cm (standard)** 7 cm
 0 cm (biamputated) 11 cm
 3 cm 13 cm
 5 cm

Backrest tilt ⑦



- 89°
 91° (standard)
 93°

Separate aluminium profile tilt adjustable foot-plates: ⑧



- forward
 backward
 Foot-holding belts

Rear wheel camber ⑨


- 1° 2° 3° 4°
 Quick rel. axles Ø 12mm with cambered bushes length 42 mm (default)
 Quick rel. axles Ø 1/2" (12,7mm) w/42mm bushes
 Quick rel. axles Ø 1/2" (12,7mm) w/50mm bushes

Side-guards ⑩



- Aluminium with carbon fibre wing Aluminium with straight carbon fibre

Height adjustable from 22 to 28cm "L" arm-rest, detachable and tip-up ⑪



- Elbow-rest**
 Tubular
 Straight (26 cm long)
 Classic (23 cm long)
 Classic (30 cm long)


Brakes ⑫



- Scissor brake
 Composite

- Standard
 Pull brake

Anti-tip device ⑬



- Right
 Left

Notes

