

Dealer


Key standard option

QUASAR GT Fix axle **QUASAR GT Adj axle**

- L** Seat width 34 36 38 40 42 44
P Seat depth 36 38 40 42 44 46
HS Back-rest height 27 29 31 33 35 37 39 41
HA Front seat/floor height 46 48 50
HP Rear seat/floor height 40 42 44 46 48
DP Foot-rest/seat distance 37 39 41 43 45 47 49 cm

QUASAR GT
 Titanium rigid frame

crash tested



Frame colour ①

- Shiny titanium Glossy black
 Pearl white Matte black
 Natural titanium
 (see "Colours" catalogue)
 Neon green Neon orange

Back-rest: Upholstery ②

- Tension adj nylon
 Breathable tension adj
 Without any backrest uph.
 FAMEKH
 other seating system.....
- Seat:**
 Fix (standard, screws)
 Tension adj. (velcro)
 Carbon fibre

Backrest upholstery colour ③

- Black
 Coloured (nylon).....
colours: orange, red, green, cornflower, light blue
 Coloured (breathable).....
colours: red, grey)


Single tubular titanium footplate ④

-
- Open Carbon fibre closed Tilt adjustable Carbon fibre closed

Push handles ⑤


-
- Without pushing handles Fixed pushing handles Height adjustable pushing handles Foldable pushing handles

Distance between the hub centre and back-rest ⑥




- 9 cm (standard) 7 cm
 0 cm (biamputated) 11 cm
 3 cm 13 cm
 5 cm

Backrest tilt ⑦



- 89°
 91° (standard)
 93°

Front angle ⑧



- 83° (standard) (A = 16 cm)
 75° (A = 21 cm)

Rear wheel camber ⑨

- 1° 2°
 3° 4°
- Quick rel. axles Ø 12mm with cambered bushes length 42 mm (default)
 Quick rel. axles Ø 1/2" (12,7mm) w/42mm bushes
 Quick rel. axles Ø 1/2" (12,7mm) w/50mm bushes

Side-guards ⑩





- Aluminium with carbon fibre wing
 Aluminium with straight carbon fibre

Security belt (included)


Car fastening stickers (included)

Brakes ⑫



- Scissor brake
 Composite
 Standard
 Pull brake

Anti-tip device ⑬



- Right
 Left

Notes

Dealer _____

Rear wheels (14)

- 24" Ultralight (24 spokes)
- 24" Light (24 spokes)
- 24" Star 6x3 (18 spokes, silver or black rim)
- 24" (36 spokes)
- 24" Mountain-bike with Citybike tyre
- 24" Mountain-bike with Offroad tyre
- 24" Spinergy LX
spokes colour.....(black def.)
- 25" Spinergy LX
spokes colour.....(black def.)
- 25" Light (24 spokes) black

Tyres

- Pneum. Marathon Plus
- Pneum. high pressure Smooth
- Pneum. high pressure Plugged 24"
- Pn.Schwalbe Right Run 24"
- Solid black KIK 24"
not available for Ultralight wheels

Pushrims (15)

- Aluminium silver
- Aluminium dark
- Chromium-plated steel 24"
- Rubber coated
- Titanium
- OVAAL (silver aluminium)
- OVAALS (dark aluminium)
- OVAALG (w/sil.insert)
- CarboLife CURVE L 24" *
- CL Curve L TetraGrip 24" *
- CL Quadro rubber ctd 24" *
- CarboLife GEKKO 24" *
- Silicon covers 24"


* All OVAAL & CL pushrims: not suitable for MTB wheels
 * Only OVAAL & GEKKO suitable for Star 6x3 wheels

Spokes guards (16)

- Transparent 24"
- Black 24"
- Fantasy
code pattern

(see catalogue of spokeguards)

Lined cushion (18)
 in foam polyurethan




- 3 cm high
- 5 cm high

Front castors (17)

- 80 (3") solid
- 80 (3") solid with lights
- 100 (4"x1.1/4) black solid with plastic rim
- 100 (4"x1) solid with lights
- 100 (4"x1.1/4) solid with alum. rim
- 125 (5"x1.1/4) black solid with plastic rim
- 125 (5"x1.1/4) solid with alum. rim
- 125 (5"x1) solid with lights
- 150 (6") pneumatic
- 150 (6") solid

Pockets (19)




- Small under the fix seat or
- Small hanging from the seat
- Big under the seat
- On the back-rest

Anodized set (20)
 (Rear wheel rims & hubs if LIGHT 24", sideguard supports, forks and fittings)

- Black
- Silver
- Light blue
- Red
- Orange
- Cyclamen light purple

(only 24" Light wheels)

Backrest bar (21)




- Standard:
P(bar) 7,5cm;
H(bar) 19cm;
dia 19mm.
- Custom:
P.....cm;
H.....cm;

Options

- 23 "M+F" velcro strips (hook and loop) on seat
- 24 "M" velcro strips (hook) on seat
- 25 "F" velcro strips (loop) on seat
- 26 Padded rigid cushion h 5
- 27 Chest harness
- 28 Chest belt
- 29 Abductor belt
- 30 Catheter holder pocket
- 31 Rx Lx Crutch holder
- 32 Pump for pneumatics
- 33 Black backpack
- 34 Front protective covers

STANDARD SEAT * Carbon fibre rigid seat (only suitable for standard seat)

ERGO SEAT



"Y" SHAPED FRONT FRAME (stand.)
standard: L-13cm

"V" SHAPED FRONT FRAME
standard: L-13cm
(The tilt adjustable footplate only is suitable with the "V" shaped front frame)

STRAIGHT FRONT FRAME
(Standard front frame for seat width of 34cm)

Custom front openingcm

7,0 kg configuration (FIXED)
7,4 kg configuration (ADJ.)

- Wheelchair 38x40x31 cm
- Nylon upholstery
- Tit. frame colour: Shiny titanium
- Open foot-rest
- No pushing handles
- Aluminium side guards with carbon fibre wing
- Scissor brakes
- Ultralight wheels
- Pneum. high press. smooth tyres
- Ultralight handrims
- 80 mm front castors
- No pockets

Max load: 100 kg