

Dealer


Key  standard  option

**QUASAR Fix axle**  **QUASAR Adj axle**

- L** Seat width  34  36  38  40  42  44  
**P** Seat depth  36  38  40  42  44  46  
**HS** Back-rest height  27  29  31  33  35  37  39  41  
**HA** Front seat/floor height  46  48  50  
**HP** Rear seat/floor height  40  42  44  46  48  
**DP** Foot-rest/seat distance  37  39  41  43  45  47  49 cm

# QUASAR

Titanium rigid frame



### Frame colour

- Shiny titanium  Glossy black  
 Pearl white  Matte black  
 Natural titanium  
 ..... (see "Colours" catalogue) from:  
 Neon green  Neon orange

### Upholstery

- Back-rest:**  
 Tension adj nylon  
 Breathable tension adj  
 Without any back-rest uph.  
 **FAMEKH**  
 other seating system.....  
**Seat:**  
 Fix (standard, screws)  
 Tension adj. (velcro)  
 Carbon fibre

### Backrest upholstery colour

- Black  
 Coloured (nylon).....  
colours: orange, red, green, cornflower, light blue  
 Coloured (breathable).....  
colours: red, grey

### Single tubular titanium footplate

- Open  Carbon fibre closed  Tilt adjustable Carbon fibre closed

### Push handles

- Without pushing handles  Fixed pushing handles  Height adjustable pushing handles  Foldable pushing handles

### Distance between the hub centre and back-rest

- 9 cm (standard)  7 cm  
 0 cm (biamputated)  11 cm  
 3 cm  13 cm  
 5 cm

### Backrest tilt

- 89°  91° (standard)  93°

### Front angle

- 83° (standard) (A = 16 cm)  
 75° (A = 21 cm)

### Rear wheel camber

- 1°  2°  
 3°  4°  
 Quick rel. axles Ø 12mm with cambered bushes length 42 mm (default)  
 Quick rel. axles Ø 1/2" (12,7mm) w/42mm bushes  
 Quick rel. axles Ø 1/2" (12,7mm) w/50mm bushes

### Side-guards

- Aluminium with carbon fibre wing  Detachable with CRB wing  
 Aluminium with straight CRB  Detachable with straight CRB

### Height adjustable from 22 to 28cm "L" arm-rest, detachable and tip-up

- Elbow-rest  
 Tubular  
 Straight (26 cm long)  
 Classic (23 cm long)  
 Classic (30 cm long)

### Brakes

- Scissor brake  Composite  
 Standard  Pull brake

### Anti-tip device

- Right  Left

Notes

Dealer

**Rear wheels** (14)

- 24" Ultralight (24 spokes)
- 24" Light (24 spokes)
- 24" Star 6x3 (18 spokes, silver or black)
- 24" (36 spokes)
- 24" Mountain-bike with Citybike tyre
- 24" Mountain-bike with Offroad tyre
- 24" Spinergy LX  
spokes colour.....(black def.)
- 25" Spinergy LX  
spokes colour.....(black def.)
- 25" Light (24 spokes) black

**Tyres**

- Pneum. Marathon Plus
- Pneum. high pressure Smooth
- Pneum. high pressure Plugged 24"
- Pn.Schwalbe Right Run 24"
- Solid black KIK 24"  
not available for Ultralight wheels

**Pushrims** (15)

- Aluminium silver
- Aluminium dark
- Chromium-plated steel 24"
- Rubber coated
- Titanium
- OVAAL (silver aluminium)
- OVAALS (dark aluminium)
- OVAALG (w/sil.insert)
- CarboLife CURVE L 24" \*
- CL CurveL TetraGrip rubb.c.24" \*
- CL Quadro rubber ctd 24" \*
- CarboLife GEKKO 24" \*
- Silicon covers 24"  
OVAAL/S/G only in size 24"  
Only OVAAL/GEKKO available for Star6x3  
OVAAL & Carbolife not available for MTB

**Spokes guards** (16)

- Transparent 24"
  - Black 24"
  - Fantasy  
code pattern .....
- (see catalogue of spokeguards)

**Lined cushion**  
 in foam polyurethan (18)




- 3 cm high
- 5 cm high

**Front castors** (17)

- 80 (3") solid
- 80 (3") solid with lights
- 100 (4"x1.1/4) solid with plastic rim
- 100 (4"x1) solid with lights
- 100 (4"x1.1/4) solid with alum. rim
- 125 (5"x1.1/4) solid with plastic rim
- 125 (5"x1.1/4) solid with alum. rim
- 125 (5"x1) solid with lights
- 150 (6") pneumatic
- 150 (6") solid

**Pockets** (19)



- Small under the fix seat or
- Small hanging from the seat
- Big under the seat
- On the back-rest


**Anodized set** (20)

(Rear wheel rims & hubs if LIGHT 24", sideguard supports, forks and fittings)

- Black
- Silver
- Light blue
- Red
- Orange
- Cyclamen light purple

(only 24" Light wheels)

**Backrest bar** (21)




- Standard:  
P(bar) 7,5cm;  
H(bar) 19cm;  
dia 19mm.
- Custom:  
P.....cm;  
H.....cm;

**Options**

- 23  "M+F" velcro strips (hook and loop) on seat
- 24  "M" velcro strips (hook) on seat
- 25  "F" velcro strips (loop) on seat
- 26  Padded rigid seat cushion h 5
- 27  Chest harness
- 28  Chest belt
- 29  Abductor belt
- 30  Catheter holder pocket
- 31  Rx  Lx Crutch holder
- 32  Pump for pneumatics
- 33  Suitcase holder
- 34  45° pelvic belt
- 35  45° padded pelvic belt
- 36  Black backpack
- 38  Front protective covers

**STANDARD SEAT**


\* Carbon fibre rigid seat (only suitable for standard seat)



**ERGO SEAT**




**"Y" SHAPED FRONT FRAME (stand.)**




standard:  
L-13cm

**"V" SHAPED FRONT FRAME**



standard:  
L-13cm

**STRAIGHT FRONT FRAME**



(Standard front frame for seat width of 34cm)

(The tilt adjustable footplate only is suitable with the "V" shaped front frame)

- Custom front opening .....cm

- Reinforced frame (max load 100-120 Kg)

**7,0 kg configuration (FIXED)**  
**7,4 kg configuration (ADJ.)**

- Wheelchair 38x40x31 cm
- Nylon upholstery
- Tit. frame colour: Shiny titanium
- Open foot-rest
- No pushing handles
- Aluminium side guards with carbon fibre wing
- Scissor brakes
- Ultralight wheels
- Pneum. high press. smooth tyres
- Ultralight handrims
- 80 mm front castors
- No pockets
- Max load: 100 kg