


Dealer

**Key**  standard  option

- L** Seat width  34  35  36  37  38  39  40  41  42  43  44  45  46  47  48
- P** Seat depth  34  36  38  40  42  44  46  48  50
- HS** Backrest height  27  29  31  35  37  39
- HA** Front seat/floor H  40  41  42  43  44  45  46  47  48  49  50  51  52
- HP** Rear seta/floor H  37  38  39  40  41  42  43  44  45  46  47  48  49  50
- DP** Footrest/seat dist.  35  36  37  38  39  40  41  42  43  44  45  46  47 cm

**THEMIS**  
**PRESTIGE**  
 Ergonomic Aluminium  
 custom rigid frame



- Colour:**
- Glossy black
  - Met. blue
  - Pearl white
  - Tiffany
  - Champagne
  - RAL/Changing col.....
- Upholstery**
- Tension adj. nylon backrest
  - Tension adj. breathable backrest
  - Upholst.colour: Black (default) .....
- Seat**
- Nylon black (default)
  - Carbon fibre
- Matte black  Matte Blue Klein  
 Matte purple  Matte London smoke  
 Colour.....  
(excluded camouflage, neon & special colours)  
(see "Colours" catalogue)

- Front frame**
- V shaped (allows only tilt adj. footplate)
  - Straight
- Footplate**
- Tilt adjustable Carbon fib.closed (available with V shaped or straight front frame)
  - Open (available only with straight front frame)
  - Carbon fibre closed (available only with straight front frame)

- Push handles**
- Without any pushing handles
  - Fixed pushing handles
  - Foldable pushing handles
  - Height adjustable pushing handles

- COG Distance between the hub centre and back-rest**
- 
- 0 cm
  - 1 cm
  - 2 cm
  - 3 cm
  - 4 cm
  - 5 cm
  - 6 cm
  - 7 cm
  - 8 cm
  - 9 cm
  - 10 cm
  - 11 cm
  - 12 cm
  - 13 cm
  - 14 cm
  - 15 cm

- Backrest tilt**
- 
- 89°
  - 91° (def.)
  - 93°
- Foldable and tilt adjustable backrest (default)
  - Fixed, welded backrest

- Front frame angle**
- 
- 83° (default) A = 16 cm
  - Custom A= .....cm (approx measurements)

- Carbon fiber side-guards with foldable backrest:**
- fixed with CRB wing (def.)
  - fixed with straight CRB
  - detachable with CRB wing
  - detachable with straight CRB
- Carbon fiber side-guards with welded backrest:**
- fixed with CRB wing (def.)
  - fixed with straight CRB

- Backrest bar**
- Default: P(bar depth) 7,5cm; H(bar height) 19cm; dia 21mm.
  - Custom: P.....cm; H.....cm.
- 

- Brakes**
- Composite brake
  - Scissor brake
  - Classic push brake

- Revolving anti-tip**
- Right
  - Left

- Rear wheel camber**  1°  2°  3°  4°
- Quick release axles Ø 12mm with cambered bushes 42mm long (def.)
  - Quick release axles Ø 1/2" (12,7mm) with cambered bushes 42mm long
  - Quick release axles Ø 1/2" (12,7mm) with cambered bushes 50mm long

