


Dealer \_\_\_\_\_

**Key**       standard       option

**VENUS ELITE**  
 Titanium rigid frame

- L** Seat width      34 36 38 40 42 44 46  
**P** Seat depth      36 38 40 42 44  
**HS** Back-rest height      27 29 31 33 35 37 39  
**HA** Front seat/floor height      46 48 50  
**HP** Rear seat/floor height      39 41 43 45 47  
**DP** Foot-rest/seat distance      31 33 35 37 39 41 43 45 47 49  
 cm



- Frame colour**
- Shiny titanium       Glossy black  
 Pearl white       Matte black
- Natural titanium  
 ..... (see "Colours" catalogue)  
 Neon green       Neon orange

- Upholstery**
- Tension adjustable nylon backrest  
 Tension adj. breathable backrest  
 Without any back-rest upholstery
- SAMEKH** (separate order form)  
 other seating system:.....

- Upholstery colour**
- Black  
 Coloured (nylon).....  
colours: orange, red, green, cornflower, light blue  
 Coloured (breathable).....  
colours: red, grey

- Single tubular titanium footplate**  ③
- 
- Open       Carbon fibre closed  
 Tilt adjustable Carbon fibre closed

- Push handles**  ④
- 
- Without any pushing handles       Fixed pushing handles       Height adjustable pushing handles       Foldable pushing handles

- Distance between the hub centre and back-rest**  ⑤
- 
- 4 cm       9 cm  
 5 cm       11 cm  
 7 cm (standard)

- Backrest tilt**  ⑥
- 
- 89°  
 91° (standard)  
 93°

- Frame length**  ⑦
- 
- 42,5 (stand.)  
 46,5  
 .....  
 cm

- Seat with shock absorbers**  ⑧
- 
- User's weight : ..... kg  
 It is not possible to have two anti-tip devices

- Side-guards**  ⑨
- 
- Aluminium with carbon fibre wing       Detachable  
 Aluminium with straight carbon fibre       Detachable

- "L" detachable arm-rest - h 22 cm**  ⑩
- 
- Tip-up
- Elbow rest**
- Tubular  
 Straight (26 cm long)  
 Classic (23 cm long)  
 Classic (30 cm long)

- Brakes**  ⑪
- Push brake  
 Pull brake
- Composite brake  
 Scissor brake
- 

- Anti-tip device**  ⑫
- 
- Right  
 Left

Notes

