

☐ Order N° \_\_\_\_\_ ☐ Quotation  
 Date \_\_\_\_\_ Signature \_\_\_\_\_  
 Customer Ref. \_\_\_\_\_ User's weight Kg \_\_\_\_\_


PAGE 1/2  
 May 2024

Dealer

- L** ☐ RANGE 1 W P > 36cm ☐ 34 ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46  
☐ RANGE 2 W P > 36cm ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48 ☐ 50 ☐ 52  
☐ RANGE 3 W P > 36cm (compulsory expandable stretched bar to be added) ☐ 46 ☐ 48 ☐ 50 ☐ 52 ☐ 54 ☐ 56 ☐ 58
- P** Seat depth ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46  
 Adjustable up to 6 cm(max. 50 cm) ☐ 48 ☐ 50
- HS** Back-rest height ☐ 29 ☐ 31 ☐ 33 ☐ 35 ☐ 37 ☐ 39 ☐ 41 ☐ 43 ☐ 45  
☐ 47 ☐ 49 ☐ 51 ☐ 55 ☐ 60 (stretched bar - that can grow - included)
- HA** Front seat/floor height ☐ 42 ☐ 43 ☐ 44 ☐ 46 ☐ 48 ☐ 50 ☐ 52  
☐ 54
- HP** Rear seat/floor height ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48 ☐ 50
- DP** Foot-rest/seat distance ☐ 34 ☐ 35 ☐ 37 ☐ 39 ☐ 41 ☐ 43 ☐ 45 ☐ 47 ☐ 49 ☐ 51  
 (from 17 to 32 cm foot-rest mounted on the frame, not possible with plastic footrests)
- Seat width (adjustable within the selected range)

# ALTHEA PLUS

Width & depth adjust.  
 fold.alum.frame



- Frame colour** ☐ Matte black ☐ Metallic blue ☐ Pearl red ☐ Metallic grey  
☐ Pearl white ☐ Sh.black ☐ Neon green ☐ Neon orange  
☐ ..... (see "Colours" catalogue)

- Backrest upholstery** ☐ Tension adj. nylon. back  
☐ Tension adj. breath. back  
☐ Without any back-rest upholstery  
☐ With divided seat (joint by velcro, def.)  
 Seat upholstery can grow +6cm max

- Upholstery colour** ★ (Coloured uphol. only with nylon)  
☐ Black  
☐ Red★ ☐ Orange★  
☐ Light blue★ ☐ Green★  
☐ Cornflower blue★

- ☐ **Not Detachable foot-rest**
- ☐ **Detachable foot-rest 82° angle**  
 (DP 41 cm)  
 (82° footrests do not allow any 175 front castors)
- ☐ **Detachable foot-rest 68° angle**  
 (DP 41cm)
- ☐ **Elevating foot-rest (detachable)**

- ☐ **Plastic separate**
- ☐ Cover for footplates (recommended for seat width >52cm)  
☐ Foot-holding belts (cover plates N.A.)  
☐ Heel stops (cover plates N.A.)

- Separate aluminium profile tilt adjustable footplates:**
- ☐ forward ☐ backward
- ☐ Cover for footplates (recommended for seat width >52cm)  
☐ Foot-holding belts (cover plates N.A.)

- Single aluminium profile tilt adjustable footplate:**  
 (limited range of width adjustment)  
 (not available with elevating foot-rest)
- ☐ forward ☐ backward
- ☐ Foot-holding belts  
☐ Heel stops

- Side-guards**
- ☐ Plastic ☐ Aluminium (same colour)
- Detachable DESK armrest - tip-up**
- ☐ OPEN H22 with side-guard (only if HP ≥44 cm)  
☐ OPEN H27 with side-guard (if HP≤42cm)  
☐ CLOSED H22 (from HP40/R24"x1 only if without any table fastenings)  
☐ CLOSED H27 cm from seat (compulsory for HP=42cm&R24")
- Elbow-rest**
- ☐ Classic (23 cm long)  
☐ Straight (26 cm long)  
☐ Tubular (only if open armrest)  
☐ Hemiplegic RX  
☐ Hemiplegic LX
- note: see also Hemi modular armrest ord.form

- Detachable SPORT armrest - tip-up**
- ☐ OPEN H22 with side-guard (only if HP ≥44 cm)  
☐ OPEN H27 with side-guard (if HP≤42cm)  
☐ CLOSED H22 (from HP40/R24"x1 only if without any table fastenings)  
☐ CLOSED H27 cm from seat (compulsory for HP=42cm&R24")
- Elbow-rest**
- ☐ Classic (23 cm long)  
☐ Straight (26 cm long)  
☐ Tubular (only if open armrest)  
☐ Hemiplegic RX  
☐ Hemiplegic LX
- note: see also Hemi mod.armrest ord.form

- Detachable "U" armrest - tip-up**
- ☐ OPEN H22 with side-guard (only if HP ≥44 cm)  
☐ OPEN H27 with side-guard (if HP≤42cm)  
☐ CLOSED H22 (from HP40/R24"x1 only if without any table fastenings)  
☐ CLOSED H27 cm from seat (compulsory for HP=42cm&R24")
- Elbow-rest**
- ☐ Straight (26 cm long)  
☐ Straight (36 cm long)  
☐ Tubular (only if open armrest)  
☐ Hemiplegic RX  
☐ Hemiplegic LX
- note: see also Hemi mod.armrest ord.form

- Height adjustable tip-up arm-rest (only closed)**
- ☐ Desk model  
☐ Sport model  
☐ "U" model
- Elbow-rest**
- ☐ Classic (30 cm long)  
☐ Classic (36 cm long)  
☐ Straight (26 cm long)  
☐ Straight (36 cm long)  
☐ Hemiplegic RX  
☐ Hemiplegic LX
- note: see also Hemi mod.armrest ord.form

☐ Order N° \_\_\_\_\_ ☐ Quotation  
 Date \_\_\_\_\_ Signature \_\_\_\_\_  
 Customer Ref. \_\_\_\_\_

PAGE 2/2  
 May 2024

# ALTHEA PLUS

Dealer \_\_\_\_\_

## Rear wheels

- ☐ 24" 36 spokes ☐ 12" solid\*  
☐ 22" 36 spokes ☐ 12" pneumatic\*  
☐ 20" 36 spokes (N.A. if HP=42 or HP≥ 46  
 if HP44: only curved rear frame)  
☐ 24" MTB citybike tyre  
☐ 24" MTB offroad tyre  
☐ 24" Light (24 spokes)

- ☐ Pneum. Plugged (x 1"3/8) **Tyres**  
 Pneum. high press. ☐ Plugged ☐ Smooth  
☐ Solid grey x 1" 3/8  
☐ Solid black KIK 24"  
☐ 24"/22" Pneum. Marathon Plus

## Pushrims

- ☐ Aluminium silver  
☐ Chromium-plated steel 24"  
☐ Aluminium dark  
☐ Rubber coated  
☐ OVAAL al.silver (\*)  
☐ OVAALS al.dark (\*)  
☐ OVAALG al.dark w/sil. insert (\*)  
☐ Carbolife CURVE L 24"(\*)  
☐ Carbolife GEKKO 24" (\*)  
☐ Silicon covers 24"  
 (\*) not available with MTB wheels

## Spoke guards


- ☐ Transparent  
☐ Black  
☐ Fantasy  
 code pattern: .....  
 (see catalogue of spokeguards)

## Front castors


- ☐ 100 (4"x1.1/4) solid - plastic rim  
☐ 100 (4"x1.1/4) solid - alum. rim  
☐ 125 (5"x1.1/4) solid - plastic rim  
☐ 125 (5"x1.1/4) solid - alum. rim  
☐ 150 (6") pneumatic  
☐ 150 (6") solid  
☐ 175 (7") solid \*\*  
☐ 175 (7") pneumatic \*\*  
 (\*\*) only if HA 48cm as minimum and 68°detach.footr.

## Rear frame


- ☐ STRAIGHT rear frame with small wheels for narrow passages (default)




- ☐ STRAIGHT rear frame WITHOUT any small wheels f/narrow passages



- ☐ STRAIGHT rear frame with LEVER-activated narrow passage wheels  
 € 134,00



- ☐ CURVED rear frame



(Only if rear seat/floor height HP is between 44 and 48 cm)

## Unbalancing system

- ☐ Right ☐ Right  
☐ Left ☐ Left  
 (Only with straight rear frame) (Only with curved rear frame)

## Anti-tip device

- (Only with straight rear frame)  
☐ Right  
☐ Left


## Revolving anti-tip device

- (Only with curved rear frame)  
☐ Right  
☐ Left


## Single drive

- ☐ Right  
☐ Left

(camber can be only 0°)  
 (rear wheels can be only 22" or 24")  
 (the axle can grow +2cm or +4cm)




## Extended wheel plate




## Push handles

- ☐ Standard ☐ Height adjustable(\*) ☐ Hinged back-rest (no stretched bar)




## Distance between the hub centre and backrest

- ☐ -2 cm (ext. wheel plate)  
☐ 3,5 cm  
☐ 5,5 cm  
☐ 7,5 cm



(with HP=44 cm)


## Adjustable swing away Lateral supports



- ☐ Right  
☐ Left  
 Specify the height "h" from the center of the swing away to the seat  
 ..... cm

## Head-rest

**Padded**  
 (cannot grow in width)




- ☐ Fixed  
☐ Height adjustable

## Orbital


- ☐ Anatomic



- ☐ Moulded



- ☐ With lateral supports



Specify "h" & "d" where requested  
 Standard: h = 20 cm d = 15 cm

h = ..... cm d = ..... cm

## Table

- ☐ Transparent  
☐ RxOLx Hemiplegic transp.  
☐ Plastic grey  
☐ Soft padded

Specify the fastening:  
☐ single ☐ double



(Only with arm-rest) OXL if L>50cm

## Lined cushion in foam polyurethan

(cannot grow in width)



- ☐ 3 cm high  
☐ 5 cm high


## Brakes

- ☐ Classic - short handle  
☐ Classic - long handle  
☐ Drum  
☐ Right single brake  
☐ Left single brake  
☐ Drum+classic

## Options

- 29 ☐ "M+F" velcro stripon sea  
 30 ☐ "M" velcro strip (hook) on seat  
 31 ☐ "F" velcro strip (loop) on seat  
 32 ☐ Seat/backrest connection flap  
 33 ☐ Chest harness  
 34 ☐ Chest belt  
 35 ☐ Abductor belt  
 36 ☐ Front protective covers  
 37 ☐ Stretched bar (width adj.)  
 38 ☐ R OL Crutch hold./single  
 39 ☐ R OL Crutch hold./double  
 40 ☐ Pump for pneumatics

## Abductor



- ☐ 1 ☐ 3 ☐ 5 ☐ 7 cm  
☐ ..... cm

Choose or specify the height from seat cover

- 41 ☐ 45° pelvic belt  
 42 ☐ 45° padded pelvic belt  
 43 ☐ Black backpack  
 44 ☐ IV bar + poles  
 45 ☐ IV poles  
 46 ☐ Q.rel.axles Ø12mm w/42mm bushes (def.)  
 47 ☐ Q.rel.a. Ø1/2" (12,7mm) w/42mm bushes  
 48 ☐ Q.rel.a. Ø1/2" (12,7mm) w/50mm bushes