

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____ User's weight Kg _____

Key


☐ standard☐ option

L Seat width	<input type="checkbox"/> 34 <input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46 <input type="checkbox"/> 48 <input type="checkbox"/> 50 <input type="radio"/> 52 <input type="radio"/> 56 <input type="radio"/> 60 (including eventual stretched bar, tension adj. upholstery & cover for footplate according to setup)
P Seat depth	<input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46 <input type="radio"/> 48 <input type="radio"/> 50 (Each depth is adjustable +2, +4, +6 cm up to max. 50 cm)
HS Back-rest height	<input type="checkbox"/> 29 <input type="checkbox"/> 31 <input type="checkbox"/> 33 <input type="checkbox"/> 35 <input type="checkbox"/> 37 <input type="checkbox"/> 39 <input type="checkbox"/> 41 <input type="checkbox"/> 43 <input type="checkbox"/> 45 <input type="radio"/> 47 <input type="radio"/> 49 <input type="radio"/> 51 <input type="radio"/> 55 <input type="radio"/> 60 (stretched bar included)
HA Front seat/floor height	<input type="checkbox"/> 42 <input type="checkbox"/> 43 <input type="checkbox"/> 44 <input type="checkbox"/> 46 <input type="checkbox"/> 48 <input type="checkbox"/> 50 <input type="checkbox"/> 52 <input type="radio"/> 54
HP Rear seat/floor height	<input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46 <input type="checkbox"/> 48 <input type="checkbox"/> 50
DP Foot-rest/seat distance	<input type="checkbox"/> 34 <input type="checkbox"/> 35 <input type="checkbox"/> 37 <input type="checkbox"/> 39 <input type="checkbox"/> 41 <input type="checkbox"/> 43 <input type="checkbox"/> 45 <input type="checkbox"/> 47 <input type="checkbox"/> 49 <input type="checkbox"/> 51 (from 17 to 32 cm foot-rest mounted on the frame, not possible with plastic footrests)

cm

ALTHEA

Alum. foldable frame
with double crossbar



Frame colour

- ☐ Metallic blue ☐ Pearl red ☐ Metallic grey ☐ Matte black
☐ Pearl white ☐ Shiny black ☐ Neon green ☐ Neon orange
☐ (see "Colours" catalogue)

Upholstery

- ☐ Backrest Nylon upholstery (max up to L 50)
☐ Tension adj. NYLON backrest uphol.
☐ Tens.adj. BREATHABL. backrest uphol.
☐ Without any back-rest upholstery
☐ With divided seat (joint by velcro - std)
☐ With one entire seat (not depth adjustable)

Upholstery colour


- ★ (Coloured uphol. only with nylon)
☐ Black
☐ Red★ ☐ Orange★
☐ Light blue★ ☐ Green★
☐ Cornflower blue★

☐ Not Detachable foot-rest☐ Detachable foot-rest 82° angle


(DP 41 cm)
(82° footrests do not allow any 175 front castors)

☐ Detachable foot-rest 68° angle

(DP 41cm)

☐ Elevating foot-rest (detachable)☐ Plastic separate


(only L = from 34 to 56 cm)



- ☐ Cover f/footplates ★
☐ Foot-holding belts (W50max)
☐ Heel stops (W50max)
 (* eventual cover included for W 52 to 60)


Single aluminium profile tilt adjustable footplate:

(not available with elevating foot-rest)



- ☐ Foot-holding belts
☐ Heel stops


Separate aluminium profile tilt adjustable footplates:



- ☐ Foot-holding belts (W50max)
 (* eventual cover included for W 52 to 60)


☐ Automatic closure footplate in carbon

(not available with detachable or elevating foot-rest; only W= from 36 to 46cm)



- ☐ Foot-holding belts

Side-guards



☐ Plastic


☐ Aluminium
(same colour)

Detachable DESK armrest - tip-up

- ☐ OPEN H22 with side-guard (only if HP ≥44 cm)
☐ OPEN H27 with side-guard (if HP≤42cm)



- ☐ CLOSED H22 (from HP40/R24"x1 only if without any table fastenings)
☐ CLOSED H27 cm from seat (compulsory for HP=42cm&R24")




Elbow-rest

- ☐ Classic (23 cm long)
☐ Straight (26 cm long)
☐ Tubular (only if open armrest)
☐ Hemiplegic RX
☐ Hemiplegic LX


note: see also Hemi modular armrest ord.form

Detachable SPORT armrest - tip-up

- ☐ OPEN H22 with side-guard (only if HP ≥44 cm)
☐ OPEN H27 with side-guard (if HP≤42cm)



- ☐ CLOSED H22 (from HP40/R24"x1 only if without any table fastenings)
☐ CLOSED H27 cm from seat (compulsory for HP=42cm&R24")




Elbow-rest

- ☐ Classic (23 cm long)
☐ Straight (26 cm long)
☐ Tubular (only if open armrest)
☐ Hemiplegic RX
☐ Hemiplegic LX


note: see also Hemi mod.armrest ord.form

Detachable "U" armrest - tip-up

- ☐ OPEN H22 with side-guard (only if HP ≥44 cm)
☐ OPEN H27 with side-guard (if HP≤42cm)



- ☐ CLOSED H22 (from HP40/R24"x1 only if without any table fastenings)
☐ CLOSED H27 cm from seat (compulsory for HP=42cm&R24")




Elbow-rest

- ☐ Straight (26 cm long)
☐ Straight (36 cm long)
☐ Tubular (only if open armrest)
☐ Hemiplegic RX
☐ Hemiplegic LX

note: see also Hemi mod.armrest ord.form

Height adjustable tip-up armrest (only closed)



Elbow-rest

- ☐ Classic (30 cm long)
☐ Classic (36 cm long)
☐ Straight (26 cm long)
☐ Straight (36 cm long)
☐ Hemiplegic RX
☐ Hemiplegic LX

note: see also Hemi modular armrest ord.form

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____

PAGE 2/2
 May 2024

ALTHEA

Dealer _____

Rear wheels ⑭

- ☐ 24" 36 spokes (def.) ☐ 12" solid*
☐ 22" 36 spokes ☐ 12" pneumatic*
☐ 20" 36 spokes (* N.A. if HP=42 or HP≥46cm if HP44: only curved rear frame)
☐ 24" MTB citybike tyre
☐ 24" MTB offroad tyre
☐ 24" Light (24 spokes)
☐ Pneum. Plugged (x 1.3/8") **Tyres**
☐ Plugged ☐ Smooth Pneum. high press.
☐ Solid grey x 1.3/8"
☐ Solid black KIK 24"
☐ 24" Pneum. Marathon Plus

Pushrims ⑮

- ☐ Aluminium silver
☐ Chromium-plated steel 24"
☐ Aluminium dark
☐ Rubber coated
☐ Carbolife CURVE L 24"(*)
☐ OVAAL silver alum. (*)
☐ OVAALS dark alum. (*)
☐ OVAALG dark w/sil. insert (*)
☐ Carbolife GEKKO 24" (*)
☐ Silicon covers 24"
 (*) not available with MTB wheels

Spokeguards ⑯


- ☐ Transparent
☐ Black
☐ Fantasy
 code pattern:
 (see catalogue of spokeguards)

Castor wheels ⑰


- ☐ 100 (4"x1.1/4) solid - plastic rim
☐ 100 (4"x1.1/4) solid - alum. rim
☐ 125 (5"x1.1/4) solid - plastic rim
☐ 125 (5"x1.1/4) solid - alum. rim
☐ 150 (6") pneumatic
☐ 150 (6") solid
☐ 175 (7") solid **
☐ 175 (7") pneumatic **
 (**) only if HA ≥ 48cm and with 68° detach. footrests

Rear frame ⑱


- ☐ STRAIGHT rear frame with small wheels for narrow passages (default)



- ☐ STRAIGHT rear frame WITHOUT any small wheels f/narrow passages




- ☐ STRAIGHT rear frame with LEVER-activated narrow passage wheels



Max
user weight
100 kg

(Only if rear seat/floor height
HP is between 44 and 48 cm)

- ☐ CURVED rear frame



Unbalancing system ⑲

- ☐ Right ☐ Right
☐ Left ☐ Left
 (Only with straight rear frame) (Only with curved rear frame)

Anti-tip device (Only with straight rear frame)



- ☐ Right
☐ Left

Revolving anti-tip device (Only with curved rear frame)




- ☐ Right
☐ Left


Single drive ⑳

(only with 24" or 22" rear wheels)

- ☐ Right
☐ Left




Extended wheel plate ㉑



Push handles ㉒

- ☐ Standard ☐ Height adjustable(*) ☐ Hinged back-rest(*) (no stretched bar)



(*) tension adj. upholstery is compulsory

COG Distance between the hub centre and backrest ㉓


- ☐ -2 cm (ext. wheel plate)
☐ 3,5 cm
☐ 5,5 cm
☐ 7,5 cm



(with HP=44 cm)

Adjustable swing away Lateral supports ㉔

(tension adj. upholstery is compulsory)




- ☐ Right
☐ Left

Specify the height "h"
from the center of the
swing away to the seat
.....
cm

Head-rest ㉕

Padded




- ☐ Fixed
☐ Height adjustable

Specify "h" & "d" where requested
 Standard: h = 20 cm d = 15 cm


h = cm d = cm

Orbital

- ☐ Anatomic



- ☐ Moulded



- ☐ With lateral supports




Table ㉖

- ☐ Transparent
☐ Rx OLx Hemiplegic transp.
☐ Plastic grey
☐ Soft padded

Specify the fastening:
☐ single ☐ double



(Only with arm-rest) OXL if L > 50cm

Lined cushion in foam polyurethane ㉗



- ☐ 3 cm high
☐ 5 cm high


Brakes ㉘

- ☐ Classic - short handle
☐ Classic - long handle
☐ Drum
☐ Right single brake
☐ Left single brake
☐ Drum+classic

Options

- 29 ☐ "F" velcro strip (loop) on seat
 30 ☐ "M" velcro strip (hook) on seat
 31 ☐ "M" velcro for Vicair
 32 ☐ Seat/back-rest connection flap
 33 ☐ Chest harness
 34 ☐ Chest belt
 35 ☐ Abductor belt
 36 ☐ Front protective covers
 37 ☐ Stretched bar
 38 ☐ Catheter holder pocket (note: only with entire seat)
 39 ☐ R ☐ L Crutch hold./single
 40 ☐ R ☐ L Crutch hold./double
 41 ☐ Pump for pneumatics

Abductor ㉙



- ☐ 01 ☐ 03 ☐ 05 ☐ 07 cm

- ☐ cm

Choose or specify the height from seat upholstery

- 42 ☐ Suitcase holder
 43 ☐ 45° pelvic belt
 44 ☐ 45° padded pelvic belt
 45 ☐ Black backpack
 46 ☐ IV bar + poles
 47 ☐ IV poles
 48 ☐ Padded rigid cushion(W52max)
 49 ☐ Padded rigid cushion(W56-60)
 50 ☐ Q.rel.axles Ø12mm w/42mm bushes (def.)
 51 ☐ Q.rel.a. Ø1/2" (12,7mm) w/42mm bushes
 52 ☐ Q.rel.a. Ø1/2" (12,7mm) w/50mm bushes
 53 ☐ **Stickers for car fastenings** (escl.: load > 120Kg, R12", UltraL.RW, hinged bckr)