

☐ Order N° \_\_\_\_\_ ☐ Quotation  
 Date \_\_\_\_\_ Signature \_\_\_\_\_  
 Customer Ref. \_\_\_\_\_ User's weight Kg \_\_\_\_\_

**Key** ☐ standard  
☐ option

**PAGE 1/2**  
 May 2024

Dealer \_\_\_\_\_

**L** ☐ RANGE 1 W ☐ 28 ☐ 30 ☐ 32 ☐ 34 ☐ 36 ☐ 38  
☐ RANGE 2 W P > 36cm ☐ 32 ☐ 34 ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44  
☐ RANGE 3 W P > 36cm ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48 ☐ 50 ☐ 52

Seat width (adjustable within the selected range)

**P** Seat depth ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46  
 Adjustable up to 6cm (max.50cm) ☐ 48 ☐ 50

**HS** Back-rest height ☐ 29 ☐ 31 ☐ 33 ☐ 35 ☐ 37 ☐ 39 ☐ 41 ☐ 45  
☐ 47 ☐ 49 ☐ 51 ☐ 55 ☐ 60 (included stretched bar that can be enlarged)

**HA** Front seat/floor height ☐ 38(fc80max) ☐ 39(fc100max) ☐ 40(fc125max) ☐ 42(fc150max) ☐ 45(fc175max)  
 with MAX size of front castors (fc)  
 (NOTE: 175 front castors do not allow 82° detachable footrests) ☐ 44(fc150max) ☐ 46(fc175max)


**HP** Rear seat/floor height ☐ 34 ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 45 ☐ 46 ☐ 48  
 NOTE: 24" rear wheels (if HP min.38/max 48) 22" (if HP min.36/max 45) 20" (if HP min.34/max 43)

**DP** Foot-rest/seat distance ☐ 34 ☐ 35 ☐ 37 ☐ 39 ☐ 41 ☐ 43 ☐ 45  
 (from 17 to 32 cm footrest directly assembled on the frame) ..... cm

# ALTHEA19PLUS

## LOWERED SET-UP

Width & depth adjust.  
 foldable alum.frame



**Frame colour** ☐ Matte black ①  
☐ Metallic blue ☐ Pearl red ☐ Metallic grey  
☐ Pearl white ☐ Glossy black  
☐ .....from: Neon colour .....  
 (see "Colours" catalogue)

**Backrest upholstery** ②  
☐ Tension adj. nylon. back  
☐ Tension adj. breath. back  
☐ Without any backrest upholstery  
☐ With divided seat (joint by velcro, def.)  
 Seat upholstery can grow +6cm max

**Upholstery colour** ③  
 ★ (Coloured uphol. only with nylon)  
☐ Black  
☐ Red★ ☐ Orange★  
☐ Light blue★ ☐ Green★  
☐ Cornflower blue★

☐ **Not Detachable foot-rest** ④

☐ **Detachable foot-rest 82° angle**  
 (DP 41 cm)  
 (STANDARD)

☐ **Detachable foot-rest 68° angle**  
 (DP 41 cm)

☐ **Elevating foot-rest (detachable)**

☐ **Plastic separate** ⑤

☐ Cover for footplates  
 (recommended for seat width >52cm)  
☐ Foot-holding belts (cover plates N.A.)  
☐ Heel stops (cover plates N.A.)

**Separate aluminium profile tilt adjustable footplates:** ⑥

☐ forward ☐ backward

**Single aluminium profile tilt adjustable footplate:** ⑦  
 (limited range of width adjustment)  
 (not available with elevating foot-rest)

☐ forward ☐ backward

☐ Foot-holding belts  
☐ Heel stops

**Side-guards** ⑧

☐ Plastic

☐ Aluminium  
 (same colour)

**Detachable DESK armrest - tip-up** ⑨

☐ OPEN H22 with side-guard (only if HP ≥44 cm)  
☐ OPEN H27 with side-guard (if HP ≤42cm)

☐ CLOSED H22 (from HP40/R24"x1 only if without any table fastenings)  
☐ CLOSED H27 cm from seat (compulsory for HP=42cm&R24")

**Elbow-rest**

☐ Classic (23 cm long)  
☐ Straight (26 cm long)  
☐ Tubular (only if open armrest)  
☐ Hemiplegic RX  
☐ Hemiplegic LX

note: see also Hemi modular armrest ord.form

**Detachable SPORT armrest - tip-up** ⑩

☐ OPEN H22 with side-guard (only if HP ≥44 cm)  
☐ OPEN H27 with side-guard (if HP ≤42cm)

☐ CLOSED H22 (from HP40/R24"x1 only if without any table fastenings)  
☐ CLOSED H27 cm from seat (compulsory for HP=42cm&R24")

**Elbow-rest**

☐ Classic (23 cm long)  
☐ Straight (26 cm long)  
☐ Tubular (only if open armrest)  
☐ Hemiplegic RX  
☐ Hemiplegic LX

note: see also Hemi mod.armrest ord.form

**Detachable "U" armrest - tip-up** ⑪

☐ OPEN H22 with side-guard (only if HP ≥44 cm)  
☐ OPEN H27 with side-guard (if HP ≤42cm)

☐ CLOSED H22 (from HP40/R24"x1 only if without any table fastenings)  
☐ CLOSED H27 cm from seat (compulsory for HP=42cm&R24")

**Elbow-rest**

☐ Straight (26 cm long)  
☐ Straight (36 cm long)  
☐ Tubular (only if open armrest)  
☐ Hemiplegic RX  
☐ Hemiplegic LX

note: see also Hemi mod.armrest ord.form

**Height adjustable tip-up arm-rest (only closed)** ⑫

☐ Desk model  
☐ Sport model  
☐ "U" model

**Elbow-rest**

☐ Classic (30 cm long)  
☐ Classic (36 cm long)  
☐ Straight (26 cm long)  
☐ Straight (36 cm long)  
☐ Hemiplegic RX  
☐ Hemiplegic LX

note: see also Hemi mod.armrest ord.form

☐ Order N° \_\_\_\_\_ ☐ Quotation  
 Date \_\_\_\_\_ Signature \_\_\_\_\_  
 Customer Ref. \_\_\_\_\_

PAGE 2/2  
 May 2024

# ALTHEA19PLUS

## ALTHEA LOWERED SET-UP

Dealer \_\_\_\_\_

### Rear wheels ⑬

- ☐ 24" 36 spokes ☐ 12" solid\*  
☐ 22" 36 spokes ☐ 12" pneumatic\*  
☐ 20" 36 spokes (N.A. if HP=42 if HP44: only curved rear frame)  
☐ 24" MTB citybike tyre  
☐ 24" MTB offroad tyre  
☐ 24" Light (24 spokes)

- ☐ Pneum. Plugged (x 1" 3/8) **Tyres**  
 Pneum. high press. ☐ Plugged ☐ Smooth  
☐ Solid grey x 1" 3/8  
☐ Solid black KIK 24"  
☐ 24" Pneum. Marathon Plus

### Pushrims ⑭

- ☐ Aluminium silver  
☐ Chromium-plated steel 24"  
☐ Aluminium dark  
☐ Rubber coated  
☐ OVAAL al.silver (\*)  
☐ OVAALS al.dark (\*)  
☐ OVAALG al.dark w/sil. insert (\*)  
☐ Carbolife CURVE L 24"(\*)  
☐ Carbolife GEKKO 24" (\*)  
☐ Silicon covers 24"  
 (\*) not available with MTB wheels

### Spoke guards ⑮

- ☐ Transparent  
☐ Black  
☐ Fantasy  
 code pattern: .....  
 (see catalogue of spokeguards)

### Front castors ⑯

- ☐ 100 (4"x1.1/4) solid - plastic rim  
☐ 100 (4"x1.1/4) solid - alum. rim  
☐ 125 (5"x1.1/4) solid - plastic rim  
☐ 125 (5"x1.1/4) solid - alum. rim  
☐ 150 (6") pneumatic  
☐ 150 (6") solid  
☐ 175 (7") solid  
☐ 175 (7") pneumatic

### Rear frame ⑰

- ☐ STRAIGHT rear frame with small wheels for narrow passages (default)  
☐ STRAIGHT rear frame WITHOUT any small wheels f/narrow passages  
☐ STRAIGHT rear frame with LEVER-activated narrow passage wheels  
☐ CURVED rear frame  
 Max user weight 100 kg  
 (Only if rear seat/floor height HP is between 44 and 48 cm)  
 (It is not possible to assemble the small wheels for narrow passages)

### Unbalancing system ⑱


- ☐ Right ☐ Right  
☐ Left ☐ Left  
 (Only with straight rear frame) (Only with curved rear frame)

- Anti-tip device** (Only with straight rear frame)  
☐ Right  
☐ Left


- Revolving anti-tip device**  
 (Only with curved rear frame)  
☐ Right  
☐ Left

### Single drive ⑲

- ☐ Right  
☐ Left  
 (camber can be only 0°)  
 (rear wheels can be only 22" or 24")  
 (the axle can grow +2cm or +4cm)




### Extended wheel plate ⑳



### Push handles ㉑

- ☐ Standard ☐ Height adjustable  
☐ Hinged backrest (no stretched bar) (max L=50)  
 (with HP=44 cm)




### Distance between the hub centre and back-rest ㉒

- ☐ -2 cm (ext. wheel plate)  
☐ 3,5 cm  
☐ 5,5 cm  
☐ 7,5 cm  
 (with HP=44 cm)




### Adjustable swing away Lateral supports ㉓

- ☐ Right ☐ Left  
 Specify the height "h" from the center of the swing away to the seat  
 ..... cm




### Head-rest ㉔

- Padded** (cannot grow in width)  
☐ Fixed  
☐ Height adjustable  
**Orbital**  
☐ Anatomic  
☐ Moulded  
☐ With lateral supports  
 Specify "h" & "d" where requested  
 Standard: h = 20 cm d = 15 cm  
 h = ..... cm d = ..... cm




### Table ㉕

- ☐ Transparent  
☐ RxOLx Hemiplegic transp.  
☐ Plastic grey  
☐ Soft padded  
 Specify the fastening:  
☐ single ☐ double



### Lined cushion in foam polyurethane ㉖



- ☐ 3 cm high  
☐ 5 cm high


### Brakes ㉗

- ☐ Classic - short handle  
☐ Classic - long handle  
☐ Drum  
☐ Right single brake  
☐ Left single brake  
☐ Drum+classic

### Options ㉘

- 29 ☐ "M+F" velcro on seat  
 30 ☐ "M" velcro strip (hook) on seat  
 31 ☐ "F" velcro strip (loop) on seat  
 32 ☐ Seat/back-rest connection flap  
 33 ☐ Chest harness  
 34 ☐ Chest belt  
 35 ☐ Abductor belt  
 36 ☐ Front protective covers  
 37 ☐ Stretched bar (width adjust.)  
 38 ☐ R OL Crutch hold./single  
 39 ☐ R OL Crutch hold./double  
 40 ☐ Pump for pneumatics

### Abductor ㉙



- ☐ 1 ☐ 3 ☐ 5 ☐ 7 cm  
☐ ..... cm  
 Choose or specify the height from seat upholstery

Choose or specify the height from seat upholstery

- 41 ☐ 45° pelvic belt  
 42 ☐ 45° padded pelvic belt  
 43 ☐ Black backpack  
 44 ☐ IV bar + poles  
 45 ☐ IV poles

- ☐ Quick rel. axles Ø 12 mm with 42mm cambered bushes (default)  
☐ Quick rel. axles Ø 1/2" (12,7 mm) with 42mm cambered bushes  
☐ Quick rel. axles Ø 1/2" (12,7 mm) with 50mm cambered bushes