

Dealer _____

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____ User's weight Kg _____

Key ☐ standard
☐ option


PAGE 1/2
 May 2024

- L** Seat width ☐34 ☐36 ☐38 ☐40 ☐42 ☐44 ☐46 ☐48 ☐50
☐52* ☐56* ☐60 (including eventual stretched bar, tension adj. upholstery & cover for footplate according to setup)
- P** Seat depth ☐36 ☐38 ☐40 ☐42 ☐44 ☐46
 Adjustable up to 6 cm (max.50cm) ☐48 ☐50
- HS** Back-rest height ☐29 ☐31 ☐33 ☐35 ☐37 ☐39 ☐41 ☐43 ☐45
☐47 ☐49 ☐51 ☐55 ☐60 (included stretched bar)
- HA** Front seat/floor height ☐38(fc80max) ☐39(fc100max) ☐40(fc125max) ☐42(fc150max) ☐45(fc175max)
 (with MAX size of front castors)
 (NOTE: 175 front castors do not allow 82° detachable footrests) ☐44(fc150max) ☐46(fc175max)
- HP** Rear seat/floor height ☐34 ☐36 ☐38 ☐40 ☐42 ☐44 ☐45 ☐46 ☐48
 NOTE: 24" rear wheels (if HP min.38/max 48) 22" (if HP min.36/max 45) 20" (if HP min.34/max 43)
- DP** Foot-rest/seat distance ☐34 ☐35 ☐37 ☐39 ☐41 ☐43 ☐45
 (from 17 to 32 cm footrest directly assembled on the frame) cm

ALTHEA.19

LOWERED SET-UP

alum.foldable frame



- Frame colour** ☐ Matte black ☐ Metallic blue ☐ Pearl red ☐ Metallic grey
☐ Pearl white ☐ Glossy black
☐ Neon colour
 (see "Colours" catalogue)

- Upholstery** ☐ Nylon (up to L50max)
☐ Tension adj. nylon. backrest up.
☐ Tension adj. breath. backrest up.
☐ Without any backrest upholstery
☐ With divided seat (joint by velcro, def.)
☐ With one entire seat (not depth adjustable)

- Upholstery colour** ☐ Black
 * (Coloured uphol. only with nylon)
☐ Red* ☐ Orange*
☐ Light blue* ☐ Green*
☐ Cornflower blue *

- ☐ **Not Detachable foot-rest** ④
- ☐ **Detachable foot-rest 82° angle**
 (DP 41 cm)
 (STANDARD) ④
- ☐ **Detachable foot-rest 68° angle**
 (DP 41 cm) ④
- ☐ **Elevating foot-rest (detachable)** ④

- ☐ **Plastic separate** ⑤
- ☐ Cover for footplates*
☐ Foot-holding belts (up to W50)
☐ Heel stops (up to L50)
 * Eventual cover f/footplates included for W52-60

- Separate aluminium profile tilt adjustable foot-plates:** ⑥
- ☐ forward ☐ backward
- ☐ Foot-holding belts (up to W50)
 * Eventual cover f/footplates included for W52-60

- Single aluminium profile tilt adjustable foot-plate:** ⑦
- (not available with elevating foot-rest)
- ☐ forward ☐ backward
- ☐ Foot-holding belts
☐ Heel stops

- ☐ **Automatic closure foot-plate in carbon** ⑧
- (not available with detachable or elevating foot-rest; only L= from 36 to 46 cm)
- ☐ Foot-holding belts

- Side guards** ⑨
- ☐ Plastic ☐ Aluminium (same colour) ⑨

- Detachable DESK armrest - tip-up** ⑩
- ☐ OPEN H22 with side-guard (only if HP ≥44 cm)
☐ OPEN H27 with side-guard (if HP ≤42cm)
☐ CLOSED H22 (from HP40/R24"x1 only if without any table fastenings)
☐ CLOSED H27 27 cm from seat (compulsory for HP=42cm&R24")
- Elbow-rest**
☐ Classic (23 cm long)
☐ Straight (26 cm long)
☐ Tubular (only if open armrest)
☐ Hemiplegic RX
☐ Hemiplegic LX
 note: see also Hemi modular armrest ord.form

- Detachable SPORT armrest - tip-up** ⑪
- ☐ OPEN H22 with side-guard (only if HP ≥44 cm)
☐ OPEN H27 with side-guard (if HP ≤42cm)
☐ CLOSED H22 (from HP40/R24"x1 only if without any table fastenings)
☐ CLOSED H27 27 cm from seat (compulsory for HP=42cm&R24")
- Elbow-rest**
☐ Classic (23 cm long)
☐ Straight (26 cm long)
☐ Tubular (only if open armrest)
☐ Hemiplegic RX
☐ Hemiplegic LX
 note: see also Hemi mod.armrest ord.form

- Detachable "U" armrest - tip-up** ⑫
- ☐ OPEN H22 with side-guard (only if HP ≥44 cm)
☐ OPEN H27 with side-guard (if HP ≤42cm)
☐ CLOSED H22 (from HP40/R24"x1 only if without any table fastenings)
☐ CLOSED H27 27 cm from seat (compulsory for HP=42cm&R24")
- Elbow-rest**
☐ Straight (26 cm long)
☐ Straight (36 cm long)
☐ Tubular (only if open armrest)
☐ Hemiplegic RX
☐ Hemiplegic LX
 note: see also Hemi mod.armrest ord.form

- Height adjustable tip-up armrest (only closed)** ⑬
- ☐ Desk model ☐ Sport model ☐ "U" model ⑬
- Elbow-rest**
☐ Classic (30 cm long)
☐ Classic (36 cm long)
☐ Straight (26 cm long)
☐ Straight (36 cm long)
☐ Hemiplegic RX
☐ Hemiplegic LX
 note: see also Hemi modular armrest ord.form

☐ Order N° _____ ☐ Quotation
 Date _____ Signature _____
 Customer Ref. _____

PAGE 2/2
 May 2024

ALTHEA19 LOWERED SET-UP

Dealer

Rear wheels 14

- ☐ 24" 36 spokes ☐ 12" solid*
☐ 22" 36 spokes ☐ 12" pneumatic*
☐ 20" 36 spokes (N.A. if HP=42 if HP44: only curved rear frame)
☐ 24" MTB citybike tyre
☐ 24" MTB offroad tyre
☐ 24" Light (24 spokes)
- ☐ Pneum. Plugged (x 1"3/8) **Tyres**
 Pneum. high press. ☐ Plugged ☐ Smooth
☐ Solid grey x 1" 3/8
☐ Solid black KIK 24"
☐ 24" Pneum. Marathon Plus

Pushrims 15

- ☐ Aluminium silver
☐ Chromium-plated steel 24"
☐ Aluminium dark
☐ Rubber coated
☐ Carbolife CURVE L 24"(*)
☐ OVAAL silver alum. (*)
☐ OVAALS dark alum. (*)
☐ OVAALG dark w/sil. insert (*)
☐ Carbolife GEKKO 24" (*)
☐ Silicon covers 24"
 (*) not available with MTB wheels

Spoke guards 16

- ☐ Transparent
☐ Black
☐ Fantasy
 code pattern:
 (see catalogue of spokeguards)

Front castors 17

- ☐ 100 (4"x1.1/4) solid - plastic rim
☐ 100 (4"x1.1/4) solid - alum. rim
☐ 125 (5"x1.1/4) solid - plastic rim
☐ 125 (5"x1.1/4) solid - alum. rim
☐ 150 (6") pneumatic
☐ 150 (6") solid
☐ 175 (7") solid
☐ 175 (7") pneumatic

Rear frame 18

- ☐ STRAIGHT rear frame with small wheels for narrow passages (default)
☐ STRAIGHT rear frame WITHOUT any small wheels f/narrow passages
☐ STRAIGHT rear frame with LEVER-activated narrow passage wheels
☐ CURVED rear frame
- Max user weight 100 kg
 (Only if rear seat/floor height HP is between 44 and 48 cm)
 (It is not possible to assemble the small wheels for narrow passages)

Unbalancing system 19

- ☐ Right ☐ Right
☐ Left ☐ Left
 (Only with straight rear frame) (Only with curved rear frame)

Anti-tip device (Only with straight rear frame)


- ☐ Right
☐ Left

Revolving anti-tip device (Only with curved rear frame)


- ☐ Right
☐ Left

Single drive 20

- ☐ Right
☐ Left




Extended wheel plate 21



Push handles 22

- ☐ Standard ☐ Height adjustable(*) ☐ Hinged backrest(*)
 (no stretched bar)



(*) tension adj. upholstery is compulsory


Distance between the hub centre and back-rest 23

- ☐ -2 cm (ext. wheel plate)
☐ 3,5 cm
☐ 5,5 cm
☐ 7,5 cm

(with HP=44 cm)

Adjustable swing away Lateral supports 23

(tension adj. upholstery is compulsory)




- ☐ Right
☐ Left

Specify the height "h" from the center of the swing away to the seat
 cm

Head-rest 24

Padded




- ☐ Fixed
☐ Height adjustable

Orbital

- ☐ Anatomic



- ☐ Moulded



- ☐ With lateral supports



Specify "h" & "d" where requested
 Standard: h = 20 cm d = 15 cm
 h = cm d = cm

Table 25

- ☐ Transparent
☐ RxOLx Hemiplegic transp.
☐ Plastic grey
☐ Soft padded

Specify the fastening:
☐ single ☐ double



(Only with arm-rest)

Lined cushion in foam polyurethane 28



- ☐ 3 cm high
☐ 5 cm high


Brakes 26

- ☐ Classic - short handle
☐ Classic - long handle
☐ Drum
☐ Right single brake
☐ Left single brake
☐ Drum+classic

Options

- 29 ☐ "M" velcro strip (hook) on seat
 30 ☐ "M" velcro strip for Vicair
 31 ☐ "F" velcro strip (loop) on seat
 32 ☐ Seat/back-rest connection flap
 33 ☐ Chest harness
 34 ☐ Chest belt
 35 ☐ Abductor belt
 36 ☐ Front protective covers
 37 ☐ Stretched bar
 38 ☐ Catheter holder pocket
 39 ☐ R ☐ L Crutch hold./single
 40 ☐ R ☐ L Crutch hold./double

Abductor 27



- ☐ 1 ☐ 3 ☐ 5 ☐ 7 cm
☐ cm

Choose or specify the height from seat cover

- 41 ☐ Pump for pneumatics
 42 ☐ Suitcase holder
 43 ☐ 45° pelvic belt
 44 ☐ 45° padded pelvic belt
 45 ☐ Black backpack
 46 ☐ IV bar + poles
 47 ☐ IV poles
 48 ☐ Padded rigid cushion (L52max)
 49 ☐ Padded rigid cushion (L56-60)
 50 ☐ Q.rel.axles Ø12mm w/42mm bushes (def.)
 51 ☐ Q.rel.a. Ø1/2" (12,7mm) w/42mm bushes
 52 ☐ Q.rel.a. Ø1/2" (12,7mm) w/50mm bushes
 53 ☐ Stickers for car fastenings
 (escl.: load>120Kg, R12", UltraL.RW, hinged bckr)